




Product Manual
I3 Control Box Membrane Chromatography



| | |
|---------------------------|-------------------------|
| Customer | i3 Membrane GmbH |
| Product name | i3 Membrane Control Box |
| Item Number Elektrosil | MU-05.0008 (C sample) |
| Project Number Elektrosil | 22-00111 |
| Date of release | 2025-01-28 |
| Version of document | 04 |
| Customer | i3 Membrane GmbH |

| Rev. | Date | Author | Description |
|------|------------|--------|---|
| 04 | 2025-01-28 | JBO | Added: application area, altitude, relative humidity, overvoltage category, pollution degree of environment, suitability for damp locations |
| 03 | 2024-11-20 | YL | Adding technical drawing |
| 02 | 2024-10-28 | YL | Updating screen shots for graphical user interface |
| 01 | 2024-10-17 | YL | Document creation |

Content

| | | |
|----|--------------------------------|---|
| 1. | PRODUCT SAFETY | 3 |
| 2. | FUNCTIONAL DESCRIPTION | 3 |
| 3. | ELECTRICAL PROPERTIES | 3 |
| 4. | MECHANISCHE EIGENSCHAFTEN | 4 |
| 5. | FEATURES | 5 |
| | BUTTON | 5 |
| | LED-INDICATOR | 5 |
| | GRAPHICAL USER INTERFACE (GUI) | 5 |
| | MICRO SD-CARD | 7 |
| | CONNECTOR | 7 |

1. Product safety

1. Elektrosil assumes no warranty for this product if it is used under conditions other than those specified in the specification.
2. Please handle the device with care. It can be damaged by hard impacts.
3. Connection of electricity beyond the defined interfaces, dust, water, dampness, or other erosion elements can lead to safety issues or product failure.
4. Please ensure that the hardware is stored according to the specified storage temperature.
5. Please ensure that the hardware is used according to the specified operation temperature.
6. Do not turn on the device directly if it was stored in a temperature condition beyond the range for operation.
7. Please ensure that the hardware is only operated indoors in the laboratory environment (Pollution Degree 2).
8. For power supply, only use the power supply defined by Elektrosil, such as EA1019HVR5(T09)
9. For operating the device, only use the connection cables delivered by Elektrosil
10. To remove the power supply: Do not pull on the cable, always pull on the plug.
11. Do not open the casing. Opening the casing voids any warranty claimed regarding the i3 control box.
12. Do not clean the device with a corrosive cleaning agent.

2. Functional description

The i3 Membrane Control Box is used to control the applied voltages (max. +/- 5V) to the membrane electrodes. Voltage profiles are set via a 7-inch touch display. The voltage sequence, voltage/current measurements, and time are graphically displayed in real-time.

For operations of the control box, a micro SD card is required, where the measurement data is stored.

The electrical supply of the i3 Membrane Control Box is provided by a separate power supply unit (EA1019HVR5(T09)).

The control box is equipped with a button that allows controlled shutdown and startup of the device.


The power supply must only be disconnected from power when the device is shut down.

3. Electrical Properties

| Nr. | Designation | Specification MU-05.0007 |
|--------|--|--------------------------|
| 3 – 1 | Rated voltage | DC 12V |
| 3 – 2 | Operating voltage | DC 12V to 12.4V |
| 3 – 3 | Operating current max. | 2 A |
| 3 – 4 | Total input power i3 Control Box max. | 24 W |
| 3 – 5 | Permissible ambient temperature during operation | 0°C to 40°C |
| 3 – 6 | Permissible storage temperature | -10°C to 50°C |
| 3 – 7 | Permissible relative humidity | 15% to 90% |
| 3 – 8 | Overvoltage category | 1 |
| 3 – 9 | Operating altitude | -30,5m to 2000m |
| 3 – 10 | Storage and shipping altitude | -30,5m to 10.000m |

4. Mechanical Properties

| Nr. | Designation | Specification MU-05.0008 |
|-------|-----------------------------|---|
| 4 – 1 | Dimensions | Control Box Height: 68 mm Length: 204 mm Width: 180 mm |
| 4 – 2 | Weight | Control Box 1140g Power supply in the box 150g Cable 140g Sum: 1430g |
| 4 – 3 | Colour of casing | silver |
| 4 – 4 | Material of casing | Extruded profile: Al Mg Si 0.5 anodized; Lids: Al Mg Si 0.5 anodized |
| 4 – 5 | Power Socket i3 Control Box | Socket 2/5,5mm |



5. Features




Button

The control box will automatically start up once connected to the power. Pressing the button (for at least approx. 1s) will initiate a controlled shutdown of the device.

LED-Indicator

Once the control box starts up, the LED indicates green. After the controlled shutdown, the LED switches to red.

Graphical User Interface (GUI)

| | |
|---|---|
|  | <p>The graphical user interface (Version 3.2.1) of the control box starts automatically after the device starting up. The start up display indicates to connect the cable with i3 Membrane Adsorber. After connecting the adsorber. Tipp "Continue" to proceed.</p> |
|  | <p>After a successful resistance check, the process can be continued.</p> |
|  | <p>Following options are available in the overview menu:</p> <ul style="list-style-type: none"> • „Run Check“: Verify the electrical resistance of the membrane • „Settings“: Change date and time in the user interface. (After rebooting the GUI-time will be reset to system time) • „Single Voltage“: Enable a single DC voltage at the output. • „Presets“: Configure the pre-defined voltage profiles. • „Custom Sequence“: Specified voltage sequence by user up to 99 data points. • „Journal“: Storage of the last 200 enabled voltage sequences. |
| | <p>In order to enable a single DC voltage the voltage level and duration are required.</p> |

| | |
|---|--|
| <p>Single Voltage</p> <p>Voltage -1.96 V</p> <p>Duration 0:51 min</p> <p>Duration [mm:ss]</p> <p>Preview</p> <p>Start</p> <p>10:25</p> | |
| <p>Presets</p> <p>Square Wave</p> <p>Triangle</p> <p>Sawtooth</p> <p>Descending Triangle</p> <p>10:26</p> | <p>The menu „Presets“ contains four different voltage profiles:</p> <ul style="list-style-type: none"> • The square waveform • The triangle waveform • The sawtooth waveform • The descending triangle waveform |
| <p>Square Wave</p> <p>Voltage 1.48 V</p> <p>Duration 1:00 min</p> <p>Delay 0:20 min</p> <p>Start Delay 0:05 min</p> <p>Cycles 3 x</p> <p>Repetitions 3</p> <p>Custom Continuous</p> <p>Preview</p> <p>Start</p> <p>10:27</p> | <p>Following parameters must be set for each of the four waveforms:</p> <ul style="list-style-type: none"> • „Voltage“: voltage level • „Duration“: Duration of the enabled voltage output. (On Time) • „Delay“: Duration of the disabled voltage output. (Off Time) • „Start Delay“: Delayed time, untill the first voltage output in the first cycle is enabled. <p>Cycles: Number of the repetition</p> |
| <p>Sequence 1</p> <p>Cycles 5 x</p> <p>① 0:07 min 0.21 V</p> <p>Δ Time [mm:ss] 0:07</p> <p>Voltage [V] 0.21</p> <p>Slope Increment</p> <p>Remove Add</p> <p>Delete</p> <p>Preview</p> <p>Save</p> <p>10:27</p> | <p>In „Custom Sequence“ a specific voltage profile with maximal 99 data points can be configured. Each data point contains the parameters such as voltage level, a delta-time and the type of the transition corresponding to the previous data point.</p> |

Journal

< Back

Sequence 1

executed 2024/11/17 10:28

0.3

0.2

0.1

0

-0.1

-0.2

-0.3

0.00

0.01

0.02

0.03

0.04

0.05

0.06

0.07

0.08

0.09

0.10

0.11

0.12

0.13

[min]

Edit

Start

Square Wave

executed 2024/11/09 13:24

Single Voltage

executed 2024/11/09 13:23

10 30


The submenu „Journal” provides an overview of the enabled voltage sequences. With this the stored sequence can be reactivated again.

Micro SD-Card

With the current status of the samples, a micro-SD card is provided for storing the measurement data and is already inserted at delivery. The SD card must always be inserted before powering on or starting up the control box and must not be removed throughout the entire active controlling process.

Connector

To control voltage on the membrane, the control box has an eight-pole socket. For a secure connection, the nut should be tightened solidly. The following illustration and table contain the defined pin assignment.



| Pin out to connector | Pin out on the membrane |
|----------------------|-------------------------|
| 1; 7 | Female electrode 2 |
| 2; 8 | Male electrode 1 |
| 3; 4 | Male electrode 2 |
| 5; 6 | Female electrode 1 |