

## Download for Apple users

Step 1: Find the APP Store



Step 2: Search for "Intelligent life" and download the APP



Step 3: Open the "Intelligent life" APP



Step 4: Enter the APP interface



## Connection Flow

Long press the power button/Open the ear picking stick



Open "APP" and find "Sweep"



Scan the QR code in the bottom left corner of the box



Select Wifi Channel



## Download for Android users

Step 1: Find the Google Play



Step 2: Search "Intelligent life" APP



Step 3: Download "Intelligent life" APP



Step 4: Open the APP



## Connection Flow

Connect device "SmartLife-XXXX"




Open the APP interface



If you encounter switching WIFI, please do not switch WIFI continue to use



Function button



## Install ear scoop cover



Installation of ear spoon cover: insert the silicone ear spoon cover into the integrated ear spoon and install it in place before use.

## Product parameters

Product Name:	Visual otoscope
Product Model:	CT01
Image Sensor:	CMOS
Frame rate of image transmission:	25Fps
WIFI frequency:	2.4GHz
Network standard:	IEEE 802.11b/g/n
On-off mode:	Key switch
Battery Capacity:	600mAh
Duration of use:	About 50min
Charging time:	About 30min
Incoming current:	5V 2.0A
Charging interface:	Type-C
Appearance material:	PC+ABS
Working temperature:	-10~50 C
Lens enhancement:	6 warm LED lights

## Precautions

- 1、 When cleaning the lens, please wipe it carefully with professional Alcohol cotton Swabs
- 2、 Be sure to rotate the ear spoon to the top before use to prevent the ear spoon from falling
- 3、 Please use this product in an environment where no one is running or fighting, to avoid impact
- 4、 Please do not put the product in liquid or water to avoid damage
- 5、 It is forbidden for Children to use this product alone in order to avoid accidental harm
- 6、 This product has built-in rechargeable lithium battery, if not used frequently, Must also be recharged once a month
- 7、 Avoid product exposure to sunshine directly, especially the accessories, in order to avoid softening
- 8、 During using it, the temperature of the product will rise slightly (up to 37°C), Please feel free to use
- 9、 The ear scoop accessories should not be used for a long time. To avoid any breakages caused by aging, it is suggested that the ear spoon be replaced every 2 months
- 10、 Ear spoons are used as an auxiliary tool for ear cleaning, provided the health of the ear canal is fully noticed



## FCC Statement

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions: 1) this device may not cause harmful interference, and 2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

### FCC RF Radiation Exposure Statement

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.