# LU&BA

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

RESIN NANO COLLAGEN SUPERCONDUCTOR

## STATEMENT

Thank you for purchasing our Nano Resin Microneedling Pen. Before use, please carefully read this instruction and follow the precautions for safety herein to prevent that the equipment is damaged or the apparatus cannot be normally used.

LU&BA owns the copyright of all contents herein. Without the written permission of us, anyone is not allowed to copy, duplicate or translate this instruction to other languages.

Proprietary information protected by copyright law contained herein, shall be kept confidential by users and shall not be disclosed unrelated third parties.

The Company is entitled to revise the instruction and change the product description at any time, and the final interpretation right is owned by the Company.

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## FOREWORD

Thank you for purchasing this product.

This product is carefully made by us, combining unique technologies and designs, and meeting high quality standards.

We are proud to meet the demands of our customers as much as possible. We welcome your comments and suggestions and believe that keeping in touch with our customers is very important for our future product line. We encourage you to read through the instructions to fully understand the operating characteristics of the instrument.

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## INSTRUCTION

- 1) Please read this instruction carefully before use;
- 2) Before use of this apparatus, please clean skin, wash hands and disinfect;
- 3) During charging, DO NOT power on for operation or touch the charger with wet hands;
- 4) During placement of apparatus, it is required to remove the needles;
- 5) After use of the apparatus, please wipe and place the apparatus in a dry and ventilated place;
- 6) DO NOT open the device yourself, otherwise LU&BA will not be responsible for subsequent operation and any warranty claims.

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
## IoT STARTUP PROCESS

- 1) Mobile phone download APP- LU&BA
- 2) Enter Mobile phone number to log in
- 3) Add the customer's name in the upper left corner
- 4) Enter the operation page - Select the corresponding instrument LU&BA
- 5) Press the Start button to scan the product consumable code
- 6) One customer, one yard

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## CAUTIONS

- 1) The nano resin needles are disposable, secondary use is not allowed to prevent skin problems caused by cross infection.
- 2) Please place this apparatus in the places out of the reach of children.

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## OPERATING INSTRUCTIONS

- 1) Install the needles to the device and adjust the length of the needles;
- 2) Clean skin and apply the selected primary liquid on skin;
- 3) Long press to power on, short press to adjust the gears;
- 4) After finish of the operation, long press to power off and remove the needles.

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WARRANTY

The warranty is valid within two years from the date of purchase if the product has material, manufacturing or functional defects not caused by the user.

The warranty is invalid in the following cases:


- \* The device is damaged due to improper operation or modification or maintenance by unauthorized third parties.
- \* The device is damaged caused by transportation, such as impact or improper handling.
- \* Damage caused by moisture or other liquids, foreign objects entering the device or the use of unapproved and discouraged accessories contrary to the instruction.
- \* Use of accessories which do not comply with or match with the contents wherein during use of the device.
- \* Damage or loss of partial or whole device due to force majeure or accidents caused by natural disasters (such as earthquake, flood, fire and lightning, etc.).
- \* Surface damage, including but not limited to scratches and dents. For the above situations, if maintenance is needed, the users are required to bear corresponding expenses.

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This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:  
(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
  - Increase the separation between the equipment and receiver.
  - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
  - Consult the dealer or an experienced radio/TV technician for help.
- Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.