PHONAK

Important Safety Information Declaration of Conformity

Devices (options 02 and 03):

Roger MyLink, Roger X, Roger 10, Roger 11, Roger 13, Roger 14, Roger 15, Roger 16, Roger 17, Roger 18, Roger 19, Roger 20

Important Safety Information

 Please read the information on the following pages before using your Phonak device.

 This symbol on the product or its packaging means that it should not be disposed of with your other household waste. It is your responsibility to dispose of your waste equipment separately from the municipal waste stream. The correct disposal of your old appliance will help prevent potential negative consequences for the environment and human health.

Hazard warnings

-  Keep this device out of reach of children aged under 3 years.
-  Opening the device might damage it. If problems occur that cannot be resolved by following the guidelines in the Troubleshooting section of this User Guide, consult your hearing care professional.
-  Your device should only be repaired by an authorized service center. Changes or modifications to this device that were not explicitly approved by Phonak are not permitted.
-  Please dispose of electrical components in accordance with your local regulations.
-  Do not use your device in areas where electronic equipment is prohibited. If in doubt, please ask the person in charge.
-  Device should not be used on aircrafts unless specifically permitted by flight personnel.
-  External devices should only be connected if they have been tested in accordance with corresponding IEC standards.
-  Only use accessories that are approved by Phonak AG.

Product safety information

-  Protect the device from excessive moisture (bathing, swimming), heat (radiator, car dashboard), direct contact with skin when sweating (workout, fitness, sport) and from excessive shock and vibration.
-  X-ray radiation, CT or MRI scans may destroy or adversely affect the correct functioning of the device.
-  Never use a microwave or other heating devices to dry the device.
-  Do not expose your device to temperatures and humidity outside of the recommended conditions stated in the technical datasheet.
-  Clean the device using a damp cloth. Never use household cleaning products (washing powder, soap, etc.) or alcohol to clean the device.
-  Avoid exposing the receiver to hair spray or other cosmetic chemicals.

Other important information

- The digital signal sent from the device to a connected receiver cannot be overheard by other devices which are not in the microphone's network.
- Be aware that wireless transmission from microphone to receiver may be interrupted due to interference or moving out of range, and that words which are then spoken by the microphone wearer may not be heard by the receiver wearer.
- The sensitivity of the hearing aid microphone may be reduced by the use of your device and therefore environmental awareness may be reduced.
- Roger X may sound a loud beep at start-up, which might be annoying for cochlear implant users. Please contact your hearing care professional to discuss deactivating this beep function.

Declaration of Conformity

Roger receivers are products of Phonak AG, which have been tested and confirmed to be in compliance with the essential requirements and other relevant provisions of the Medical Devices Directive 93/42/EEC and Directive 1999/5/EC. To obtain a copy of the Declaration of conformity, please contact the manufacturer or the local Phonak representative whose address can be taken from the list on www.phonak.com (world-wide locations).

This device complies with Part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. The device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. The device generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the device does cause harmful interference to radio or tele vision reception, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the device and receiver.
- Connect the device into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This Class B digital apparatus complies with Canadian ICES-003.

FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the device.

Japanese Radio Law and Japanese Telecommunications Business Law Compliance.

This device is granted pursuant to the Japanese Radio Law (電波法)

This device should not be modified (otherwise the granted designation number will become invalid).

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- ⚠ N'utilisez que des accessoires approuvés par Phonak AG.