WiTouch[®] Pro

User Manual

Before Using Your WiTouch Pro Device

Sync the Remote Control and the WiTouch Pro Device

1. Using the provided screwdriver, remove the back cover from the WiTouch Pro device.
 2. Remove the clear plastic tab from the remote control.
 3. On the remote control, press and hold the (+) and (-) buttons simultaneously until the LED appears solid green (~ 10 seconds).
 4. Quickly insert the batteries into the WiTouch Pro device.
- Note:** The LED on the remote control will flash several times indicating syncing is complete.
5. Replace the back cover on the device.

Note: You may need to re-sync your remote control and WiTouch Pro device when you replace the battery in the remote control. This will ensure communication between your remote control and WiTouch Pro device. If your remote control and WiTouch Pro device ever lose communication, please repeat these syncing instructions.

Hollywog

2830 Amnicola Hwy

Chattanooga, TN 37406 USA

+1 423 305 7778

info@hollywog.com

www.hollywog.com/witouchpro


® and © 2015 Hollwog, LLC. All rights reserved.

U.S. Patents 8,972,016, D712,052, D701,610,

D716,958. Other U.S. & Foreign Patents Pending

11S.1004J 09/30/2015

Introduction

Congratulations on your purchase of the WiTouch Pro device. WiTouch Pro is a unique wireless remote controlled pain relief device incorporating TENS technology to specifically target back pain. The thin and flexible design perfectly contours the back for maximum surface contact. The advanced electronics design maximizes energy use, providing over 150 30-minute treatment sessions per battery life.

This innovative device is safe, drug-free, easy to use, discreet and comfortable to wear, and most importantly allows you to control your pain to maintain an active lifestyle.

Indications for Use: For the symptomatic relief and management of chronic intractable back pain and relief of pain of the upper and lower back associated with arthritis. It is also used for adjunctive treatment for post-surgical and post-trauma acute back pain.

Safety

CONTRAINDICATIONS

- **Do not** use this device if you have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device. Such use could cause electric shock, burns, or electrical interference or death.

WARNING!

- The device may cause rhythmic disturbances to the heart. **Do not** use the device across or through your chest. If you are susceptible to rhythm disturbances of the heart, use of the device must be done under the direction of a physician.
- Use of this device over your neck could cause muscle spasms resulting in airway closure, difficulty in breathing, adverse effects on heart rhythm or blood pressure. **Do not** use this device over your neck
- Use of the device when you are in the care of a physician or have had medical or physical treatment for your pain. **Consult** with your physician before using this device.
- Continued use of the device when pain does not improve, becomes more severe, or lasts more than five days may indicate a severe condition. **Stop** using the device and consult with your physician.
- Use of the device on the following skin conditions may cause a condition to become worse. **Do not** use the device over, or in proximity to, these skin conditions: abnormal skin, skin that is not intact, uncleaned, unhealthy, open wounds, rashes, swollen, red, infected, inflamed areas, skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins), or cancerous lesions.
- Electrical stimulation during common activities may increase the risk of injury. **Do not** use the device when in the bath or shower, sleeping, driving, operating of machinery or any activity in which electrical stimulation can put you at risk of injury.
- Using the device around electronic monitoring equipment (e.g., cardiac monitors, ECG alarms) may cause equipment malfunction. **Do not** use this device around electronic monitoring equipment.
- The effect, and safety, of using the device on children, during pregnancy, or use across the head has not been evaluated or established and is unknown. **Do not** use the device on children or let children handle the device. **Do not** use this device if you are pregnant, or suspect that you are pregnant, unless under the direction of your physician. **Do not** apply the device across your head.

PRECAUTION

- Using the device when you have suspected or diagnosed epilepsy or heart disease may cause unexpected reactions. **Always** consult your physician before using the device.
- Modifications not expressly approved by the manufacturer may void a user's authority to operate the device.

Contents

Not shown: User manual.

WiTouch Pro device (1)

Remote control with (1) CR2032 battery (comes installed)

AAA alkaline batteries (2)


Gel pads (1 pack)

#0 Phillips head screwdriver (1)

Features and Functions

1. **Electrode Area:** Surface where the gel pads are placed
2. **On/Off Button:** Press to turn the device On/Off
3. **Start/Stop Button:** Press to start or stop treatment
4. **+/- Buttons:** Press to increase or decrease stimulation intensity



Preparation for Use

Before using your WiTouch Pro you will need to sync the remote control, apply the gel pads to the device, and prepare the device for a treatment.

Applying the Gel Pads to the WiTouch Pro Device

The gel pads are intended for single person use. They will last, depending on skin type, oils, and pH levels, approximately two to five applications. Replace the gel pads when they no longer adhere completely. Follow these steps to apply the gel pads:

1. Separate two gel pads.
2. Remove the blue liner from the side being applied to the electrode area. Do not remove the green protective liner.
3. Align the shape of the first gel pad with the electrode area. Apply the gel pad onto the electrode area and firmly press across the entire surface to ensure good adhesion.
4. Repeat steps 1 – 3 for the second gel pad.

Skin Preparation

- Trim, not shave, excessive hair on the treatment area.
- Wash the skin and dry completely.

Use

The WiTouch Pro exclusive, patent pending 3-stage waveform incorporates both clinical theories of TENS to provide pain relief. The 30-minute stimulation treatment is delivered as follows:

Stage 1: 5 minutes of high frequency stimulation that initiates a feeling of pain relief by suppressing the transmission of pain signals in nerves. This stage provides a high sensory sensation and allows you to establish a comfortable intensity setting for the entire treatment.

Stage 2: 20 minutes of low frequency stimulation that initiates an increased endorphin release in the body to reduce the sensitivity to pain for an extended period of time following the 30-minute treatment. This stage provides a low sensory sensation often described as a gentle tapping sensation.

Stage 3: 5 minutes of high frequency stimulation, providing the same high sensory sensation experienced in Stage 1. This stage allows you to maintain the feeling of pain relief and complete the overall treatment with a comfortable high sensory sensation.


Treatment Recommendations

- You can leave the device in position for multiple treatments during the day. It will automatically turn Off after two (2) hours of inactivity.
- It is recommended you wait a minimum of 30 minutes between treatments.


Conducting a Treatment

Always read the safety warnings before conducting a treatment. Follow these steps to conduct a treatment:


1. Remove the green liners from the gel pads by slowly peeling the liner diagonally from an inside corner to the opposite outside corner. **Important!** Avoid contact of the gel pad with other objects. Contact with other objects may affect the pads adhesion properties. Save the green liners for storage of the device.

2. Press, and hold, the **On/Off** button on the WiTouch Pro device for one (1) second. The LED will begin flashing indicating the device is ready for use. **Note:** The WiTouch Pro device must always be ON prior to application.

3. Align the center of the device over the spine and place the device on your back in the area of pain. If you cannot place the device properly, ask another person for assistance.


Important! Do not apply the gel pads/electrodes directly over the spine.

- Press the **Start/Stop** button on the remote control to begin the treatment.

Note: The LED on the WiTouch Pro device will flash rapidly during treatment.

- Press the (+/-) buttons to increase or decrease the intensity of the stimulation until it is at a comfortable level. **Note:** The stimulation intensity can be controlled on the WiTouch Pro device with either the remote control or the buttons located on the device. This allows for adjustment by either you or your healthcare professional.

- Press the **Start/Stop** button on the remote control to stop the treatment at any time. Treatment will automatically stop after 30 minutes. **Note:** If a treatment is stopped and restarted the treatment will restart at the first stage.

Removing the WiTouch Pro

Important! Do not remove the device until the treatment has stopped.

- After treatment, or when you want to remove the device, grasp the edge of the device and gel pad to ensure the gel pad does not stay on the skin. Slowly peel the device away from the skin.
- Align, and place the green protective liners back on the gel pads. Ensure the pads are completely covered.

Battery Replacement


When battery replacement is needed the LED will flash as follows:

Item	Description	Battery
WiTouch Pro device	LED flashes yellow every two (2) seconds when the unit is On	(2) AAA alkaline batteries
	LED flashes yellow rapidly during treatment	
Remote Control	LED will flash yellow once when the Start/Stop button is pressed or no LED flash at all.	(1) CR2032 lithium battery

Note: Fully depleted batteries will have no flashing LED.

To replace the WiTouch Pro Device Batteries:

- Using the included #0 Phillips head screwdriver remove the cover.
- Remove the old batteries, and place the new batteries in the correct polarity.
- Replace the battery compartment cover.

To replace the remote control battery:

- Place the provided screwdriver, coin or flat tool in the groove on the side of the remote control. Turn the coin and lift the cover off.
- Slide the old battery out, and slide a new CR2032 battery into the slot. **Note:** The (+) is facing upward.
- Align the cover with the remote control base. Holding the handle portion of the remote control, snap the cover down into place.

Note: Consult your local authorities for proper battery disposal.

Maintenance and Storage

- Clean the device using a damp cloth and mild soap. Gently wipe clean. **Note:** This device is manufactured with water detection technology. Do not submerge the device in water or other liquids. Water damage to the electronics will void the warranty.
- The device should be operated, transported and stored at temperatures between 50° F and 104° F (10° C and 40° C), with relative humidity between 30% - 85%. All values have +/- 10% tolerance.
- Store the device, remote control and gel pads in the original packaging when not in use.
- If the device is not working properly, stop using immediately. Do not disassemble or modify the device. Contact Hollywog Customer Service at: Phone: +1 423 305 7778
Email: info@hollywog.com.

Technical Specifications

Channels: Single channel

Waveform: Asymmetric biphasic square pulse

Pulse Amplitude: 0 ~110 mA = 0 ~ 55 volts, adjustable (at 500 ohm load)

Pulse Frequency: (Hz) 5-120

Pulse Width: (μ s) 120 - 240

Timer Control: (mins) 30

Power Supply: Device: Two (2) AAA batteries (Internally Powered)



Remote Control: One (1) CR2032 Lithium battery (Internally Powered)

Size (D x W x H): 0.7" x 7.5" x 3.5" (18 mm x 191 mm x 90 mm)

Weight (including battery): 4.8 oz (136 g)

Safety Standards: IEC/EN/UL 60601-1, IEC/EN 60601-1-2, IEC 60601-2-10, CAN/CSA C22.2 No. 601.1, EN 300200 / EN 301498-3

Made in the USA

Limitations of Use and Storage					
 <p>500 hPa ↔ 1060 hPa</p> <p>Atmospheric Pressure Limitations</p>	 <p>50°F / 10°C ↔ 104°F / 40°C</p> <p>Temperature Limitations</p>	 <p>30% ↔ 85%</p> <p>Humidity Limitations</p>			
 <p>CE 0120 European Conformity</p>	 <p>European Directive 2002/96/EC</p>	 <p>Type BF Applied Part</p>	<table border="1"><tr><td>EC</td><td>REP</td></tr></table> <p>MDSS GmbH Schiffgraben 41 30175 Hannover, Germany</p>	EC	REP
EC	REP				
<p>RF: 2.4GHz transceiver</p> <p>Device: FCC ID: N8L-HWOGTENS11 IC: 10346A-HWOGTENS11</p> <p>Remote Control: FCC ID: N8L-HWOGTENS12 IC: 10346A-HWOGTENS12</p>					

Warranty

COMPANY SHALL NOT BE LIABLE IN ANY EVENT FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES.

Hollywog, LLC ("Company") warrants the WiTouch Pro ("Product") is free of defects in materials and workmanship. This warranty shall remain in effect for 1 year (12 months) from the date of purchase. If these Products fail to function during the 1 year warranty period due to a defect in material or workmanship, at the Company's Option, Company or the selling retailer will replace the respective Product without charge within a period of thirty (30) days from the date on which the Product is returned to the Company or the retailer.

This Warranty Does Not Cover: WiTouch gel pads. Any defects, damage, loss, malfunction or failure in the Product caused by product misuse, alteration, disassembly or any use that is inconsistent with the Product's Instruction Manual.

Guidance and manufacturer's declaration –electromagnetic emissions

The WiTouch® and WiTouch® Pro are intended for use in the electromagnetic environment specified below. The customer or the user of the WiTouch or WiTouch Pro should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The WiTouch and WiTouch Pro use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The WiTouch and WiTouch Pro are suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration –electromagnetic immunity

The WiTouch and WiTouch Pro are intended for use in the electromagnetic environment specified below. The customer or the user of the WiTouch or WiTouch Pro should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the [EQUIPMENT or SYSTEM] requires continued operation during power mains interruptions, it is recommended that the [EQUIPMENT or SYSTEM] be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A / m	3 A / m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration –electromagnetic immunity

The WiTouch® and WiTouch® Pro are intended for use in the electromagnetic environment specified below. The customer or the user of the WiTouch or WiTouch Pro should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	Portable and mobile RF communications equipment should be used no closer to any part of the WiTouch or WiTouch Pro, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.17\sqrt{P}$ 150 kHz to 80 MHz $d = 1.17\sqrt{P}$ 80 MHz to 800 MHz $d = 2.33\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, "should be less than the compliance level in each frequency range." Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the WiTouch or WiTouch Pro is used exceeds the applicable RF compliance level above, the WiTouch or WiTouch Pro should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the WiTouch or WiTouch Pro.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the WiTouch and WiTouch Pro

The WiTouch and WiTouch Pro are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the WiTouch or WiTouch Pro can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the WiTouch or WiTouch Pro as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.70	11.70	23.30

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

This device complies with Industry Canada's license-exempt RSSs. Operation is subject to the following two conditions:

- (1) This device may not cause interference, and
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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


Hollywog
2830 Amnicola Hwy
Chattanooga, TN 37406 USA
+1 423 305 7778
info@hollywog.com
www.hollywog.com/witouchpro
® and © 2015 Hollwog, LLC. All rights reserved.
U.S. Patents 8,972,016, D712,052, D701,610,
D716,958. Other U.S. & Foreign Patents Pending

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