




USER MANUAL

smaXtec® System

© 2014 smaXtec animal care sales GmbH. All rights reserved.

smaXtec® is a trademark of smaXtec animal care sales GmbH. This trademark may not be used without permission of smaXtec animal care sales GmbH.

This user manual may not be copied. Changes regarding this document will not be announced. Before using the smaXtec® System, please check if there is a new version of this document available on the website: <http://www.smaxtec-animalcare.com>

The smaXtec pH & Temp Sensor (part of the smaXtec® System) is under protection of the following patents:

WO 2011/079338 A2 (PCT/AT2010/000490)
EP 2 519 098 (EP 10808965.7) (Pending)
AT 509 255 B1 (A 2052/2009)
BR 112012016375-2 (Pending)
CA 2 785 925 A1 (Pending)
RUS 2012 132 453 (Pending)
UA 2012/09263 (Pending)
US 2012/0277550 A1 (13/520,101) (Pending)
ZA 2012/04642

Design Patent:

"Messvorrichtung, insbesondere zur Messung von physiologischen Daten im Körper eines Lebewesens"

"Measuring Instrument"

European Union: RCD 001702598-0001 - Date: 30.04.2010 (DD/MM/YYYY)

USA: US D655,633 S – Date: 19.10.2010 (DD/MM/YYYY)

DD/MM/YYYY

D = Day

M = Month

Y = Year

Contents

1. Overview of the smaXtec System.....	4
2. System Requirements	6
3. Getting Started	7
4. smaXtec Sensors.....	8
4.1. smaXtec pH & Temp Sensor and smaXtec Temp Sensor	8
4.1.1. First Steps.....	8
4.1.2. Notes - Application	8
4.2. smaXtec Climate Sensor	10
4.2.1. First Steps.....	10
5. smaXtec Read-Out Devices.....	11
5.1. smaXtec Base Station.....	11
5.1.1. Installation and Start Up.....	13
5.2. smaXtec Repeater	15
5.2.1. Installation and Start Up.....	16
6. smaXtec Software	18
6.1. smaXtec Sensor Init Software	18
6.1.1. First Steps.....	18
6.2. smaXtec Messenger.....	19
6.2.1. First Steps.....	19
7. Technical Details	20

1. Overview of the smaXtec System



WARNING: For your own safety, please read all instructions in the User Manual as well as the **Important Safety and Product Information** before using the smaXtec® System! Both documents are available for download at www.smaxtec-animalcare.com.

The **smaXtec Sensor Init software**, **smaXtec Messenger** and further information about the smaXtec System, including informational videos, are available in the smaXtec Support Center (<http://support.smaxtec-animalcare.com>) and on our website (<http://www.smaxtec-animalcare.com>).

The main components of the **smaXtec System** are:



1. The **smaXtec pH & Temp Sensor** is designed to provide continuous measurement of pH and temperature in the rumen. The data recorded in the smaXtec pH & Temp Sensor is automatically transmitted via the wireless network to the smaXtec Base Station or to a smaXtec Repeater. Due to the dimensions, use of the smaXtec pH & Temp Sensor is not possible for cows under the age of 18 months and weighing less than 450 kg.
2. The **smaXtec Temp Sensor** is designed to provide continuous measurement of rumen temperature in cows. The data recorded in the smaXtec Temp Sensor is automatically transmitted via the wireless network to the smaXtec Base Station or to a smaXtec Repeater. Due to the dimensions, use of the smaXtec Temp Sensor is also not possible for cows under the age of 18 months and weighing less than 450 kg.
3. The **smaXtec Climate Sensor** is designed to provide continuous measurement of environmental factors such as ambient temperature and humidity levels. This data is automatically transmitted via the wireless network to the smaXtec Base Station or to a

smaXtec Repeater. The number of smaXtec Climate Sensors required is dependent on the size of the cow shed.

4. The **smaXtec Base Station** automatically receives the pH and temperature data in real-time. The data is transmitted to the server and stored in the online database. The smaXtec Base Station automatically receives the data from the various smaXtec sensors. smaXtec Repeaters can be used to extend the network range.
5. The **smaXtec Repeater** is designed to extend the range of the smaXtec Base Station and is recommended for use in certain configurations in order to cover remote areas of the cow shed. It provides an easy-to-use and cost-efficient extension of the smaXtec wireless network and communication with the smaXtec Sensors.
6. The **smaXtec Messenger** is an automatic alert messaging tool that keeps you informed about changes in a cow's state of health. Collection of data in real-time means you are kept fully informed at all times about your herd's state of health. The smaXtec Messenger provides you with up to the minute information confirming if rumen pH value and/or temperature or conditions in the cow shed are within or outside the normal range. The smaXtec Messenger can be accessed via PC/laptop, tablet or smartphone (Google Play Store, Apple App Store) (<http://messenger.smaxtec.com/> or via app). Before using the **smaXtec pH & Temp** and **smaXtec Temp Sensors** with the smaXtec System, these devices must first be initialised using the **smaXtec Sensor Init** software and then activated in smaXtec Messenger.

2. System Requirements

The following items are required for operation of the smaXtec System:

- a PC/laptop with a USB connection and one of the following operating systems:
 - Windows 7 Home, Professional
 - Windows 8 & Windows 8.1
- access to smaXtec Messenger via your browser (Internet Explorer, Firefox, Chrome etc.
- smaXtec Sensor Init software (internet connection required)
- smaXtec USB antenna
- smaXtec Base Station
- possibly a smaXtec Repeater (dependent on the application and local conditions in the cow shed)
- magnet (delivered with the smaXtec Base Station)
- applicator for smaXtec pH & Temp sensors and smaXtec Temp Sensors (inside diameter 35 mm – can be ordered directly from smaXtec).

3. Getting started

1. Activation of the devices in the smaXtec wireless network (smaXtec Base Station, smaXtec Repeater):

To start using the smaXtec System follow the steps described in the Quickstart Guide. Start with activation of the smaXtec Base Station and/or the smaXtec Repeater by going to <https://activate.smaxtec.com> or directly in smaXtec Messenger (<https://messenger.smaxtec.com>). Before activation can start you will first need to complete your online registration. Next, mount the smaXtec Base Station and/or the smaXtec Repeater in areas where the cows tend to spend extended periods of time (e.g. milking stations, waiting areas in front of milking stations, feeding stations, etc.). Please see Chapter 6.1 for further information.

2. Installation of the devices in the smaXtec wireless network (smaXtec Base Station, smaXtec Repeater):

Prior to installing the devices, please refer to the instructions in the Quickstart Guide. Mount the smaXtec Base Station and/or one/several smaXtec Repeater(s) in areas where the cows tend to spend extended periods of time (e.g. milking station, waiting areas in front of milking stations, feeding stations, etc.). For further information, please see Chapter 5.1 (smaXtec Base Station) and Chapter 5.2 (smaXtec Repeater).

3.a. Configuration of the smaXtec pH & Temp Sensors/smaXtec Temp Sensors:

Using the smaXtec Sensor Init software, which you can download from our website, you now need to initialise the smaXtec pH & Temp Sensors/smaXtec Temp Sensors. Please see Chapter 6.1 for further information.

3.b. Installation of the smaXtec Climate Sensors:

The smaXtec Sensor Init software is not required in order to use smaXtec Climate Sensors. To ensure proper use, please follow the steps described in the Quickstart Guides provided. For further information, please see Chapter 6.1.

4. Read-out of measurement data

Following correct installation of the smaXtec Base Station (and a smaXtec Repeater if required) measurement data is automatically read-out from the sensors when the cows are in range (read-out range is 5 – 10 Meter/16 – 32 feet).

5. Working with smaXtec Messenger

Start your browser, go to <https://messenger.smaxtec.com> and log in with your user name and password. Before using smaXtec Messenger for the first time you need to register. Please ensure that you enter a valid email address at registration.

4. smaXtec Sensors



Caution: The Quickstart Guide included with the products provides a brief summary of the most important information and steps regarding initialisation, activation and proper use of smaXtec products. Please ensure you follow these steps carefully in order to make proper use of the smaXtec devices.

4.1. smaXtec pH & Temp Sensor and smaXtec Temp Sensor

4.1.1. First Steps

Follow the steps described in the accompanying Quickstart Guide and ensure they are carried out correctly.

- 1) Configuration and switching on of the smaXtec pH & Temp Sensor/smaXtec Temp Sensor using the **smaXtec Sensor Init** software (download at <http://www.smaxtec-animalcare.com> and installation on your PC/laptop) – please note that to follow the steps for the **smaXtec Sensor Init** software, the PC/laptop must be connected to the internet.
- 2) Activation in **smaXtec Messenger** – menu item “Activate device”: Device activation is required to ensure the device can be correctly assigned to the user and used in smaXtec Messenger.
- 3) Choice of an appropriate cow and application of the sensor (please see relevant notes in Chapter 4.1.2)

Note:

To enable your device to communicate with the **smaXtec Sensor Init** software you need to use the smaXtec USB antenna, which is delivered with the smaXtec Base Station. The magnet required for switching on purposes is also delivered with smaXtec Base Station.

4.1.2. Notes - Application



Caution: The smaXtec pH & Temp Sensor and the smaXtec Temp Sensor are exclusively for use in cows. Under no circumstances may they be used for other animals or for humans. Use of the smaXtec pH & Temp Sensors and smaXtec Temp Sensors is not possible for cows under the age of 18 months weighing under 450 kg. If you have any concerns, please consult your veterinarian.

Cows that are to be fitted with a smaXtec pH & Temp Sensor or a smaXtec Temp Sensor must have an empty stomach when the sensor is applied. This means that the cow must not have eaten anything for at least 2 hours prior to application.

Applicator for the smaXtec pH & Temp Sensor and the smaXtec Temp Sensor:

The smaXtec pH & Temp Sensor or the smaXtec Temp Sensor can be inserted into the animal's maw using an applicator. It is recommended that you use a smaXtec applicator. This can be ordered from your distributor or directly from smaXtec animal care.

Application Procedure:

smaXtec pH & Temp Sensor

Position the smaXtec pH & Temp Sensor into the applicator with the sensor tip at the head. Insert the smaXtec pH & Temp Sensor into the applicator in such a way that sensor tip comes to rest at the far end of the applicator. Please ensure that the push rod of the applicator is kept flat – the push rod should not press directly on the sensor tip. This ensures that the smaXtec pH & Temp Sensor is positioned sensor tip first in the cow's stomach.

Sensor shaft



smaXtec Temp Sensor

In contrast to the smaXtec pH & Temp Sensor, the smaXtec Temp Sensor must be positioned with the sensor shaft (the bottom end) at the head of the applicator.



Picture: application procedure

With the cow's head and neck angled in an extended and straightened position you can now administer the smaXtec pH & Temp Sensor or the smaXtec Temp Sensor. The front end of the applicator must extend beyond the base of the tongue before the smaXtec pH & Temp Sensor or the smaXtec Temp Sensors is released. Please ensure that the cow does not eat anything for at least two hours prior to administration of the sensor. It is preferable to administer the smaXtec pH & Temp Sensor or the smaXtec Temp Sensor in the morning before the first feed.



Warning: After insertion the smaXtec pH & Temp Sensor or the smaXtec Temp Sensor needs to settle on the floor of the reticulum as quickly as possible. For this reason, it is important to ensure that the fibrous material in the rumen is not too dense, i.e. the cow should not recently have eaten anything. It is therefore of utmost importance that your animal should have an empty stomach when the sensor is administered. Acting in accordance with the user instructions provided by smaXtec will ensure smooth operation of the measurement procedure.

4.2. smaXtec Climate Sensor

The smaXtec Climate Sensor is designed to provide continuous measurement of environmental factors such as ambient temperature and humidity levels. This data is wirelessly transmitted to the smaXtec Base Station or a smaXtec Repeater. The number of smaXtec Climate Sensors required depends on the size of the cow shed.



Scope of delivery

The smaXtec Climate Sensor package contains the following components:

- 1x Climate Sensor
- 4x screws and 4x wall plugs for mounting
- 1x Drill template
- 1x Quickstart-Guide

4.2.1. First Steps

Follow the steps described in the Quickstart Guide provided and please ensure they are carried out correctly.

- 1) Activation in **smaXtec Messenger** – menu option “Add device”: It is necessary to activate the device so that the device can be correctly assigned to a user and used in smaXtec Messenger.
- 2) Switch on using the magnet provided (magnet is included in delivery with the smaXtec Base Station):
 - a. Switch on by holding the magnet in front of the On/Off symbol ()
 - b. Power LED () glows yellow for 1 second
 - c. Connection (takes approx. 1 second) to the smaXtec Base Station/Repeater (wireless network)
 - d. Connection OK (connected to the smaXtec Base Station/Repeater) – LED lights up green for 3 seconds
 - e. Connection NOT OK – LED flashes red 5 times: smaXtec Climate Sensor then automatically switches off → switch on again for new connection attempt


Note – TO SWITCH OFF:

The smaXtec Climate Sensor can be switched off using the magnet. Simply hold the magnet in front of the On/Off symbol again – LED glows red for 3 seconds.

5. smaXtec Read-Out Devices

5.1. smaXtec Base Station

The smaXtec Base Station automatically receives the pH & temperature and cow shed climate data. The data can be transmitted via the in-house network (Ethernet cable required) or via GSM (technical standard for mobile communications).



Scope of delivery

The smaXtec Base Station package contains the following components:

- 1x smaXtec Base Station
- 1x smaXtec USB antenna (required for initialisation of the smaXtec pH & Temp Sensor/smaXtec Temp Sensor)
- 1x magnet with tag (required for switching the sensors on)
- 1x mains adapter (with 3 m/10 ft cable)
- 4x screws and 4x wall plugs for mounting
- 1x Drill template
- 1x optional GSM antenna
- 1x Quickstart-Guide

Should any of these components be missing or damaged, please contact your smaXtec distributor or smaXtec. Please retain the box and the original packaging materials in case you need to return the product for any reason (e.g. repairs) in the future.

The smaXtec Base Station features the following status LEDs (top left):



Power LED:

- Indicates if the smaXtec Base Station is on.



Link LED

- Indicates the status of the connection to the server.



Data LED

- Indicates if data from smaXtec sensors is being read out.

On the right side of the smaXtec Base Station there is a GSM antenna and on the left you will find the socket for the network cable and the power supply. This is where the mains adapter (included in delivery) is connected during installation. The smaXtec Base Station must also be connected to the network via a network cable or an optional GSM antenna used.

smaXtec Network Cable (not included in delivery scope)

In locations where the smaXtec Base Station will be exposed to wet conditions or moisture, a special network cable is required for the connection between the smaXtec Base Station and the network socket. This ensures that the network cable remains waterproof and dustproof when plugged in. The cable must comply with norm TIA/EIA-568-B (with regard to length, bending radius etc.).



Image: Network cable (5/10/20 m | 16/32/65 ft)

The network cable can be ordered directly from smaXtec animal care or from your distributor in the following lengths: 5/10/20 m (16/32/65 ft).

Power Supply – included mains adapter

The smaXtec Base Station is powered via the mains adapter provided. Only use this mains adapter for installation purposes. At the end of the cable there is a specially designed plug to protect the smaXtec Base Station from water and dust. Connect the mains adapter to the socket on the left-hand side of the smaXtec Base Station and tighten the screw caps. The mains adapter itself does not conform to an IP standard with regard to protection against water, i.e. the mains adapter must not be exposed to rain or moisture.




Image: Mains adapter (3 m/10 ft cable)

Specifications of the adapter:

- Manufacturer: SUNNY Computer Technology Europe, s.r.o.
- Input: 100-240 V ~, 0.5 A MAX, 50-60 Hz
- Output: 5 V --- 2.0 A
- Output power: 10 W MAX.

5.1.1. Installation and Start Up



Caution: Please observe all safety instructions concerning electrical devices and installations.

Note:

Following correct installation of the smaXtec Base Station including the smaXtec Repeater as required (see Quickstart Guide and in this Chapter 5.2 in this User Manual), measurement data is automatically read out from the smaXtec Sensors whenever the cows are in range (read-out range is 5 – 10 Meter/ 16 – 32 feet).



Warning: If the smaXtec Base Station is mounted in the cow shed it is recommended that you use a special network cable to ensure proper functionality in conditions typically found in the cow shed (waterproofness of the smaXtec Base Station). This cable is available from smaXtec or from your distributor in the lengths 5/10/20 m (16/32/65 ft). It is, however, also possible to connect the smaXtec Base Station to your network with a standard network cable if the smaXtec Base Station is not located in the wet area/outdoor area of the cow shed.

The most important steps for successful installation of the smaXtec Base Station are described below:

1. Before mounting the smaXtec Base Station – Activation of the smaXtec Base Station:

Activation in **smaXtec Messenger** – menu item “Add devices”: device activation is required so that the device can be correctly assigned to a user and used in smaXtec Messenger.

2. Choose the most suitable location for the smaXtec Base Station

Before mounting the smaXtec Base Station, please select the most suitable location as described in the Quickstart Guide. The following points must be observed and implemented:

- Suitable position: Waiting area, milking area, feeding stations, massage brushes (areas where the cows tend to spend more time)
- Power and Ethernet connection or mobile communications reception required using the optional GSM.
- Power outlet may be up to a maximum of 3 m/10 ft away as the power supply unit / adapter cable provided is 3 m/10 ft long.
- Important for the Ethernet option:
 - Waterproof network cable for the connection to the smaXtec Base Station and the network socket are available from smaXtec or from your distributor in the lengths 5/10/20 m (16/32/65 ft).
 - If using a network cable you have installed yourself, please ensure a maximum length of 90 m/295 ft between the smaXtec Base Station and the network socket is not exceeded.
 - If using a network cable you have installed yourself, please ensure you only use network cabling with double shielding (SSTP, SFTP)
- Mount at a minimum height of 2 m/6 ft (not within reach of cows)
- Do not mount in a closed storeroom (transmission problems)
- Close to cows (the smaXtec Base Station can read out smaXtec Sensors within a radius of 5 – 10 m/16 – 32 ft; otherwise a smaXtec Repeater is required)

- Remember to consider locations for a smaXtec Repeater (if planned)
- An unobstructed view enables a greater range for planned smaXtec Repeaters
- Do not expose the smaXtec Base Station and mains adapter to direct sunlight or rain.

3. Mounting of the smaXtec Base Station:

Please follow the following steps:

- Use the screws and wall plugs provided
- Lift up plastic panel
- Drill mounting holes (use the drill template provided, see Quickstart Guide)
- Wood: screws can be used on their own (2 mm (1/12") holes can be pre-drilled)
- Walls: use wall plugs (drill a 6 mm (1/4") hole and insert wall plug)

4. Connecting the mains adapter and network cable/GSM antenna:

When the smaXtec Base Station has been mounted in accordance with instructions, please connect the mains adapter and the network cable/GSM antenna. Please follow the steps described below:

- To ensure dustproof and waterproof use of the smaXtec Base Station, please use the smaXtec network cable 5/10/20 m (16/32/65 ft).
- Connect the power supply unit provided with the smaXtec Base Station and plug in to a wall socket. If there is no wall socket within 3 m/10 ft then it is best to install a new one to the left of the smaXtec Base Station.
- Please note (applies to network cable and mains adapter cable):
 - It is not permitted to loop cables around the smaXtec Base Station (transmission problems).
 - If possible, the cable should lead to the smaXtec Base Station from above.

5. Verification of correct functionality with the help of status LEDs:

Following successful completion of all steps, the status of the smaXtec Base Station can be verified with the help of the LEDs. As soon as the mains adapter is plugged in the smaXtec Base Station starts up. Following successful installation, please check the status of the LEDs (connection can take up to 10 minutes):

- smaXtec Base Station has started up: Power LED glows steady green
- connected to the internet: Link LED glows steady green

The LEDs in detail:



Power LED:

- Flashing yellow: Starting up
- Glowing green: Started up/OK
- Flashing red/red: Error



Link LED:

- Flashing yellow: Connecting (depending on use of Ethernet or GSM this can take up to 10 minutes)
- Glowing green: Connected to the internet
- Flashing red/red: No connection to the server (via Ethernet/GSM)



Data LED:

- Flashing green: Data is being read out from the smaXtec Sensors

5.2. smaXtec Repeater

In certain configurations we recommend that you use a smaXtec Repeater (to ensure coverage of various areas of the cow shed). This enables an easy and low-cost extension of the smaXtec wireless network and communication with the smaXtec Sensors.



Scope of Delivery

The smaXtec Repeater packaging contains the following items:

- 1x smaXtec Repeater
- 1x mains adapter (3 m/10 ft cable length)
- 4x screws and 4x wall plugs for mounting
- 1x Drill template
- 1x Quickstart Guide

Should one of these items be missing or damaged, please contact your smaXtec distributor or contact smaXtec directly. Please retain the box and the original packaging material in case you need to return the product for repairs at a later date.

The smaXtec Repeater has the following status LEDs (top left):



Power LED:

- Indicates if the smaXtec Repeater is switched on.



Link LED

- Indicates the status of the connection in the wireless network (either to another smaXtec Repeater or to the smaXtec Base Station).



Data LED

- Indicates if data is being read out from smaXtec Sensors.

Power supply – mains adapter supplied

The smaXtec Repeater is powered by the mains adapter provided. For installation purposes please only use this mains adapter. At the end of the cable there is a specially constructed plug, which protects the smaXtec Repeater from water and dust. Plug the connector into the socket on the left-hand side of the smaXtec Repeater and tighten the screw cap. The mains adapter itself does not conform to an IP standard with regard to protection against water, i.e. the mains adapter must not be exposed to rain or moisture.




Image: Mains adapter (3 m/10 ft cable)

Specifications of the mains adapter:

- Manufacturer: SUNNY Computer Technology Europe, s.r.o.
- Input: 100-240 V ~, 0.5 A MAX, 50-60 Hz
- Output: 5 V --- 2.0 A
- Output power: 10 W MAX.

5.2.1. Installation and Activation



Caution: Please observe all safety instructions concerning electrical devices and installations.

Note:

Following correct installation of the smaXtec Repeater (see Quickstart Guide and Chapter 5.1 in this User Manual), the measurement data are automatically read out from the smaXtec Sensors whenever the cows are within range (read out range is 5 – 10 meters/16 – 32 feet).

The most important steps for successful installation of the smaXtec Repeater are described below:

1. Before installation and activation of the smaXtec Repeater:

Activation in **smaXtec Messenger** – menu item “Add device”: Device activation is required to ensure the device can be correctly assigned to the user and used in smaXtec Messenger.

2. Choose the most suitable location for the smaXtec Repeater

Before mounting the smaXtec Repeater, please select the most suitable location as described in the Quickstart Guide. The following points must be observed and implemented:

- Suitable position: Waiting area, milking area, feeding stations, massage brushes (areas in which the cows spend longer periods of time)
- Power outlet must be no further than 3 m/10 ft away as the mains adapter cable provided is 3 m/10 ft long.
- Mount at a minimum height of 2 m/6 ft (not within reach of cows)
- Do not mount in a closed storeroom (transmission problems)
- Position close to the cows (the smaXtec Repeater can read out smaXtec Sensors within a radius of 5 – 10 m/16 – 32 ft)

- An unobstructed view enables greater range for additional smaXtec Repeaters and/or to the smaXtec Base Station
- Do not expose the smaXtec Repeater and the mains adapter to direct sunlight or rain.

3. Mounting of the smaXtec Repeater:

Please follow the steps below:

- Please use the screws and wall plugs provided
- Lift up plastic panel
- Drill mounting holes (please use the drill template provided, see Quickstart Guide)
- Wood: screws are sufficient (pre-drill a 2 mm (1/12“) hole)
- Walls: Use wall plugs (6 mm (1/4“) drill a hole and push in wall plug)

4. Connection to mains adapter:

Following proper mounting of the smaXtec Repeater please connect the mains adapter. Please follow the steps below:

- Connect the power supply provided to the smaXtec Repeater and plug in. If there is no plug socket within 3 m/10 ft then it is best to install a new one to the left of the smaXtec Repeater.
- Please note:
 - It is not permitted to loop cables around the smaXtec Repeater (transmission problems).
 - If possible, the cable should lead to the smaXtec Repeater from above

5. Verification of correct functionality with the help of status LEDs:

Following successful completion of all steps, the status of the smaXtec Repeater can be verified with the help of the LEDs. As soon as the mains adapter is connected the smaXtec Repeater starts up. Following successful installation, please check the status of the LEDs:

- smaXtec Repeater has started up: Power LED glows steady green
- connected to the smaXtec Base Station: Link LED glows steady green

The LEDs in detail:



Power LED:

- Glowing green: Started up/OK
- Flashing red/red: Error



Link LED:

- Flashing yellow: Connecting
- Glowing green: Connected to the smaXtec Base Station (either direct or via a smaXtec Repeater)
- Flashing red/red: No connection to the smaXtec Base Station (either direct or via other smaXtec Repeaters)



Data LED:

- Flashing green: Data being read out from smaXtec Sensors

6. smaXtec Software

6.1. smaXtec Sensor Init software

To use the **smaXtec pH & Temp Sensor** and the **smaXtec Temp Sensor** these sensors must first be correctly initialised using the **smaXtec Sensor Init** software.

6.1.1. Initial steps



WARNING: Before using the magnets provided please consult the **Safety and Product Information** (available for download from <http://www.smaxtec-animalcare.com>)!

1. Downloading the smaXtec Sensor Init software:

Go to <http://www.smaxtec-animalcare.com>, click on the Download icon, download the **smaXtec Sensor Init** software and install the software. Should you have any problems with the driver installation please go to <http://support.smaxtec-animalcare.com> or follow the steps in the Driver Handbook for the smaXtec USB antenna.

2. Starting the software:

The smaXtec Sensors must first be initialised so that you can use the functions in the smaXtec System. Please ensure you have the provided magnet and the earmark number/identification/name of the cow that you want to fit with a smaXtec pH & Temp Sensor or smaXtec Temp Sensor to hand. Please first connect the **smaXtec USB Antenna** to your PC/laptop and then start the **smaXtec Sensor Init** software.

3. Follow the initial steps as prompted by the software

Please read the instructions in the Quickstart Guide provided. To start the initialisation process for the **smaXtec pH & Temp Sensor** as well as for the **smaXtec Temp Sensor**, switch the sensor on using the magnet provided and follow the steps as prompted by the software. For further information, please go to <http://support.smaxtec-animalcare.com>.



Following application the smaXtec pH & Temp Sensor/smaXtec Temp Sensor remains in its final position in the reticulum and can no longer be accessed. It is therefore essential that you act in strict compliance with all instructions provided by smaXtec with regard to activation, calibration and application of the smaXtec pH & Temp/smaXtec Temp Sensor (see Chapter 4.1.2).

The smaXtec USB antenna can now be disconnected from your USB interface. It is only required for initialisation of the smaXtec pH & Temp Sensor/smaXtec Temp Sensor.

6.2. smaXtec Messenger

Actionable alerts with smaXtec Messenger. Monitoring of data in real-time means that you are kept fully informed at all times about the current health status of your animals. smaXtec Messenger automatically sends you alerts detailing any abnormal changes in rumen pH and/or temperature. The suggested actions contained in the alert messages enable a rapid and targeted response – for example by adjusting feed quantities or contacting your veterinarian.

The suggested actions and physiological data can be accessed in real-time via PC, tablet or smartphone from the smaXtec server. Where necessary you automatically receive an actionable alert message. You are kept easily updated with all relevant information. This means curves no longer have to be looked at or interpreted.



Caution: For display of temperature events, measurement data must be available from over a minimum of 4 days.

6.2.1. Initial Steps

Start your browser, go to <https://messenger.smaxtec.com> and log in with your user name and password. Before using for the first time you need to register and set up an account for your farm. Please use a valid email address.

The most important steps for successful use of **smaXtec Messenger** are described below:

1. Setting up a smaXtec Messenger user account:

Start your browser, go to <https://messenger.smaxtec.com> and register with your chosen user name and password. Please use a valid email address at registration.

2. Activating a device

Activation in **smaXtec Messenger** – menu item “Add device”: Device activation is required to ensure the device can be correctly assigned to the user and used in smaXtec Messenger.

3. Working with smaXtec Messenger

Following successful activation of devices (smaXtec Base Station, Repeater, Sensors) you can administer your cows, view important events and manage alerts.

7. Technical Details

smaXtec pH & Temp Sensor:

- Dimensions: 132 x 35 mm (length x diameter)
- Internal antenna
- Internal Li-Ion battery
- Measurement interval: 10 min
- Measurement range (pH): pH 3 - 9
- Measurement range (temperature): 0°C - 80°C/32 - 176°F
- Measurement accuracy (pH): up to day 90 pH ± 0.2 up to day 150 pH ± 0.4
- Measurement accuracy (temperature at 39°C/102.2°F): $\pm 0.25^\circ\text{C}/\pm 0.45^\circ\text{F}$
- Measurement resolution (pH): pH 0.01
- Measurement resolution (temperature): 0.01°C/0.018°F
- Measurement duration pH measurement: 150 days (5 months)
- Battery life temperature measurement: up to 4 years
- Shelf life from date of manufacture: 18 months
- Guarantee: 50 days
- Internal measurement value memory: 50 days
- Ruminal fluid resistant – materials tested by the German Agricultural Society (Deutsche Landwirtschaftliche Gesellschaft DLG)

smaXtec Temp Sensor:

- Dimensions: 105 x 35 mm (length x diameter)
- Internal antenna
- Internal Li-Ion battery
- Measurement interval: 10 min
- Measurement range: 0 - 50°C/32 - 122°F
- Measurement accuracy (at 39°C/102.2°F): $\pm 0.25^\circ\text{C}/\pm 0.45^\circ\text{F}$
- Battery life: up to 4 years
- Internal measurement value memory: 60 days
- Ruminal fluid resistant - materials tested by the German Agricultural Society (Deutsche Landwirtschaftliche Gesellschaft DLG)

smaXtec Climate Sensor:

- Dimensions: 120 x 90 x 50 mm (L x B x H)
- Dustproof and waterproof plastic casing (IP 66)
- Internal antenna
- Internal Li-Ion battery

- Measurement interval: 10 min
- Measurement range (relative humidity): 0 - 100% RH
- Measurement range (temperature): -40 - +80°C/-40 - +176°F
- Measurement accuracy (temperature; 5 – 60°C/41 – 140°F): ± 0.3 °C/ ± 0.54 °F
- Measurement accuracy (relative humidity; 20 – 80% RH): $\pm 2\%$ RH
- Battery life: up to 5 years

smaXtec Base Station:

- Dimensions: 200 x 110 x 60 mm (L x B x H)
- Dustproof and waterproof plastic casing (IP 66)
- Antenna for communication with smaXtec Sensors (max. distance 5-10 m/16-32 ft)
- Antenna for communication with smaXtec Repeater (max. distance between read-out stations 100 m/328 ft)
- Device available with GSM/without GSM
- Ethernet socket: for internet connection
- Supplied mains adapter: 100-240V, 2A
- Easy to mount on a wall or wooden post etc.

smaXtec Repeater:

- Dimensions: 160 x 90 x 50 mm (L x B x H)
- Dustproof and waterproof plastic casing (IP 66)
- Antenna for communication with smaXtec Sensors (max. distance 5-10 m/16-32 ft)
- Antenna for communication with smaXtec Base Station/Repeater (max. distance between read-out stations 100 m/328 ft)
- Supplied mains adapter: 100-240V, 2A
- Easy to mount on a wall or wooden post etc.

smaXtec USB Antenna:

- Dimensions: 21 x 85 x 12.4 mm (B x L x H)
- Internal antenna
- USB Type A plug
- USB 2.0