

## **USING THE DRY4LIFE BEDWETTING ALARM**

### **SET UP:**

- 1) Insert 1 size C battery into the alarm clock.
- 2) Press and release the silver button on top of the alarm. It will ring to indicate it is functional.
- 3) Leave the button in the downright position until ready for bedtime. It is now deactivated. (The clock will continue to tell time, but the wake up alarm is permanently disabled to prevent confusion.)

### **PREPARATION FOR BED:**

- 1) Cut a slit in a common, disposable, self-adhesive panty shield.
- 2) Insert the wet-sensor into the panty shield, the metal side facing the body.
- 3) Peel and stick the panty shield containing the wet sensor to the inside of close fitting underwear. (As close to the source of urine as possible.)
- 4) Release the alarm button. It will ring indicating the unit is ready to receive signals from the transmitter when the user wets.

### **WHEN THE ALARM GOES OFF:**

- 1) Depress the button. The alarm will stop sounding.
- 2) Clean the wet-sensor with soap and water.