



Empowering Caregivers

Arial® Wi-Fi Wireless Nurse Call System

Deployment & Maintenance Guide for UL1069

0981-739-000 REV B

KB #: 12472

Published: June 2022



Disclaimer

The information and know-how included in this document are the exclusive property of STANLEY Healthcare and are intended for the use of the addressee or the user alone. The addressees shall not forward to another their right of using the information, know-how or document forwarded herewith, in whole or in part in all matters relating or stemming from or involved therein, whether for consideration or without consideration, and shall not permit any third party to utilize the information, know-how or the documents forwarded herewith or copies or duplicates thereof, unless at the company's consent in advance and in writing. Any distribution, advertisement, copying or duplication in any form whatsoever is absolutely prohibited. The Company reserves the right to sue the addressee, user and/or any one on their behaves, as well as third parties, in respect to breaching its rights pertaining to the intellectual rights in particular and its rights of whatever kind or type in the information, know-how or the documents forwarded by them herewith in general, whether by act or by omission.

This document is confidential and proprietary to STANLEY Healthcare and is not to be distributed to any persons other than licensed Arial System users or other persons appointed in writing by STANLEY Healthcare.

Trademark Acknowledgements

AeroScout is a trademark of Securitas or its affiliates. Other brand products and service names are trademarks or registered trademarks of their respective holders. Below is a partial listing of other trademarks or registered trademarks referenced herein:

Cisco™ is a trademark of Cisco Systems, Inc.

Aruba Networks® is registered trademark of Hewlett Packard Enterprise

Sun, Sun Microsystems, the Sun Logo, Java, JRE and all other Sun trademarks, logos, product names, service names, program names and slogans that are referred to or displayed in this document are trademarks or registered trademarks of Sun Microsystems, Inc. in the United States and other countries. This product includes software developed by the Apache Software Foundation (<http://www.apache.org/>).

This product includes code licensed from RSA Data Security

Skype, SkypeIn, SkypeOut, Skype Me, the Skype Logo and the S logo and other marks indicated on Skype's website are trademarks of Skype Limited or other related companies.

Esper is a trademark of EsperTech, Inc.

Jboss is a trademark of Red Hat Middleware, LLC.

Oracle 10G is a registered trademark of Oracle Corporation and/or its affiliates.

Windows and MS SQL Server are registered trademarks of Microsoft Corporation in the United States and/or other countries.

JasperSoft, the JasperSoft Logo, JasperReports, the JasperReports logo, JasperIntelligence, JasperDecisions, JasperAnalysis, Scope Center, Scope Designer, and JasperServer are trademarks or registered trademarks of JasperSoft, Inc. in the United States and other countries.

NETGEAR is a registered trademark of NETGEAR, Inc.

CONTROL and Device Manager are registered trademarks of Peperl_Fuchs Control, Inc.

Table of Contents

Introduction	7
Supervised Wireless and Electrical Data Circuits	7
Arial Wi-Fi Power Connection Overview	8
Uninterruptable Power Supply	9
Arial Wi-Fi Supervised Data Connection Diagram	10
System Software Requirements	10
Wireless Infrastructure Compatibility	10
Reference Documents & Software	11
Solution Component Overview	12
Fundamental Arial Wi-Fi Nurse Call System Components	14
Server Software	14
Ethernet Switches	15
Gateway GW3X00	16
GW3X00 Accessory Power Supply	20
Gateway GW1000	21
System Sensor Siren (Chime)	24
Wi-Fi Call Stations (CSK-1000 and CSK200 Series)	27
How to Use a Call Station (CSK-1000 Series)	29
How to Use a Call Station (CSK200 Series)	40
Remote Push Button Call Cords	46
Dome Lights	48
Pendants	50
How to Use the P10 Pendant	50
Supplemental Hardware Components	52
Display Sign	53
Paging Base Station	54
Arial Pager	56
Adam IP I/O Module	57
Arial Software Components	58
Arial Software	58

Arial Mobile Application (Supplemental)	59
Solution Deployment Checklist.....	60
Assigning IP Address for PC to Communicate on Fundamental Network	61
Deploying the AeroScout Location Engine	62
Important ALE Solution Deployment Notes.....	62
ALE Deployment Checklist	62
Configuring the ALE Collection Time.....	63
SNMP Trap Interval.....	64
Location Engine Map Configurations	65
Steps for Adding Maps in AES.....	66
Arial and ALE Integration.....	72
Applying the UL1069 License in Arial.....	72
Configuring Arial for ALE.....	75
Adding and Configuring Gateways in Engine.....	78
Adding the GW1000 Gateway	78
Adding the GW3X00 Gateways.....	82
Installing Gateways	86
GW3X00 Wall Mounting Parts	86
Mounting the Gateway on a Wall	86
Installing and Configuring Sirens	90
Selecting 56108 Siren Audio	90
Mounting a 56108 Siren	91
56108 Gang Box Mounting	92
56108 Surface Mounting.....	93
Siren 56108B Mounting Instructions	95
Configuring Gateways to Activate Sirens.....	97
Prerequisites	97
Configuring HPS Server Supervision.....	97
Deploying Arial Wi-Fi Call Stations.....	99
Activating the Call Station.....	99
Call Station Alerts	99
Dismissing Call Station Alerts.....	99
Batteries and External Power	99
Preparing the Batteries	99
External Power	100

Using External Power	101
Turning the Call Station On or Off	101
Installing and Removing the Call Station	101
Mounting the Bracket	102
Attaching the Call Station	104
Adjusting the Pull Cord Length.....	105
Shortening the Pull Cord	105
Using the Strain Relief Spools	106
Removing the Call Station	108
Adding Call Stations to Arial	109
Task 1: Add a Call Station to Arial	110
Task 2: Add the Call Station to a Map	113
Configuring Gateways to Activate Dome Lights.....	116
Prerequisite	116
GW1000 to Dome Light 4 LEDs	120
ADAM 6060 to Dome Light with 4 LEDs.....	121
GW1000 to Dome Light with 2 LEDs	122
ADAM 6060 to 2 LED Dome Light	123
Setting the Default Alarm Sound	124
Call Station Escalation Alerts	124
Check-In Events (optional).....	125
Reports	126
Solution Testing.....	127
General Solution Testing Requirements	127
Call Station Testing	127
Testing Requirements	127
Testing Call Station Alert/Dismiss Events.....	127
Backing Up and Restoring.....	129
Arial	129
Location Engine	129
Service and Maintenance	130
Servicing	130
Fault Indicators Requiring Service	130
Routine Maintenance	131
Power Back-Up.....	131

Uninterruptable Power Supply Battery Replacement.....	131
Verifying Data Backups	134
Replacing a Removable Drive	135
Cleaning the PC Fan	137
Arial Software Updates.....	140
Computer Operating System Updates	140
Operation Testing	143
Coverage Testing	143
Appendix	145
Wiring Diagrams	145
GW1000 to 56108 Siren (Chime) – Server	145
GW1000 to 56108B Siren (Chime) – Server	146
2-LED Dome Lights to GW3X00	147
4-LED Dome Lights to GW3X00	147
Troubleshooting	148
ALE SNMP Communication Failure	148
Safety, FCC Warnings and Warranty	149

Introduction

The STANLEY Healthcare Arial Wi-Fi Wireless Nurse Call System, for Skilled Nursing Communities, is designed to enhance the safety of senior living residents, while also providing advanced tools for event management, staff notification and reporting.

The Solution uses fixed Call Stations, in resident rooms and bathrooms, together with Wi-Fi wireless technology, to help provide residents throughout the senior living community a convenient way to summon help.

When a resident activates a call for assistance, an alert signal is sent over the wireless network to the Arial Server and is displayed at a nurse's station PC with the alarm point's details and room information. Alarms are reset physically from the device that generated the alert.

Reports

Reports are available to show alarm activity, caregiver response times, as well as other maintenance and system information.

This guide covers the deployment steps, best practices, and maintenance procedures for the Arial Wi-Fi Nurse Call System.

Supervised Wireless and Electrical Data Circuits




Figure 1. Supervised Wireless and Electrical Data Circuits

Arial Wi-Fi Power Connection Overview




Figure 2. Arial Wi-Fi Power Connection Overview

Devices that connect to Battery Backup side of an Uninterruptable Power Supply (UPS):

- PC Running Arial Server and ALE Software
- Main Network Switch

IMPORTANT: The UPS must be connected to a generator backed up outlet. The battery backup will all keep the server PC and these critical pieces of equipment running for at least five minutes while the generator starts and takes load, and once again when power from the grid is restored and the generator backed-up outlet transfers back to grid power.

NOTE: The GW1000 and Siren are backed up as well since they are powered from the POE port of the Main Network Switch, which is backed up by the UPS.




Figure 3. Backup and Surge outlet locations on the back of the UPS

Uninterruptable Power Supply

STANLEY Healthcare supplies an Uninterruptable Power Supply (sold as item number 0768-004) with systems it installs.

Manufacturer & Model	APC BK500
Output	10A Max 120VAC 60Hz +/- 3Hz
Input	98-140VAC 60Hz +/- 3Hz
Backup Time	95.5 Min @ 30W, 3.4 min @ 300W (Max Load)
Battery	Lead-acid RBC2, 106VA, 4-6 Yr. Life
Physical	3.58" x 11.18" x 6.5" (28.4 cm x 9.1 cm x 16.5cm) 13.82 lbs (6.27Kg)
Environmental	32-104 °F (0-40 °C), 0-95% RH
Compliance	UL1778, Energy Star V2.0, FCC Part 15, RoHS

Arial Wi-Fi Supervised Data Connection Diagram




Figure 4. Arial Wi-Fi Supervised Data Connection Diagram

System Software Requirements

This section describes required software requirements for the Arial® Wireless Nurse Call System.

Software	Minimum Requirements & Notes
Arial	10.6.2 and above
AeroScout Location Engine (ALE)	5.5 (5.5.5.168) and above <i>Refer to the AeroScout Location Engine 5.5 Deployment Guide for full ALE system requirements.</i>

Wireless Infrastructure Compatibility

The Arial Wireless Nurse Call system uses a Wi-Fi wireless infrastructure that uses a network of gateways to receive transmissions and relay them to the Arial® software. The Arial Wireless Nurse Call system's software and fundamental components must be installed on a separate and dedicated network that connects to the Gateways. Other network connections can be made to the Arial server, provided it is not part of the UL1069 network.

Reference Documents & Software

The following articles can be accessed by logging into the STANLEY Healthcare Support Community site at the following

URL: <https://stanleyhealthcare.force.com/HomePage>.

KB #	Documents
12472	Arial Wi-Fi Nurse Call System Deployment and Maintenance Guide for UL1069
12138	Arial 10.6 Quick Reference Guide
12137	Arial 10. 6 Wireless Emergency Call Solution Admin Guide
12139	Arial 10.6 Software Installation Guide
11837	Gateway GW3000 and GW3100 Data Sheet
11715	Gateway GW3000/GW3100 Installation & Configuration Guide
10334	Arial 56108 Siren Datasheet
12473	Arial 56108B Siren Datasheet
12271	AeroScout Location Engine 5.5 Deployment Guide
1100	Product Versions Compatibility Matrix Reference Guide
7870	Arial Paging Transmitter Installation Guide
12141	Arial 10.6 Mobile Application User Guide
10336	Arial 2-LED Datasheet
10337	Arial 4-LED Datasheet

KB #	Software
12322	Arial 10.6 Software Download
12274	Location Engine 5.5 for Non-Cisco Environment Download

Solution Component Overview

This Chapter describes the fundamental and supplemental hardware and software components used for the Arial Wi-Fi Nurse Call System for UL1069 Deployments.

Fundamental Components

The fundamental components are the required components for the Arial Nurse Call System.

NOTE: It is required that one of the CSK-1000, CSK-1000-CD, or CSK-1000MR call stations be assigned to each resident apartment and be permanently located in the bathroom of that apartment. At least one 1069-54372 call notification station serving the property must be permanently placed and continually operating at a fixed location. Dome lights shall be located in the corridor directly outside of each resident apartment.

Supplemental Components

The supplemental components may already exist in a facility from a previous Arial deployment or can be optionally added. They are not required for the Arial Wi-Fi Nurse Call System.

NOTE: The nurse call system must not share its nurse call communication network with any other network or system. All transmitters using the same frequency or channel and that are within range of a receiver of a wireless system, must be compatible with the nurse call system and evaluated for the intended purpose. Supplementary devices not evaluated for the intended purpose are permitted only if their connection to the network utilizes a compatible device evaluated for the purpose.




Figure 5. Arial Wi-Fi Nurse Call System Overview

Fundamental Arial Wi-Fi Nurse Call System Components

This section describes the following fundamental software and hardware components.

- Server Software
- Ethernet Switches
- Gateway GW3000 & GW3100
- Gateway GW1000
- System Sensor Siren (Chime)
- Wi-Fi Call Stations
- Remote Push Button Cords
- Pendants

Server Software

Model	Comments
SKU: 54302 Arial Software	Latest version of Arial Server Software
SKU: 1069-54304 Arial UL1069 License	License that grants customer the right to use one instance of the Arial Server Software in their community. Application of license automatically adjusts key settings for compatibility with UL1069 devices.
SKU: ENG4100 Arial Wi-Fi and Location Engine Integration License	License that grants customer the right to use one instance of the AeroScout Location Engine software in their community. License is required to support GW1000 and Trouble Siren that are required for UL1069

The Arial Server software runs on a Windows PC. The Arial Server receives alarm and status signals from call stations, processes this information, and relays it to Client computers and messaging devices such as IP and mobile phones, pagers and display signs. The Server also stores the database of system activity. Note that the Server can also act as a Client PC, displaying the Arial software for alarm annunciation and other day-to-day tasks.

Arial 10.6 supports the integration of the Arial software with the AeroScout Location Engine (ALE).

Ethernet Switches

Models	Comments
SKU: SWC-POEP8-1000 8 Port Switch	Unmanaged 8-Port Ethernet Switch with PoE Outputs UL 60950 Listed
SKU: SWC-POEP16-1000 16 Port Switch	Unmanaged 16-Port Ethernet Switch with PoE Outputs UL 60950 Listed

The Ethernet Switch is the ‘central hub’ for the solution’s components and must conform to UL60950 or UL62368 standards.




Figure 6. 8-Port Switch Front View




Figure 7. 8-Port Switch Back View




Figure 8. 16-Port Front View




Figure 9. 16-Port Back View

Gateway GW3X00

Models	Comments / Firmware Version
SKU: GW3000 and GW3100	Wi-Fi receiver and dome light controller Firmware Version: 415.15 or higher

Gateways GW3100 and GW3000 serve as receivers for messages from Wi-Fi and BLE Tags. Received Tag messages are relayed to the AeroScout Location Engine (ALE).

The Gateways also include Wi-Fi and BLE transceivers, which monitor the Wi-Fi and BLE receivers. The Gateways report an error if an issue with the receivers is detected.

The LED indicator of the GW3000/GW3100 Gateways displays color based on the device status as follows:

Status	LED
Connected to the network via LAN	Green
Connected to the network via Wi-Fi	Ivory
Not Connected to the Network	Blinking Green
Malfunction	Red
Blink	Blink 5 times in Ivory for a total duration of 3 seconds
Firmware Update	Blink in Ivory (similar to phase 4)
IP Reset	Upon Reset, single LED blink



LED GW3000/GW3100 Gateway

Figure 10. Gateway GW3000/GW3100




Figure 11. GW3000 & GW3100 Connections and Controls – Viewed from Back

Specifications: GW3X00

Physical and Mechanical	<p>Dimensions: 5.5 Length x 5.5 Width x 1.19 Thickness inches (140 x 140 x 30 mm)</p> <p>Weight: Including mounting plate:</p> <p>GW3000 243.0 g (8.57 ounces)</p> <p>GW3100 245.5 g (8.65 ounces)</p> <p>Housing: Polycarbonate and ABS</p>
Power	<p>Input Voltage: 24-48 VDC</p> <p>Power Consumption: 10.2 Watts @ 24 VDC, 30 Watts @ 48VDC</p> <p>PoE (802.3af)</p> <p>PoE+ (802.3at)</p>
Network Interface	Ethernet 802.3 100Base-T (RJ-45) Wi-Fi 802.11 b/g/n
Output Ratings	<p>Up to 30mA @ 24VDC per output using 24VDC output of the Gateway. 120mA total @ 24VDC if all four outputs used at maximum can be supplied by the GW3X00.</p> <p>If an external power supply is used to power remote devices, each output can switch up to .5A @ 30VDC.</p>
Environmental	<p>Operating temperature: 0°C to 49°C (32°F to 120°F)</p> <p>Humidity: 0 to 93%, non-condensing</p>
Certification and Regulatory Compliance	<p>UL2560, UL1069</p> <p>Safety: USA/Canada: UL 60950-1 2nd edition; CAN/CSA-C22.2 No. 60950-22:17; UL 62368-1 2nd edition; IEC 62368-1:14</p> <p>Europe: EN 60950-1:06+A11:09+A1:10+A12:11+A2:13; EN 62368-1:14+A11:17</p> <p>International: IEC 60950-1:05+A1:09+A2:13; IEC 62368-1:14</p>

	<p>EMC: EU/USA/Canada: Draft EN 301 489-1 V2.2.0; Draft EN 301 489-17 V3.2.0; 47 CFR Part 15 Sub B; ICES-003</p> <p>Radio: Europe: EN 300 328 V2.1.1:16 Canada: RSS-247 USA: FCC-certified</p>
--	--

GW3X00 Accessory Power Supply

STANLEY Healthcare sells a separate power supply (Item WGB-ADP-047-U) that can be used with Gateways GW3100 and GW3000 if they are not powered from a PoE network switch.

Manufacturer & Model	Adapter Tech. ATS050T-P480
Input Power	100-240VAC 50/60Hz 1.2A Max
Output Power	48VDC 1.05A
Environmental	32- 104°F (0-40°C), 20-80% RH
Size and Weight	115mm x 53mm x 38mm, 310g
Compliance	UL2468, UL94V-0, FCC Part 15 Class B, RoHS

Gateway GW1000

Models	Comments / Firmware Version
SKU: GW-1000 (GW-1000)	Wi-Fi receiver and dome light and siren controller Firmware Version: 415.15 or higher

The GW1000 monitors the software on the PC, acting as a watchdog. If the PC or software stops running, the GW1000 activates an output to turn on a siren that is controlled by the output. This siren serves as a trouble signal and will activate within 90 seconds should the Engine software or PC stop executing.



Figure 12. Gateway GW1000




Figure 13. GW1000 Connections and Controls – Viewed from Back




Figure 14. GW1000 I/O Connectors

IO Connector Descriptions

Letter	Description
a	24-48 IN: Power connection accepts 24-48 V Direct Current
b	GND IN: Ground/Earth
c	24V OUT 1: 24V output
d	R1 NC: Relay 1 normally closed connection
e	R1 COM: Relay 1 common connection
f	R1 NO: Relay 1 normally open connection
g	R2 NC: Relay 2 normally closed connection
h	OUT1: Open drain output 1
i	OUT 2: Open drain output 2
j	OUT 3: Open drain output 3
k	OUT 4: Open drain output 4
l	OUT 5/IN 1: Open drain output 5 / Input 1
m	OUT 6/IN 2: Open drain output 6 / Input 2
n	OUT 7: Open drain output 7
o	OUT 8: Open drain output 8
p	R2 COM: Relay 2 common connection
q	R2 NO: Relay 2 normally open connection
r	GND: Ground Input connection/Earth
s	OUT 9/IN 3: Open drain output 9 / Input 3
t	OUT 10: Open drain output 10
u	24V OUT 2: 24V output
v	OUT 13/IN 4: Open drain output 13 / Input 4
w	GND: Ground Input connection/ Earth
x	OUT 11/IN 5: Open drain output 11 / Input 5
y	OUT 12/IN 6: Open drain output 12 / Input 6

Specifications: GW1000

Physical and Mechanical	<p>Dimensions: 245mm X 200mm X 60mm (9.6in x 7.9in x 2.4in)</p> <p>Weight: 865g (31oz)</p> <p>Housing: Polycarbonate and ABS</p>
Network Interface	Ethernet (RJ-45) Wi-Fi 802.11b/g/n
Firmware	DSP: 415415 or higher, BOOT: 61515 or higher
Power	<p>Input voltage: 24-48VDC PoE (802.3af) 48VDC</p> <p>Maximum power consumption: 10W.</p>
Environmental	<p>Operating temperature: 0°C to 49°C (32°F to 120°F)</p> <p>Humidity: 0 to 93%, non-condensing</p>
Relays	<p>Two: 220 VDC/ 250 ACD, 1 Amp, NO (Normally Open) or NC (Normally Closed)</p> <p>24VDC Output, up to 300mA</p> <p>13 Open-Drain Outputs (up to 100mA each)</p>
Certification and Regulatory Compliance	UL2560, UL1069

System Sensor Siren (Chime)

Model	Comments
<p>SKUs: 56108 (System Sensor CHWL) and 56108B (System Sensor MHW)</p>	<p>Chime with selectable chime tones and volume settings. 24V, UL approved.</p>

Sirens (Chimes) produce a distinctive chime tone for alerting trained personnel to investigate Server Software or Hardware failures.

Specifications: 56108 and 56108B Sensor Sirens

Physical and Mechanical	<p>56108 Dimensions: 5.6 in L × 4.7 in W × 1.25 in D (143 mm L × 119 mm W × 32 mm D)</p> <p>Wall Surface Mount Black Box (SBBWL):</p> <p>5.6 in L × 4.7 in W × 4.3 in D (142 mm L × 119 mm W × 109 mm D)</p> <p>56108B Dimensions: 4.6 in L x 2.9 in W x .45 in D (117 mm L x 74 mm W x 11 mm D)</p>
Environmental	<p>Standard Operating Temperature: 32°F to 120°F (0°C to 49°C)</p> <p>Humidity: 10 to 93% Non-condensing</p>
Voltage	Nominal Voltage: Regulated 12VDC or regulated 24DC/FWR
Input terminal wire gauge	12 to 18 AWG
Certification and Regulatory Compliance	UL2560, UL1069, UL464




Figure 15. 56108 Siren




Figure 16. 56108 Siren Wiring Illustration

The siren only requires two wires for power and supervision.




Figure 16B 56108B Siren shown (left) with wire connections and sound select switch shown (right).

Use Position 4 (normal) or 3 (loud) for sound select switch. Positions 1 and 2 shall not be used.

Wi-Fi Call Stations (CSK-1000 and CSK200 Series)

It is required that one of the Call Stations assigned to each resident apartment should be permanently located in the bathroom of that apartment.

CSK-1000 Series

Models	Comments
SKU: CSK-1000	Standard Wi-Fi Call Station
SKU: CSK-1000-MR	Moisture Resistant Call Station
SKU: CSK-1000-CD	Code Call Station


The CSK-1000, CSK-1000-MR (Moisture Resistant) & CSK-1000-CD (Code Call) Wi-Fi Call Stations are components of STANLEY Healthcare's Wireless Call Solutions. The wireless, battery-operated fixed-location devices are designed to trigger alarms when the call button is pressed, the pull cord is pulled, or when an optional remote push-button (EAC-1000) is used.

The Wi-Fi Call Stations are equipped with a large call button, a durable antibacterial pull cord and a dual-purpose Dismiss/Check-In button for dismissing alarms and daily resident check-ins. Triggered alarm information is shown in Arial, and includes the device's fixed location, resident details and event times. Additionally, response time and resident check-in reports can be generated by the Arial software.

CSK-1000 Series Key features:

- Wi-Fi Compatible
- Large and Tactile Call Button
- Durable Antibacterial Pull Cord
- Battery and Optional AC/DC Power
- Back-Up Power
- Visual and Audible indications
- Resident Check-Ins
- Night Illumination
- Fully Supervised by Arial
- Easy Mounting
- Moisture Resistant (CS100MR)

Front View



Back View





Figure 17. Call Station Features and Controls

Battery Compartment



Bottom View




Figure 18. Call Station Features and Controls



Figure 19. Code Call Station

How to Use a Call Station (CSK-1000 Series)


The Call Station turns on automatically and is ready for use once power is supplied (battery or power adaptor) and a confirmation tone is heard. To turn the Call Station Off, simply remove the batteries or un-plug the power adaptor.

Activating a Call Station Alarm

The Call Station alarm is activated by performing either of the following:

- Pressing the large call button (1)
- Pulling the pull cord (2)
- Optional: Pushing the remote push-button

A confirmation tone is heard when an alarm is activated. An alarm message is sent to Arial. The Call Station remains in alarm mode (flashing red LED) until the alarm is physically dismissed.



Dismissing a Call Station Alarm

IMPORTANT! Before dismissing an alarm, a staff member must attend to the resident and dismiss the alarm on the Call Station.

Press on the Dismiss/Check-in button (1) to dismiss the alarm. A confirmation tone is heard and LED stops flashing red. The alarm is automatically dismissed in Arial.



Volume Settings

The Wi-Fi Call Station has eight (8) volume levels ranging from silent (mute) to high. The volume buttons are located inside the battery compartment. Voice commands ('minimum volume', 'maximum volume') are heard when minimum or maximum volume is reached.

1. Unscrew (1) and open (2) the battery compartment.



2. Press the volume down (3) or volume up (4) buttons repeatedly to select the desired volume. A tone is heard as you adjust the volume. A voice command is heard once minimum or maximum volume is reached.
3. You can mute the volume by pressing the volume down (3) button repeatedly until the 'Minimum Volume' voice command is heard.




CSK-1000 Series LED and Audio Indications

Action	LED Indication	Audio Indication
Power On	Green	Confirmation tone
Alarm Call	Flashing Red	Confirmation tone
Alarm Reset	Green	Confirmation tone
Check-In	Green	Confirmation tone
Low Battery	N/A	Single negative tone. <i>Low battery levels are displayed in Arial.</i>
Volume Up/Down	N/A	Tone on presses and voice feedback for minimum or maximum volume

Night Illumination

The Wi-Fi Call Stations have two (2) green dim night illumination LEDs on the front panel for use in darkness. The LEDs are On by default, and provide enough light for comfort and don't illuminate the room. The LEDs flash dimly every 3 seconds if the Call Station is powered by batteries. When the device is powered by the power adaptor, the LEDs are On constantly.



Daily Check-in


A check-in is performed by pressing the Dismiss/Check-In button (1) once. A confirmation tone is heard and a check-in alert is sent to Arial.



Turning Night Illumination Off/ON

The night illumination LEDs can be turned off at any time, if needed. The LEDs are turned Off by pressing the volume up and volume down buttons simultaneously and holding for 3 seconds. The same procedure is used to turn the LEDs back On.

1. Unscrew (1) and open (2) the battery compartment.
2. Press the volume down (3) and volume up (4) buttons simultaneously and hold for 3 seconds. The night LEDs are turned Off/On.



Specifications (CSK-1000, CSK-1000-MR, CSK-1000-CD)

Physical and Mechanical	<p>Dimensions: 415 mm x 84 mm x 40 mm (5.70 in x 3.30 in x 1.57 in)</p> <p>Weight: 200 g (7.05 oz)</p>
Radio	<p>802.11 radio (2.4 GHz); b/g/n compliant Low frequency receiver for chokepoint detection (125kHz) Transmission power: up to +19dBm (~81mW) Patented clear channel sensing avoids interference with wireless networks</p>
Security Protocol	<p>WPA2 security with AES-256 encryption</p>
Environmental	<p>Operating temperature: 0°C to 49°C (32°F to 120°F)</p> <p>Humidity: 0 to 93%, non-condensing (CSK-1000 & CSK1000-CD)</p> <p>Humidity: Condensing (CSK-1000-MR)</p> <p>Ingress Protection Rating: IP 20</p>
Electrical	<p>Battery: 3x 3.0V CR123A Lithium batteries (replaceable) AC/DC Power input: 5V 2A</p>
Audio	<p>Volume Level: Max 96dBA</p>
Certifications	<p>Radio:</p> <p>FCC Part 15, sub-part C class B, sub-part B EN 300-328, EN 300- 330, EN 301-489, RSS 210 (Canada)</p> <p>Safety:</p> <p>US: UL 60950, IEC60601-1 and applicable, IEC60601-1-X standards.</p> <p>Europe: CE mark - EN 60950, IEC/EN 60601-1 and applicable IEC 60601-1-X standards</p> <p>UL1069, UL2560</p>

CSK200 Series Call Stations

Models	Comments
CSK200-1069	Standard Wi-Fi Call Station for UL1069 Systems
CSK200-1069MR	Moisture Resistant Wi-Fi Call Station for UL1069 Systems

It is required that one of the Call Stations assigned to each resident apartment should be permanently located in the bathroom of that apartment.


The CSK200 and CSK200-MR (Moisture Resistant) Wi-Fi Call Stations are components of STANLEY Healthcare's Wireless Call Solutions. The wireless, battery-operated fixed-location devices are designed to trigger alarms when the call button is pressed, the pull cord is pulled, or when an optional remote push-button (56115 and 56116) is used.

The Call Stations are equipped with a large call button, a durable antibacterial pull cord and a dual-purpose Dismiss/Check-In button for dismissing alarms and daily resident check-ins. Triggered alarm information is shown in Arial and includes the device's fixed location and event times. Additionally, response time and resident check-in reports can be generated by the Arial software.


CSK200 Series Key features:

- Wi-Fi Compatible
- Large and Tactile Call Button
- Durable Antibacterial Pull Cord
- Battery Power
- Visual Status Indication
- Resident Check-Ins
- Fully Supervised by Arial
- Easy Mounting
- Moisture Resistant Option


Front View



Back View



Bottom View




Antimicrobial Pull Cord


Case Retaining Screw

1/4" Accessory
Pushbutton Cord Jack

Battery Location (Inovonics)



Battery Location (Wi-Fi)



How to Use a Call Station (CSK200 Series)

The Call Station turns on automatically and is ready for use once power is supplied by inserting the battery. To turn the Call Station Off, simply remove the battery.

Activating a Call Station Alarm

The Call Station alarm is activated by performing any of the following:

- Pressing the large call button (1)
- Pulling the pull cord (2)
- Optional: Pushing the remote push-button

The alarm is confirmed by a flashing red LED indicator (3). An alarm message is sent to Arial. The Call Station remains in alarm mode (flashing red LED) until the alarm is physically dismissed by pressing the reset button (4).



Dismissing a Call Station Alarm

IMPORTANT! Before dismissing an alarm, a staff member must attend to the resident and dismiss the alarm on the Call Station.

Press on the Dismiss/Check-in button (4) to dismiss the alarm. The LED stops flashing red. The alarm is automatically dismissed in Arial.


IMPORTANT! If an optional remote push button cord is used, verify the cords are securely connected. Failure to correct a loosened connection may prevent the remote cord from placing an alarm. A cord out notice will be active in Arial if the cord is disconnected.

CSK200 Series LED Indications

Action	LED Indication
Alarm Active	Red LED – “Fast” 3Hz Flash
Standby (No Active Alarm)	LED Off
Staff Presence	Red LED – “Slow” .5Hz Flash Note: LED Stops after 10 minutes
Resident Check-in	Red LED – Single Flash

Daily Check-in

A check-in is performed by pressing and releasing the Dismiss/Check-In button (1) once. The Red LED (2) will flash once to provide feedback that a check-in event was sent to Arial.



Staff Presence

Staff presence can optionally be used by facilities to indicate when staff enter a resident room to respond to an active alarm. A caregiver presses and holds the Dismiss/Check-In button (1) for at least two seconds on the Call Station that the resident used to place an alarm. The LED (3) changes from the fast-flashing rate of an alarm to a slower flashing rate indicating a staff presence. A staff presence message is

sent to Arial. Once the caregiver is finished assisting the resident, they press and release the Dismiss/Check-In button to end the alarm and staff presence events.

Note: After 10 minutes, the Call Station LED indicator stops flashing for staff presence to extend the battery life, but the alarm and staff presence events remain active in Arial until the Call Station is reset using the Dismiss/Check-In button.

CSK200 Series Specification Table

Models	CSK200 / CSK200-1069 Series Specifications (Wi-Fi)
Physical and Mechanical	Dimensions: 132mm x 85mm x 33mm (5.2in x 3.3in x 1.3in)
	Weight: 145 g (5.1 oz)
Radio	802.11a/b/g/n 2.4 GHz Transmission power: up to +19dBm (~81mW)
Visual Indicator	Red LED Indicates alarm (3Hz), Staff Presence (.5HZ) and 2 second blink for Check-In
Features Activation	<p>Primary Alarm Activation: Push button, pull cord, or remote button</p> <p>Staff Presence: >2 second press of check/reset pushbutton while in Primary Alarm</p> <p>Resident Check-in Activation: Press check/reset pushbutton on call station</p> <p>Primary Alarm Cancel: Press check/reset pushbutton on call station</p> <p>Secondary Alarm: Closed condition across secondary alarm input terminals (Not for use on UL call points)</p> <p>Secondary Alarm Cancel: Open condition across secondary alarm input terminals</p>
Field Connections	<p>One 1/4 inch mono jack for pushbutton or accessory cord Primary Alarm input</p> <p>Two position terminal block for Secondary Alarm input</p>
Environmental	Operating temperature: 0°C to 49°C (32°F to 120°F)
	Humidity: 0 to 93%, non-condensing (CSK200 & CSK200-1069)
	Humidity: Condensing (CSK200MR and CSK200-1069MR)
	IP 20
Electrical	Battery: 1x 3.0V CR123A Lithium batteries (replaceable)
Certifications	
	FCC Part 15 sub part C, EN 401 489-1, EN 301 489-17, EN 300 328, EN 300 220*, RSS 247
	UL 1069* (CSK200-1069 and CSK-1069MR)
	UL 2560* (CSK200, CSK200MR)
	UL/IEC 62368

*Depends on device model/SKU.

Accessories

Models	Comments
CSK200-PLATE	Plate to cover paint/wallpaper lines exposed by removing old call stations when displacing older system or call points.
CSK200-AH	Cover that fits around call station for use in rooms where ligature resistance is required.

CSK200-PLATE Specifications

Characteristic	Specification
Dimension	137mm x 90mm x 3mm (5.4in x 3.5in x .1in)
Weight	8.5g (.3 oz)



Figure 20. CSK200-PLATE

CSK200-AH Specifications

Characteristic	Specification
Dimension	145mm x 90mm x 22mm (5.7in x 3.5in x .9in)
Weight	28 g (1 oz)



Figure 21. CSK200-AH Cover

Remote Push Button Call Cords

Model	Comments / Version
SKU: 56115	Single Call Station Remote Push Button Call Cord – 7' (2.1 m)
SKU: 56116	Single Call Station Remote Push Button Call Cord – 15' (4.6 m)

This optional Remote Push Button Call cord can be used for Residents to quickly call for help by pushing the call button to trigger an alarm.



Figure 22. Remote Push Button Call Cord




Figure 23. Remote Push Button Call Cord Component Locations

Specifications: Remote Push Button Call Cord

Physical and Mechanical	Length: Approximately 7 feet/2.1 m (56115) or 15 feet/4.6 m (56116) Diameter: 1.25" (3.18 cm) at widest point Button Life: 100,000 cycles Switch Type: Normally Open SP/ST
Electrical:	Max Voltage Rating: 30 VDC
Certification and Regulatory Compliance	RoHS UL2560, UL1069

Dome Lights

Models	Comments
SKU: 0900-685	2 LED Dome Light
SKU: 0900-686	4 LED Dome Light

It is required that one of the Dome Lights shall be installed in the corridor outside of the resident room and that the light shall illuminate when call station(s) in the room are in alarm.

The LED Dome Lights are used to indicate a call for assistance from a resident's room, common area, or as a zone indication.

Key features:

- Long-life LED
- Modern appearance
- Bright white color with frosted lens
- Mounts to standard single or double-gang electrical boxed

Specifications: LED Dome Lights

Dimensions:	5.0 x 4.8 x 2.755" (12.7 x 12.2 x 7.0 cm)
Electrical:	24VDC, 40 mA per indicated LED color
Weight:	10 oz. (280 g)
Environmental:	Indoor 32°–120° F (0°–49° C) Up to 90% relative humidity (non-condensing)
Certification and Regulatory Compliance	UL 1069





Figure 24. 0900-685 2 LED Dome Light (right) and 0900-686 Dome Light (left)

Pendants

Model	Comments
SKU: TAG-P1000	P10- Wi-Fi Resident Pendant Tag Basic model includes lanyard

The P10 Pendant is a component of the STANLEY Healthcare Arial Wi-Fi Nurse Call Solution targeted for the Senior Living market. The battery powered P10 Pendant is worn by residents, and enables them to move about the facility, be accurately located and, by pressing a single button, alert staff for assistance.



Key Features

- Leverages standard Wi-Fi infra-structure for low total cost of ownership and ease of deployment
- Large and tactile call button with LED and buzzer indications
- Designed for durability against significant impacts and for water and dust resistance
- Bi-directional communication for over-the-air updates
- Lanyard and belt clip attachment options
- Colored clip-on covers provide personalization for residents
- Rugged performance

How to Use the P10 Pendant

Activating the Pendant for the First Time

1. After attaching the lanyard, activate the pendant by pressing the call button for 3 seconds.
2. The LED lights green and an audible tone is played to indicate that the pendant is activated.