**OPERATION MANUAL**  
**FOR**  
**I-ASSIST KIT**  
**MODEL UIIA-771A**  
**FOR**  
**DIAGNOSTIC ULTRASOUND SYSTEM**  
**MODELS SSA-700A/SSA-770A**  
( 2 B 7 0 8 – 1 1 5 E )

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This operation manual is used at sites in Japan and overseas.  
The English pages should be used in areas outside Japan.

**TOSHIBA MEDICAL SYSTEMS CORPORATION**

**IMPORTANT!**

1. No part of this manual may be copied or reprinted, in whole or in part, without written permission.
2. The contents of this manual are subject to change without prior notice and without our legal obligation.

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# 1. Outline

## 1.1 Introduction

The IASSIST controller (hereinafter referred to as the "controller") can be used when the UIIA-771A option has been installed.

## 1.2 Outline of the Controller

IASSIST functions (selection and execution of activities) can be performed by remote control without using the operating panel. In addition, the functions of the operating panel switches such as FREEZE and GAIN can be assigned to the controller, allowing them to be remote-controlled.

This is convenient when the transducer positioning region is far from the operating panel or if you want to concentrate on the monitor during examination.

- NOTE:**
1. This controller is designed exclusively for IASSIST software that runs on the Toshiba diagnostic ultrasound system, Aplio, and cannot be used for other purposes.
  2. Standby mode cannot be used if the USIA-770A (IASSIST option) is installed. For details of the USIA-770A, refer to the operation manual for Aplio, Application Volume.

## 1.3 Definitions of Terms

- IASSIST controller : A remote controller which wirelessly controls the IASSIST functions.
- Bluetooth-USB adaptor : A receiver which is connected to the USB port of Aplio and receives radiowave control signals from the IASSIST controller.
- Charger : A battery charger for the IASSIST controller.
- Bluetooth™ : A low-power and short-range wireless standard.

### <Specifications>

- Standard version : V1.1
- Frequency band : 2.4 GHz band (2400 to 2483.5 MHz)
- Output classification : Class 2
- Maximum output : 1.0 mW (0 dBm)
- Communication distance: 20 m

## 2. Safety Precautions

- ⚠ WARNING:**
1. The controller may cause cardiac pacemakers that are in the vicinity to malfunction.  
Do not use this controller within 25 cm of persons who use a cardiac pacemaker.
  2. Do not use this controller in locations where precision instruments that may affect human life are installed, such as an operating room. Malfunction in the instruments may occur.

**⚠ CAUTION:** The transmission output power of this controller uses weak radiowaves(\*) that may affect other electronic equipment. Avoid using this controller within 25 cm of other electronic equipment.  
(\*Reference) PHS: 80 mW > Wireless LAN: 20 mW >> Bluetooth: 1.0 mW

**CAUTION:** CHANGES OR MODIFICATIONS NOT EXPRESSLY APPROVED BY THE MANUFACTURER FOR COMPLIANCE COULD VOID THE USER'S AUTHORITY TO OPERATE THE EQUIPMENT

- CAUTION:**
1. This controller is a piece of high-precision equipment. If the controller is subjected to strong shock such as falling on the floor, failure may result. Be sure to attach the provided strap and loop it over the user's neck to prevent it from falling.
  2. Do not modify or disassemble the controller. Failure may result.
  3. Do not immerse the controller in water or spill water onto the controller. Failure may result.
  4. The provided charger is dedicated to the controller. Do not use this charger for charging other equipment.
  5. Do not use the charger for purposes other than charging.
  6. Do not move the system while the controller is being charged. In addition, do not move the system while the controller is being inserted into the charger. Otherwise, the controller may come off, resulting in failure.

**NOTE:** About a use frequency band of Bluetooth

The frequency band used in this controller is also used by local area wireless stations for mobile unit identification systems that are used in industrial equipment (microwave etc.), scientific equipment, and medical equipment, and also in factory production lines and some small wireless stations (wireless LANs). Therefore, when such equipment is used in the vicinity of this controller, interference due to radiowaves (decrease in communication speed, etc.) may occur.

1. Before using this controller, confirm that there is no local area wireless station for a mobile unit identification system or a small wireless station being operated within the assumed interference distance of this controller.
2. In the event of interference to a local area wireless station for a mobile unit identification system or a small wireless station due to radiowaves from this controller, immediately change the frequency band of the local area wireless station for the mobile unit identification system or stop transmitting radiowaves.
  - \* Frequency band: 2.4 GHz
  - \* Modulation technique: FH-SS technique (frequency hopping technique)
  - \* Assumed interference distance (for a local area wireless station for a mobile unit identification system and a small wireless station): 10 m or less

**NOTE:** Battery life, charge time, and service life

A rechargeable lithium-ion (Li-ion) battery is used in this controller. The battery needs to be replaced approximately every 1 year on a for-charge basis, but this replacement period differs depending on the frequency and conditions of use. When a new battery is fully charged, the controller can be used for approximately 8 hours (depending on its usage). If the time between full charge and low charge (blinking) becomes extremely short, contact your Toshiba representative. When the controller needs to be disposed of, also contact your Toshiba representative.

(Reference value)

- Battery life : Approximately 8 hours  
This is the battery life for a fully charged new battery. The battery life varies depending on battery age and memory effect.
- Battery charge time : Approximately 2 hours and 30 minutes  
This is the charge time to fully charge a new battery after being fully discharged. The charge time varies depending on battery age and memory effect.

**NOTE:** In order to comply with FCC radio-frequency radiation exposure guidelines for an uncontrolled exposure, this device and its antenna must not be co-located or operating in conjunction with any other antenna or transmitter.

**NOTE:** The term "IC:" before the radio certification number only signifies that Industry Canada technical specifications were met.



"Operation of this device is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device."


"The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website [www.hc-sc.gc.ca/rpb](http://www.hc-sc.gc.ca/rpb)"

### 3. Features of the Controller

- (1) Aplio can be remotely controlled using Bluetooth wireless communication. (Depending on the surrounding conditions, the radiowaves have a reach of 20 m when there is no shielding.)
- (2) Each controller button is for a user function. Some functions of the operating panel and the touch panel can be assigned.
- (3) This controller can operate Aplio independently of the trackball. In addition, when controller communication is turned OFF, the GUI can be operated using the trackball, SET/NEXT switch, and F12 key of the keyboard. However, the body mark B activity cannot be operated with the controller OFF.  
For details of the body mark B, refer to the list of activities shown in the operation manual of Aplio, Application Volume.

## 4. External Appearance of the Controller and Buttons That Can Be Assigned

Function assigned to each controller button at the time of shipment

| No.  | Item            | Function equivalent to that of the Aplio operating panel                                                                                                        | Remarks                                                                                                                                         |
|------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| <1>  | Joystick        | Trackball                                                                                                                                                       | Assignment cannot be changed.                                                                                                                   |
| <2>  | Joystick button | SET                                                                                                                                                             | Other functions can be assigned.                                                                                                                |
| <3>  | Dial (+), (–)   | Palm dial                                                                                                                                                       | Other functions can be assigned.<br>(+) and (–) can be changed separately.                                                                      |
| <4>  | Button          | SET                                                                                                                                                             | Other functions can be assigned.                                                                                                                |
| <5>  | Button          | NEXT                                                                                                                                                            | Other functions can be assigned.                                                                                                                |
| <6>  | Button          | FREEZE/UNFREEZE                                                                                                                                                 | Other functions can be assigned.                                                                                                                |
| <7>  | Button          | HDD store                                                                                                                                                       | Other functions can be assigned.                                                                                                                |
| <8>  | Dial (+), (–)   | DEPTH/ZOOM                                                                                                                                                      | Other functions can be assigned.<br>The polarity of the field depth complies with the Aplio settings.<br>(+) and (–) can be changed separately. |
| <9>  | Dial button     | Switching between field depth and zoom                                                                                                                          | Other functions can be assigned.                                                                                                                |
| <10> | Button          | F12 key of the keyboard<br>Display/non-display of the menu                                                                                                      | Other functions can be assigned.                                                                                                                |
| <11> | Button          | Display/non-display of the popup menu<br>IASSIST controller dedicated function                                                                                  | Other functions can be assigned.                                                                                                                |
| <12> | Button          | Auto print                                                                                                                                                      | Other functions can be assigned.                                                                                                                |
| <13> | Wheel (+), (–)  | IASSIST menu operation (selection of activities)<br>Rotating upward (–) : Menu up<br>Rotating downward (+) : Menu down<br>IASSIST controller dedicated function | Other functions can be assigned.<br>(+) and (–) can be changed separately.                                                                      |
| <14> | Wheel button    | Execution of activities<br>IASSIST controller dedicated function                                                                                                | Other functions can be assigned.                                                                                                                |
| <15> | Power switch    | Controller power ON/OFF                                                                                                                                         | Assignment cannot be changed.                                                                                                                   |

## 5. Controller Operation Procedures

### 5.1 Checks Before Operation

Confirm that the Bluetooth-USB adaptor is inserted into the upper row of the USB ports at the back of the system.

Confirm that the controller battery is fully charged. Follow the procedures described in subsection 5.4 "LED Indicating Communication Status" for confirmation.

### 5.2 Usage Procedures

The functions assigned at the time of shipment are used in the descriptions below. (Refer to the figures shown in section 4 for the button numbers.)

| Step | Operation                                                       | Explanations                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|------|-----------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | Start of the controller<br>Press the switch <15>.               | <p>The LED indicating the communication status blinks.</p> <p>When communication with Aplio is established after approximately 12 seconds, this LED lights without blinking (refer to the NOTE on page 27). Simultaneously, the highlighted icon shown below is displayed at the left of the status bar on the screen.</p>  <p>After the switch is pressed there is a delay before the icon is displayed.</p> |
| 2    | Display/non-display of the GUI<br>Press the button <10>.        | Display and non-display can be switched every time the button is pressed.                                                                                                                                                                                                                                                                                                                                                                                                                          |
| 3    | Selection of activities<br>Operate the wheel <13>.              | The focus moves upward and downward on the menu.                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 4    | Execution of activities<br>Press the wheel button<14>.          | The selected activity is executed.                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| 5    | Display/non-display of the popup menu<br>Press the button <11>. | Display and non-display can be toggled by pressing the button.                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| 6    | Change of 2D gain<br>Rotate the dial <3>.                       | Operation is the same as for the palm dial on the operating panel.                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| 7    | Image freeze<br>Press the button <6>.                           | Operation is the same as for the FREEZE switch on the operating panel. Image freeze/unfreeze can be toggled by pressing the button.                                                                                                                                                                                                                                                                                                                                                                |
| 8    | HDD Store<br>Press the button <7>.                              | Operation is the same as for the HDD Store switch on the operating panel.                                                                                                                                                                                                                                                                                                                                                                                                                          |

| Step | Operation                                                                                                                                                                                                                             | Explanations                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 9    | Field depth/zoom<br>Change the field depth using the dial <8>.<br>Use the dial button <9> to switch to zoom.                                                                                                                          | Operation is the same as for the DEPTH/ZOOM dial on the operating panel.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 10   | Cursor operation<br>Operate the free cursor using the joystick <1>.                                                                                                                                                                   | Operation is the same as for the trackball on the operating panel.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| 11   | Setting of the cursor<br>Press button <2> or <4>.                                                                                                                                                                                     | Operation is the same as for the SET switch on the operating panel.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 12   | Measurement CH switching<br>Press the button <5>.                                                                                                                                                                                     | Operation is the same as for the NEXT switch on the operating panel.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| 13   | Output to the printer<br>Press the button <12>.                                                                                                                                                                                       | Operation is the same as for the AutoPrint switch on the operating panel.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| 14   | Communication OFF from the operating panel<br>If the controller cannot be operated for some reason while the communication is being established, the communication can be cut off by Aplio.<br><br>Press the F10 key on the keyboard. | <p>Communication between the controller and the operating panel is cut off. The IASSIST function can be operated using the trackball and the SET/NEXT switch.</p> <p>The icon at the left end of the status bar is displayed as shown below.</p>  <p>When the F10 key is pressed again, the grayed-out icon is displayed temporarily, communication is reset, and the icon is highlighted again (refer to the NOTE on page 27).</p> <p>After the F10 key is pressed there is a delay before the icon is highlighted (resumption of communication).</p> |
| 15   | Stopping of the controller<br>Press the power switch <15>.                                                                                                                                                                            | <p>The communication between the controller and the operating panel is cut off. The IASSIST function can be operated using the trackball and the SET/NEXT switch.</p> <p>The grayed-out icon is displayed at the left end of the status bar as shown below.</p>                                                                                                                                                                                                                                                                                       |

**NOTE:** The communication establishment time varies depending on the surrounding radiowave environment.

For example, if other equipment that uses Bluetooth is being used in the general vicinity, it is likely to take longer time to establish the communication. When the power supply of the controller is turned OFF while in the communication status and then turned ON immediately, it is likely to take some time for the next communication to be established. After the power supply is turned OFF, wait until the communication icon of Aplio is grayed out and the communication is cut off completely. Then, turn ON the power supply again.

To resume communication that was cut off using the F10 key, wait until the LED indicating the controller communication status starts blinking. Then, press the F10 key again.

### 5.3 Icon Indicating Communication Status (Displayed at the left end of the status bar on the screen.)

| Icon                                                                                | Operation                                                                                | Description                                                                                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|   | Communication establishment<br>Press the power switch <15> of the controller (power ON). | When communication is established, the LED indicating the communication status of the controller stops blinking and stays lit.                                                                                                                     |
|  | Stopping (communication cut-off from Aplio)<br>Press the F10 key on the keyboard.        | Communication stopping and resumption can be toggled by pressing the F10 key. When communication stopping is set, the LED indicating the communication status starts blinking. Pressing the F10 key again stops it blinking and makes it stay lit. |
|  | Communication termination<br>Press the power switch <15> of the controller (power OFF).  | The LED indicating the communication status of the controller goes out.                                                                                                                                                                            |

## 5.4 LED Indicating Communication Status (At the top-right of the controller)

|                                               | Status                                         | Description                                                                                                                                                                 |
|-----------------------------------------------|------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immediately after turning ON the power supply | Quick blinking approximately every 0.5 seconds | After turning ON the power supply, the LED blinks until communication starts.                                                                                               |
| Communication establishment                   | Continuously lit                               | The LED changes from quick blinking to being continuously lit in approximately 12 seconds.                                                                                  |
| Going out                                     | Going out                                      | When the power supply is turned OFF and the communication stops, the LED goes out.                                                                                          |
| Blinking for battery                          | Slow blinking approximately every 2 seconds    | When the LED enters this status, the power supply is automatically turned OFF in approximately 30 minutes.<br>Immediately charge the battery using the charger provided.(*) |

(\*) Put the controller in the charger after every use and charge it to prevent the battery from suddenly running down.

## 5.5 Change of Button Assignment






| Step | Operation                                                                                                                                                                                                                                                                                                                                                                                    | Description                                                                                                                                                                        |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | <p>Display the popup menu and perform either of the operations described below.</p> <p>(1) Press the NEXT switch on the operating panel while the controller is OFF.</p> <p>(2) Press the button &lt;11&gt; while the controller is ON.</p>                                                                                                                                                  | The sub menu (pull-down menu) is displayed.                                                                                                                                        |
| 2    | <p>Selection of the controller button to be assigned</p> <p>Open the setting screen by performing either of the operations described below.</p> <p>(1) Select the "DevAssign" switch using the trackball and the SET switch while the controller is OFF.</p> <p>(2) Select the "DevAssign" switch using the wheel &lt;13&gt; and the wheel button &lt;14&gt; while the controller is ON.</p> | <p>The controller assignment screen opens.</p> <p>The cursor is displayed on the screen.</p>     |
| 3    | <p>Selection of the button to be changed</p> <p>Select the button to be changed using the trackball and the SET switch.</p> <p>For the list of the buttons that can be changed, refer to the next subsection.</p>                                                                                                                                                                            | <p>The currently assigned type, group, and function of the selected button are displayed.</p>  |

Figure 2 Type, Group, and Function

| Step | Operation                                                                                                                                                                                                                                                                                                                                                                | Description                                                                                                                                                                                                                                                     |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 4    | <p>Selection of the function to be assigned</p> <p>Select the desired item from the Type, Group, and Function categories.</p>                                                                                                                                                                                                                                            | <p>When the function is selected, the function display on the controller assignment screen is automatically changed.</p> <p>At this point, only the display is changed.</p> <p>The function assignment is not yet set.</p>                                      |
| 5    | <p>Applying the changes</p> <p>Press the Apply button using the trackball and the SET switch.</p> <p>To return the changes to their original settings, press the Quit button.</p>                                                                                                                                                                                        | <p>The Apply button is grayed out.</p> <p>Pressing the Quit button cancels the changes and closes the screen.</p>  <p style="text-align: center;">Figure 3 Apply and Quit</p> |
| 6    | <p>Termination of the controller assignment screen</p> <p>Press the Quit button at the lower right of the screen.</p> <p>When Quit is pressed without setting the changes, a dialog box is displayed asking whether the changes should be saved or the screen closed without saving (figure 4).</p> <p>When the screen is closed, the popup menu is displayed again.</p> |  <p style="text-align: center;">Figure 4 Confirmation of saving/canceling of the changes</p>                                                                                |

| Step | Operation                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 7    | <p>&lt;Exception processing&gt;</p> <p>When any of the functions given below is removed from the controller, a confirmation message is displayed (figure 5).</p> <ol style="list-style-type: none"> <li>1.Selection of activities (Roll Up)</li> <li>2.Selection of activities (Roll Down)</li> <li>3.Execution of activities</li> <li>4.F12 (PAS GUI display)</li> <li>5.Popup menu display</li> <li>6.SET</li> <li>7.NEXT</li> </ol> <p>(Important)</p> <p>To assign the function again, turn OFF the power supply of the controller to enable the trackball and the SET/NEXT switch on the operating panel.</p> <p>During communication with the controller, IASSIST operation can not be performed using the trackball and the SET/NEXT switch.</p> | <p>When the OK button is selected on the dialog box and any of the functions shown at left is removed, IASSIST operation cannot be performed using the controller.</p>  <p>Figure 5 Confirmation of deletion/change of important functions</p> <p>When the Cancel button is selected without applying the changes, the dialog box is closed and the controller assignment screen is displayed again.</p> |

\* The same function can be assigned to multiple switches.

## 5.6 List of the Buttons of the Controller

Refer to the figures in section 4 for the button No.

| Button No. | Controller button type | Location of the switch to be assigned at Aplio | Remarks                                                  |
|------------|------------------------|------------------------------------------------|----------------------------------------------------------|
| <1>        | Joystick               | Cannot be changed.                             | Fixed to the trackball.                                  |
| <2>        | Joystick button        | GUI/Panel/TCS                                  | The details of the functions are shown in the next list. |
| <3>        | Dial                   | GUI/Panel                                      |                                                          |
| <4>        | Button                 | GUI/Panel/TCS                                  |                                                          |
| <5>        | Button                 | GUI/Panel/TCS                                  |                                                          |
| <6>        | Button                 | GUI/Panel/TCS                                  |                                                          |
| <7>        | Button                 | GUI/Panel/TCS                                  |                                                          |
| <8>        | Dial                   | GUI/Panel                                      |                                                          |
| <9>        | Dial button            | GUI/Panel/TCS                                  |                                                          |
| <10>       | Button                 | GUI/Panel/TCS                                  |                                                          |
| <11>       | Button                 | GUI/Panel/TCS                                  |                                                          |
| <12>       | Button                 | GUI/Panel/TCS                                  |                                                          |
| <13>       | Wheel                  | GUI/Panel/TCS                                  |                                                          |
| <14>       | Wheel button           | GUI/Panel                                      |                                                          |
| <15>       | Power switch           | Cannot be changed.                             | Fixed to the power switch.                               |

## 5.7 Functions Which Can Be Set to the Controller


[illegible]

| Type        | Group               | Function (in alphabetical order)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | Remarks                           |
|-------------|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------|
| Panel       | Nothing             | Acoustic Power (-)/(+)<br>Angle Cor. (-)/(+)<br>Angle Cor. (-5)/(+5)<br>Base Line (-)/(+)<br>Base Line (-5)/(+5)<br>Color Filter (-)/(+)<br>Color Frequency (-)/(+)<br>Color Gain (-)/(+)<br>Color Gain (-5)/(+5)<br>Color Steering (-)/(+)<br>D Gain (-)/(+)<br>D Gain (-5)/(+5)<br>D Gate Size (-)/(+)<br>D Gate Size (-4)/(+4)<br>Depth/Zoom (-)/(+)<br>DOP Filter (-)/(+)<br>Meas. Edit (-)/(+)<br>PW Steering (-)/(+)<br>R-Delay1 (-)/(+)<br>R-Delay1 (-10)/(+10)<br>R-Delay1 (-30)/(+30)<br>R-Delay2 (-)/(+)<br>R-Delay2 (-10)/(+10)<br>R-Delay2 (-30)/(+30)<br>Scale (-)/(+) | Can be assigned to a dial only.   |
|             | ↓                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |                                   |
| Touch panel | Base (Common)       | CH1 Delay<br>CH2 Delay<br>OTHER<br>PIMS<br>PRESET<br>SPECIAL_APPLI                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | Can be assigned to a button only. |
|             | PIMS                | Auto Print<br>B/W Print<br>Color Print<br>DICOM Print<br>Print Size<br>Still<br>Still Destinat.<br>Still Size                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | Can be assigned to a button only. |
|             | ↓                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |                                   |
|             | Special Application | Dynamic Flow<br>Panoramic View<br>TDI                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | Can be assigned to a button only. |

(\*) Operation differs depending on the output destination set in the operating panel.

## 6. Charger

### 6.1 External Appearance

### 6.2 Charging Time

It takes approximately 2 hours and 30 minutes to charge a new battery after it has been fully discharged. The charging time varies depending on the battery age and memory effect.

### 6.3 LED Indicating Charging Status (At the lower middle of the controller)

When the controller is inserted in the charger, the LED indicating the charging status at the lower middle of the controller lights in green. When charging is finished, this LED goes out. When the controller is removed after charging and then inserted again, the LED lights but charging is completed immediately and the LED goes out.

### 6.4 Power Switch

Turning ON the power switch at the back of the battery enables charging. Leave this switch turned ON normally.

## 7. Bluetooth-USB Adaptor

- CAUTION:**
1. Do not connect/disconnect the Bluetooth-USB adaptor during system startup/shutdown. The system may not start up correctly.
  2. Do not remove the Bluetooth-USB adaptor during communication. Communication cannot be established. If so, turn OFF the power supply of the controller, insert the adaptor, and then perform connection again.
  3. Do not apply excessive force to the Bluetooth-USB adaptor from the above when it is inserted into the USB port. The Bluetooth-USB adaptor may be damaged.

### 7.1 External Appearance

### 7.2 Operating Procedures for the Bluetooth-USB Adaptor

Insert the Bluetooth-USB adaptor into the upper row of the USB ports.

Also, if the diagnostic ultrasound system is ON, the power switch of the controller must be set to OFF when the Bluetooth-USB adaptor is connected. Similarly, the power switch of the controller must be set to OFF when the adaptor is disconnected. The adaptor should be connected/disconnected only when necessary.