

## **READ BEFORE USING THESE INSTRUCTIONS!**

1. The pants can be machine washed and dried. **Do not wring!**
2. The receiver works on electricity **or** batteries. **DO NOT USE IN COMBINATION.**
3. It is advisable to place the receiver on some distance from the child, so that he or she **has** to get out of bed to turn off the receiver.
4. Other nearby transmitters can disturb the bedwetting alarm system. (For example: wireless baby-monitor). Change the frequencies of the other devices or move them if needed.
5. Always remove the transmitter immediately from a wet pant. This extends the life of the battery. Leaving the transmitter on a wet pant might lead to slight discoloration of the pants.
6. Make sure to remove the transmitter from the pants before laundering. **The warranty is voided if this occurs.**
7. The transmitter and receiver are programmed to be used together. If either part is not working, both must be replaced.
8. Remove the batteries when the bedwetting alarm system is not going to be used for more than several weeks.
9. Do a system test before using for the first time to acquaint yourself with the operation and the sound of the unit.
10. The socket-outlet shall be installed near the equipment and shall be easily accessible.

### **System Test**

- Put the receiver into an electrical outlet and turn it on.
- The green light should shine continuously and the red light should flash once when you turn on the unit.
- Connect the two posts of the transmitter with a metal object (scissors, screwdriver). This will cause the wall unit to sound.
- The wall unit can be turned off by pushing the small button on the front.
- When the receiver does not go off, then check the battery in the transmitter.

### **Pants Test**

- Take a cup of lukewarm water (body temperature).
- Add a wiped teaspoon of salt, and stir.
- Click the transmitter on the pant.
- Put the receiver in a socket.
- Pour two teaspoons of this salty solution on the front of the pant. This should cause the wall unit to sound.
- After the pants are washed, they become more absorbent.

## **FCC Compliance**

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.

**Warning:** Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for help.

### **Responsible Party in USA**

Legal Company Name: Enuresis Associates LLC  
Address: 17737 New Hampshire Ave  
Ashton MD20861  
USA

Phone, Fax: 001301-774-1349 001301-774-7982  
Contact Person: Bryan Mercer