ABU62

User's Guide Guide d'Utilisation

Contents

Introduction	3
General Information	4
Setting up the titration system	6
Maintenance	9
Cleaning	9
Transporting the instrument	9
Servicing	9
Replacing the fuses	10
Technical specifications	11
International Standards	14

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Introduction

The ABU62 is a compact unit integrating two motor-driven burettes plus two electrode and one temperature inputs on a single platform. Two ABU62s can be connected to your existing TitrLab system providing up to six motor-driven burettes and six electrode inputs, meeting virtually all application needs. It can also be used as a sample stand or for electrode storage and reconditioning. The ABU62 offers a large screen display, giving you key information at a glance such as the reagent dispensed volume at each burette position.

The mono-block burette stand operating in conjunction with a high-resolution step-motor offers more than 18,000 steps to dispense a full burette volume in no more than 20 seconds. It also offers a resolution for volume addition down to 0.1 μ l for the 1 ml burette. The mono-block burette stands are fitted in no time and come in five sizes: 1, 5, 10, 25 and 50 ml.

General Information

Safety information

Please read this entire manual before unpacking, setting up, or operating this equipment.

Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

To ensure that the protection provided by this equipment is not impaired, do not use or install this equipment in any manner other than that specified in this manual.

Precautionary Labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	<p>Electrical and electronic equipment marked with this symbol may not be disposed of in European public disposal systems after 13 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of life equipment to the Producer for disposal at no charge to the user.</p> <p>Note: For equipment supplied or produced by "Radiometer Analytical", please contact www.hach-lange.com and select your country for instructions on how to return your equipment for proper disposal.</p>
	This symbol, when noted on the product, identifies the location of a fuse or current limiting device.

Compliance with FCC rules, part 15

Information to the user

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception which can be determined by turning the equipment off and on, the user is encouraged to try to correct interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Warning!!

Changes and modifications that are not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Setting up the titration system

Warning !


The ABU62 have been developed to meet the requirements of electrochemical applications. It is therefore aimed at experienced users who have the knowledge required to operate the instrument and implement the security instructions enclosed. Please remember that the ABU62 must not, under any circumstances, be used to perform tests on living beings.

We accept no responsibility for using the ABU62 and its peripheral devices under conditions that are not specified in this User's Guide and its associated TitraLab 9xx Reference manuals (part no. D21T073 for a TitraLab 960/965 or D21T075 for a TitraLab 980).

The ABU62 is used in association with the following Titration Managers: TIM960, TIM965 and TIM980. All Titration Managers or ABU62 in the system should have the same version of embedded software installed preferably from the Radiometer Analytical www.titration.com resource centre.


1. Set up the instrument in a properly ventilated place. The power supply connector on the rear panel must remain accessible and close to the user (2 m maximum) so that you can quickly disconnect the cables in case of emergency.
The room temperature should not exceed 40°C.

2. Connect the ABU62 **PC/TIM** socket to the Titration Manager Local socket. Use cable, part no. A95X501. Two ABU62's can be controlled by the same Titration Manager. Connect the instruments as follows:

3. Connect the mains socket of the ABU62 and the Titration Manager to the mains supply using the 3-lead power cord provided. The ABU62 and the Titration Manager must be connected to an earthed mains socket for safety reasons. Efficient grounding is essential for reliable measurements and security. In the USA or Canada, use a UL listed power cord only.
4. Switch on the ABU62 (O/I power switch set to "I") then the Titration Manager. The ABU62 display is as follows:

5. Proceed as follows to adjust the brightness of the ABU62 display from the Titration Manager:

- Display the Cell window.
- At the line Instrument, select the ABU62 (ABU1 or ABU2) you want to adjust the display brightness. ABU1 means an ABU62 connected to the **Local** socket of the Titration Manager. ABU2 means an ABU62 connected to the **Slave** socket of a second ABU62 of the ABU1 type (see figure on the previous page).
- Press 0 to increase the brightness or 7 to decrease the brightness.

The ABU62 can only be operated from the Titration Manager ("Slave" mode). It means that:

- . The Configuration parameters of the Titration Manager are imposed to the ABU62 (language, date and time in particular),
- . The electrode and reagents used on the ABU62 are created and edited from a Titration Manager screen,
- . The burettes, electrodes and reagents used on the ABU62 are installed from a Titration Manager screen,
- . Direct measurements using electrodes connected to the ABU62 are initiated on the Titration Manager,
- . Manual dosing using ABU62 burettes is initiated on the Titration Manager.

All these operations are described in the Reference Manual of the Titration Manager: consult the "Glossary" chapter at the appropriate keyword index (examples: Create electrode, Electrode calibration, Manual dosing).

Maintenance

Cleaning

The ABU62 requires minimum maintenance. The exterior surface can be cleaned with tepid water and wiped dry with a soft cloth. Never use another solvent unless you have first consulted your Radiometer Analytical representative

Transporting the instrument

Always use the packaging supplied by the manufacturer.

IMPORTANT!

Remove the metal rod before transporting the instrument.

Never pick-up or carry the instrument by the metal rod.

Servicing

DO NOT ATTEMPT TO SERVICE THIS PRODUCT YOURSELF, except as noted in the Reference Manual. For servicing, please contact your **Radiometer Analytical service representative** at the address given below:

RADIOMETER ANALYTICAL SAS
72, rue d'Alsace
69627 Villeurbanne CEDEX - France
Tel.: +33 (0) 4 78 03 38 38
Fax: +33 (0) 4 78 03 38 27
E-mail: radiometer@nalytical.com

or your local service representative:


Replacing the fuses

For continued protection replace the fuse with one of a high interrupting capacity, same type and rating:



2 x fuses, slow blow, 1.0 A (5 x 20 mm), part no. 450-020.

To replace the fuses:

1. Disconnect line cord

2. Remove the fuse holder

3. Replace the used fuses with ones of the same type and rating

4. Put the cap back in place

Technical specifications

Burette

2 detachable burette stands.

Burette volumes available: 1, 5, 10, 25, 50 ml

Burette motor: 18,000 steps.

Complies with ISO 8655-3

UV-protected encapsulated glass syringe.

Operating procedures for burette exchange, air bubble removal (Flush), Rinse, Fill, Empty function all controlled from the Titration Manager

Second burette operational during first burette refill

Inputs/outputs

2 indicator electrode inputs.

1 reference electrode input.

Selectable polarised input from -1 mA to 1 mA in 1 μ A steps, DC or AC

Differential input.

Temperature input.

Serial connections to titration workstation/PC and second ABU62 Biburette.

Measuring ranges

-9 to 23 pH.

± 2000 mV.

-10 to 100°C.

Resolution

0.001 pH.

0.1 mV.

0.1°C.

Languages

English, German, Danish, French, Italian, Spanish and Swedish.

Casing

Fully splashproof latene®.

Graphic 128x128 dot LCD (TPX®).

Dimensions (H x W x D)

380 x 230 x 450 mm (excl. tubing).

Weight

5 kg (excl. reagent bottles).

Power requirements

100-240V, 50/60Hz, 26W.

Environmental operating conditions

5 to 40°C temperature.


20 to 80% relative humidity.

Level of pollution: 2.

Transitory overvoltage: level II.

Fuses: 2 mains fuses : slow blow 1.0 A (230 V)
 1 secondary fuse: 4 A.

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CE Marking DECLARATION OF CONFORMITY

We, **Radiometer Analytical SAS**,

declare under our sole responsibility that the Potentiometric Titration Workstations **TIM9xx**, including TIM960, TIM965, TIM980 and ABU62

are in conformity with the provisions of:

Radio and Telecommunications equipment directive 99/5/EEC, dated 9/3/1999,

Low voltage directive 73/23/EEC, dated 19/2/73

EMC directive 89/336/EEC, dated 3/5/1989, 92/31/EEC, dated 28/4/1992,

CE harmonisation directive 93/68/EEC, dated 22/7/1993,

following standards and severity levels of:

EN 61010-1, 1995

EN 300330-1, 1999,

ETSI TR102 070, **2002**, and **EN 301 489**, **2001**,

CISPR16-1, 1999, and **EN 61326-1**, 1997,

EN 55022, 1998, class A,

EN 61000-3-2, 2000, class A and **EN 61000-3-3**, 1995,

EN 61000-4-2, 1995, level 2 (4kV) with contact discharges and level 3 (8kV) with air discharges,

EN 61000-4-3, 2002, level 2 (3V/m),

EN 61000-4-4, 1995, level 2 (1kV) on AC power line,

EN 61000-4-5, 1995, level 2

EN 61000-4-6, 1996, level 2 (3V),

EN 61000-4-11, 1994,

after tests and qualification of our products, regarding "Radio" and "EMC" directives, performed by:

AEMC Lab, Sassenage (France) independent laboratory, **Cofrac** accredited under the number 1-0805, with a **Radio qualification report no. R0603066R-E-C** and a **EMC qualification report no.**

R0603066C-E-C, dated June 15, 2006 which attests **the conformity** of the TIM9xx workstations,

and tests and qualification of the product, regarding "Low Voltage" directive, attested by our company, according to the self-assessment procedure (article 10.1), recommended for the CE marking, as these products are UL / CSA certified.

Year of the first CE marking: **2006**.

Villeurbanne, July 20, 2006
RADIOMETER ANALYTICAL SAS

Georges RIVOIRARD
Quality Manager

The Titration Manager complies with the following standards:

UL 61010A - 1

CAN / CSA C22 2 N° 1010.1 - 92.

The ABU62 complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

See "Compliance with FCC rules, part 15 - Information to the user" on page 5.

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