


Instructions for Use

Confidant Connector 2.0 Serial/Bluetooth Converter



Description

The Confidant Connector Serial/Bluetooth Converter collects data from monitoring devices over a serial port and then transmits that data via the Bluetooth wireless communications protocol to another Bluetooth device, such as a mobile phone.




Installing the Batteries

The battery chamber is accessed by opening the battery door on the back of the device. Gently slide the door to one side, as indicated by the Open arrow and lift upwards. The inside of the battery compartment is marked for correct battery orientation. The battery must be installed in the correct orientation for the device to operate. The device requires two AAA alkaline batteries.



Plugging in the Data Cable

The Data Cable has different sized Jacks at each end; one small (2.5mm) and one large (3.5mm). This is to accommodate a range of monitoring devices with different Socket sizes. Note that either Jack can be connected to the Confidant Connector using the corresponding Socket.



Starting the Confidant Connector

The Confidant Connector is started by sliding the ON Switch forward and releasing. The Indicator Light will come on when the device is ready to communicate. If the device is idle for 120 seconds it automatically shuts off to conserve battery life.

Note: the device will shut off after 120 seconds of its initial start up unless communications are initiated. Repeatedly pushing the ON Switch will not keep the Confidant Connector on beyond the initial 120 seconds.

Resetting the Confidant Connector

In the rare event that problems develop with communication, the Confidant Connector device can be reset by sliding the ON Switch forward and releasing. This restarts the Confidant Connector device.

Indicator Light

The Indicator Light flashes on startup and then stays on to indicate the device is ready for communications. The Indicator Light also flashes when communications are taking place. The color of the Indicator Light indicates battery strength.

Battery Strength		
GREEN = Good	ORANGE = Weak	RED = Dead

Troubleshooting

If you experience any problems with the Confidant Connector device try resetting it. It is ready to communicate when the indicator light stays on. If you continue to have difficulty, double check the batteries to make sure they are correctly oriented. Replace the batteries when they become weak.

For Technical Support, Contact:

Confidant International, LLC
2530 Meridian Parkway, Suite 300
Durham, NC 27713
Phone: (800) 886-4220
Email: support@confidantinc.com

WARNINGS

Do not immerse the Confidant Connector device in water or other liquids. Do not expose the Confidant Connector device to excessive heat or cold.

The information contained in this manual is subject to change without notice. Confidant International, LLC reserves the right to make changes to equipment design or program components, as progress in engineering, manufacturing methods, or other circumstances may warrant. No responsibility is assumed for the use of equipment or software that is not supplied by Confidant International, LLC.

Federal Communications Commission (FCC) Statement

15.21

You are cautioned that changes or modifications not expressly approved by Confidant International LLC could void the user's authority to operate the equipment.

15.105(b)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Operation is subject to the following two conditions:

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.