

ALVA BC640

User Manual English

UDC-02004 A

© 2007 Optelec, the Netherlands.

All rights reserved



Optelec

P.O. Box 399

2990 AJ Barendrecht

The Netherlands

Telephone: +31 (0)88 6783 444. Fax +31 (0)88 6783 400

Internet: www.optelec.com

Contents

1. Introduction.....	4
1.1. About this manual.....	4
1.2. ALVA BC640 package contents.....	4
2. Getting to know your ALVA BC640	6
2.1. ALVA BC640 Base Unit: Surface and front panel layout	6
2.1.1. Surface and front panel layout.....	6
2.1.2. Back panel layout	7
2.1.3. Right side layout	7
2.1.4. Left side layout.....	7
2.1.5. Front panel keys, Smartpad, eTouch navigation buttons and cursor routing keys	7
2.2. The ALVA BC640 Braille Audio Feature Pack	8
2.2.1. Braille Audio Feature Pack: Surface layout	9
2.2.2. Right Side layout.....	9
2.2.3. Left Side lay-out.....	9
3. Installation	10
3.1. Setting up the USB connection.....	10
3.2. Installing ALVA BC640 screenreader drivers.....	10
3.3. Bluetooth® connection.....	11
3.3.1. Preparing for Bluetooth® connection	11
3.3.2. Setting up the Bluetooth® USB dongle.....	11
4. The ALVA BC640 Local menu.....	15
4.1. Activating the Local Menu.....	15
4.2. Changing settings in the Local Menu	15
4.3. The Feature Pack Menus	23
4.3.1. The Braille Audio Feature Pack Bluetooth® Menu.....	23
4.3.2. The Braille Audio Feature Pack Audio Menu	24
5. Power Sources	25
5.1. USB power	25
5.2. Internal battery.....	25
5.3. External power adapter.....	25

6.	The ALVA BC6 Bluetooth® settings application	26
6.1.	Configuring Audio	26
6.2.	Using the Braille keyboard	26
7.	Using the Braille keyboard.....	28
7.1	Chord commands.....	28
7.2	Prefix chord commands	28
7.3	Contracted Braille on/off.....	29
8.	Notes	30
8.1	Starting Notes	30
8.2	Using Notes	30
8.2.1	Adding text.....	30
8.2.2	Editing text.....	30
8.2.3	Navigating in Notes.....	31
8.2.4	Defining a block	31
8.2.5	Reading in Notes.....	31
8.3	The Notes Menu	32
8.3.1	Notes main menu.....	32
8.3.1.1	File menu.....	32
8.3.1.2	Edit menu	33
8.3.1.3	Help	33
9.	Updating the Firmware	35
10.	Mounting the Braille Audio Feature Pack.....	36
	Regulatory notice.....	37
	Optelec offices worldwide	41

1. Introduction

Congratulations on choosing the ALVA BC640. The ALVA Braille Controller 6-series is the first series of wireless Braille Controllers that merges classic Braille display access with modern note taking functionality. The high quality and reliability of the ALVA BC640 gives you optimal flexibility and productivity in an ultra compact design.

The ALVA BC640 features future proof hardware standards as wireless Bluetooth® connection and USB ports. The ergonomic design with its high quality Braille cells and optimal navigation possibilities allow you to be in full control of your computer, PDA or (smart) phone. Used in combination with screen reading software the ALVA Braille Controller can be installed for particular operating systems to provide a tactile and audible user interface.

The ALVA BC640 offers an enhanced ergonomic design that extends the navigation and operation possibilities of the screen reading software. The ergonomics are designed for optimal efficiency when working with a computer, PDA or Smartphone. Note that the functionality of the ALVA BC640 depends somewhat on the Braille display options of your screenreader and your mainstream computing hardware.

The ALVA Braille Controller was designed to adhere to the Open Standard philosophy using mainstream technology. This enables compatibility with most screen reading software and operating systems. Contact your distributor or Optelec for updated information about screenreaders that fully support the ALVA BC640.

If you have any questions or suggestions concerning the use of this Braille display, please contact your distributor or Optelec's head office using the contact information at the last page of this manual. Your feedback is highly appreciated. We hope you enjoy working with your ALVA BC640.

1.1. About this manual

This manual will familiarize you with the features and basic operation of the ALVA BC640 and Braille Audio Feature Pack. Please read this manual carefully before using your device.

1.2. ALVA BC640 package contents

The ALVA BC640 package contains the following:

- A pouch
- The ALVA BC640 base unit
- Power adapter 110/240V and power cable
- Long and short USB cables
- An USB Bluetooth® dongle, Class 1, Windows XP SP2 and Windows Vista compatible
- The ALVA BC driver software for Windows XP and Windows Vista and documentation on CD-ROM
- This manual in Print

-
- Quick Start Guide in Braille and Print

The ALVA BC640 with Braille Audio Feature Pack package contains the following:

- A carrying bag
- The ALVA BC640 base unit including Braille Audio Feature Pack
- Power adapter 110/240V and power cable
- Long and short USB cables
- An USB Bluetooth® dongle, Class 2, Windows XP SP2 and Windows Vista compatible
- The ALVA BC driver software for Windows XP and Windows Vista and documentation on CD-ROM
- A headset
- This manual in Print
- Quick Start Guide in Braille and Print

The Braille Audio Feature Pack package contains the following:

- A carrying bag
- The Braille Audio Feature Pack
- The ALVA BC driver software for Windows XP and Windows Vista and documentation on CD-ROM
- A headset
- A hex key
- This manual in Print
- Quick Start Guide in Braille and Print

If any of these items are missing from your package, please contact your distributor.

2. Getting to know your ALVA BC640


The ALVA BC640 consists of an ALVA BC640 base unit containing a 40 cell Braille line and an optional Braille Audio Feature Pack. This chapter will give a detailed overview of the lay-out of the ALVA BC640 base unit and Braille Audio Feature Pack.

2.1. ALVA BC640 Base Unit: Surface and front panel layout

The unit is placed correctly in front of you when the 40 Braille cells are located at the back of the flat surface and the thumb keys are facing towards you.

The ALVA BC640 contains:

1. 40-cell Braille display
2. eTouch keys (left) for panning left and right and screenreading functions
3. eTouch keys (right) for panning left and right and screenreading functions
4. Smartpad with 5-way navigation key and 4 keys for easy access to mainstream devices
5. Thumb keys for panning, left, right, up and down and returning to the cursor position
6. 40 Cursor Routing Keys
7. Power supply and mini USB connectors
8. On/off switch



2.1.1. Surface and front panel layout

The ALVA BC640 is designed with a row of 40 Braille cells at the back of the flat surface. In front of the Braille cells one row of small buttons can be found, known as Cursor Routing Keys. Each Braille cell has a corresponding Cursor Routing Key that can be used for multiple functions. On the left and right side of the Braille line you will find two groups of two keys. These so-called eTouch keys can be used for Braille panning left and right as well as various other screen-reader functions.

In the middle, just in front of the Cursor Routing Keys, you will find a group of 9 keys. This group is called Smartpad and consists of 4 function keys and a 5-way navigation key. The Smartpad can be used for operating mainstream devices such as PC's, (smart)phones and PDA's.

2.1.2. Back panel layout

At the back of the Braille Controller, a connector can be found for the Braille Audio Feature Pack.

2.1.3. Right side layout

The power switch is located at the right side of the Braille Controller. The Braille Controller is switched on when the rectangular push button is pressed. Switching the Braille Controller on will generate an audio signal.

2.1.4. Left side layout

The left side of the Braille Controller contains three connectors. From front to back: power connector, mini USB connector and a feature connector. The feature connector has no function with the release of this manual.

2.1.5. Front panel keys, Smartpad, eTouch navigation buttons and cursor routing keys

The names of the ALVA BC640 keys listed in this section will be used in all documentation, including your screenreader documentation. This section describes the name of the key(s) and its main function, which is supported by the screenreader. The exact functionality is defined by the screenreader. Please refer to your screenreader documentation for detailed functionality of the ALVA BC640.

Smartpad

The integrated Smartpad allows for simultaneous access to mainstream devices such as PC's, PDA's and (smart)phones and navigation through menu's and dialog boxes. From left to right, the Smartpad keys are:

Smartpad key 1	First rectangular key at the left marked with 3 horizontal lines
Smartpad key 2	Second rectangular key at the left marked with 1 horizontal line
LEFT	Left key on the 5-way navigation key
UP	Top key on the 5-way navigation key
DOWN	Bottom key on the 5-way navigation key
RIGHT	Right key on the 5-way navigation key
MIDDLE	Middle key marked with a round dot of the 5-way navigation key
Smartpad key 3	Second rectangular key from the right marked with 1 vertical line
Smartpad key 4	Most right rectangular key marked with 3 vertical lines

eTouch keys

The two keys on each side of the Braille line are called the eTouch keys. They are used for panning left and right as well as various other Braille functions.

eTouch 1	Top left key marked with two vertical lines
eTouch 2	Bottom left key marked with one horizontal line
eTouch 3	Top right key marked with two vertical lines
eTouch 4	Bottom right key marked with one horizontal line

Cursor Routing Keys

The cursor routing keys correspond to the particular Braille cell immediately above each button. Typically, they are used for routing the cursor or caret to the corresponding Braille cell, generating a left or right mouse click at that point, or for obtaining additional information about the contents of the Braille cell.

Thumb Keys


At the front of the Braille Controller five front panel keys are placed, the so called thumb keys. Four bigger keys are symmetrically aligned along the front of the ALVA BC640. One smaller button is located in the middle of the front. Pressing one of the four bigger thumb keys will navigate left, right, up or down. By pressing the middle key you will return to the cursor position.

2.2. The ALVA BC640 Braille Audio Feature Pack

The Braille Audio Feature Pack is placed correctly in front of you when the Braille input keys are located at the back of the flat surface and the thumb keys of the ALVA BC640 Base Unit are facing towards you.

The ALVA BC640 Braille Audio Feature Pack contains:

1. Ctrl-key
2. Eight key Braille keyboard (From left to right Braille key 7, 3, 2, 1, 4, 5, 6, 8)
3. Windows key (to open the Windows Start Menu)
4. Enter key
5. Alt-key
6. Headphone connector
7. Microphone connector
8. Space bar
9. Volume control and mute button
10. Integrated high quality stereo speakers



2.2.1. Braille Audio Feature Pack: Surface layout

The ALVA BC640 Braille Audio Feature Pack is designed with an eight key Braille keyboard that is located directly behind the 40 cell Braille line. The Braille keyboard consist of eight Braille keys, a space bar, a Ctrl-key, a Windows key, an Enter key and an Alt-key and two integrated high quality stereo speakers that are located left and right of the Braille keyboard. The space bar is the long key that is located in the centre of the Braille Audio Feature Pack. The eight rectangular keys located left and right of the space bar are the Braille input keys. The four square keys are the Ctrl-key, Windows key, Enter key and Alt-key. From left to right the keys are labelled Ctrl-key, Braille key 7, Braille key 3, Braille key 2, Braille key 1, Windows key, Enter key, Braille key 4, Braille key 5, Braille key 6, Braille key 8 and the Alt-key.

2.2.2. Right Side layout

The volume and mute button is located at the right side of the Braille Audio Feature Pack. Moving the button towards you allows you to turn the volume down, whereas moving the button away from you will raise the volume. Pressing the button allows you to mute the audio, which will be generated with a low beep. Pressing this button again will generate a higher beep and switch the mute function off. A message will be displayed on the Braille line to indicate if the mute function is on or off.

2.2.3. Left Side lay-out

The left side of the Braille Audio Feature Pack contains two connectors. From back to front you will find a 3.5 mm headphone connector and a 3.5 mm microphone connector.

3. Installation

In this section, a general description of the installation procedure is described for the JAWS and Dolphin screenreaders for Windows XP and Windows Vista, the Braille Audio Feature Pack setup and Bluetooth® installation. Please refer to the screenreader documentation for instructions on how to setup and configure the Braille display for other screenreaders and operating systems.

In case you upgraded your ALVA BC640 base unit to an ALVA BC640 with Braille Audio Feature Pack, please make sure to upgrade the firmware to version 1.14.0 or higher and Bluetooth® version 1.6.0 and 1.4.0. To learn more about the firmware version, please refer to chapter 4. The firmware update procedure is described in chapter 8.

3.1. Setting up the USB connection

The ALVA BC640 is compatible with Windows XP and Windows Vista. Older Windows versions are not supported. The ALVA BC640 does not need any special USB drivers. Windows XP or Windows Vista will automatically detect the ALVA BC640 and install the necessary drivers.

1. Attach the USB cable to the ALVA BC640 and computer.
2. A 2GB removable drive will appear on the screen. This drive contains the ALVA BC640 screenreader drivers for JAWS and Dolphin.
3. Switch on the ALVA BC640. Windows will automatically install the necessary drivers. Generally this does not require any further user input. When the connection is established the ALVA BC640 will issue an audible signal and displays the message “USB connected”.
4. If you use an ALVA BC640 with Braille Audio Feature Pack, the audio will be rerouted to the Braille Audio Feature Pack and the Braille keyboard will be installed. Please refer to chapter 6 for configuration instructions of the USB audio and eight key Braille input.

3.2. Installing ALVA BC640 screenreader drivers

1. Insert the ALVA BC640 software CD-ROM that came with the package in your computer's disk drive.
2. The program will automatically run if the autostart feature of your CD-ROM drive is enabled.
 - a. If the autostart feature has not been enabled, open your Windows Start Menu, select 'Run' and type d:\setup.exe (assuming d: is the letter assigned to your disk drive) and press Enter.
 - b. Alternatively, locate the ALVA_BC640 folder on the ALVA BC640 removable drive and run the “ALVA BC6 Setup.exe”.

3. Select "Install screenreader drivers"
4. The setup program will guide you through the installation.
5. The setup program will automatically search your system for installations of JAWS and Dolphin screenreaders and will install all necessary drivers, manuals, firmware update software and the application to configure the ALVA BC640 Braille Audio Feature Pack Bluetooth® settings. Please note that the ALVA BC640 is only supported by Windows XP and Windows Vista.

Drivers for screenreaders other than Dolphin and JAWS are included in the products themselves and are not shipped on the ALVA BC640 CD-ROM. Please refer to your screenreader documentation for more information about the setup of the ALVA BC640. You have the choice between two setup types: a complete and a custom setup. We recommend using the "Complete" setup. The "Complete" setup will automatically install and update all the necessary components. If you would like to install the ALVA BC640 software manually you can choose the "Custom" setup option. This option is for advanced users only. The "Custom" setup will ask you what components you would like to install. For each component chosen you will need to provide the correct installation location.

Note:

The USB audio channel in Windows is configured to mix the microphone audio with the playback audio. This can result in a beeping sound when connection the Braille Audio Feature Pack to the PC via USB. To prevent Windows from mixing the microphone with the playback device mute the microphone in the playback mixer of the Windows volume control.

3.3. Bluetooth® connection

3.3.1. Preparing for Bluetooth® connection

With the ALVA BC640 a class 1 high range USB Bluetooth® dongle for Windows XP and Windows Vista is included. Depending on the Bluetooth® software included with other Bluetooth® modules, it will be possible that not all ALVA BC640 functionality will be supported. The Bluetooth® dongle and software, which is included with the ALVA BC640, is optimized for best performance and Bluetooth® range. We therefore highly recommend using the supplied dongle for best results and ease of installation.

3.3.2. Setting up the Bluetooth® USB dongle

Once the screenreader drivers have been installed, the next step is the installation of the supplied Bluetooth® USB dongle.

1. Upon Windows start-up, insert the ALVA BC640 software CD-ROM. The installation program will automatically run if the autostart feature of your CD-ROM drive is enabled. Please make sure to use the ALVA BC640 CD and not the Belkin CD.
2. Select Bluetooth® Setup
3. The setup program will guide you through the installation.

Warning!

Please note that the installation of the Bluetooth® software might take more than 20 minutes. Please do not interrupt the setup, even when the setup appears to be unresponsive. Interrupting the setup will result in an incorrect installation of the Bluetooth® software and improper functioning of especially the ALVA BC640 Audio and keyboard. In case the installation has been interrupted, please uninstall the Bluetooth® drivers before attempting to reinstall the Belkin Bluetooth® drivers again. The uninstallation procedure may take more than 20 minutes. Please make sure to not interrupt the uninstallation procedure.

After installation of the Bluetooth® dongle and its software an icon named “My Bluetooth Places” will appear on the desktop of your Windows computer. To create the Bluetooth® connection with your ALVA BC640 base unit or ALVA BC640 with Braille Audio Feature Pack, the device will need to be paired with the PC. Once the ALVA BC640 has been paired successfully, the required Bluetooth® services have to be activated.

The ALVA BC640 without Braille Audio Feature Pack has one Bluetooth® module and uses one Bluetooth® service. The ALVA BC640 with Braille Audio Feature Pack has two Bluetooth® modules and uses four Bluetooth® services.

Bluetooth® modules and services:

	Bluetooth® Services			
Model	Braille Bluetooth ® Serial Port	Bluetooth® HID Keyboard	Bluetooth® Mono Audio	Bluetooth® Stereo Audio
ALVA BC640 without Braille Audio FP	ALVA BC640 Bluetooth® module			
ALVA BC640 with Braille Audio FP	ALVA BC640 Bluetooth® module	ALVA BC640 Bluetooth® module	ALVA BC640 Bluetooth® module	ALVA BC640 FP Bluetooth module

As indicated in the figure above, the ALVA BC640 Base Unit only uses the Bluetooth® serial service. However, the ALVA BC640 with Braille Audio Feature Pack uses two modules; one module for the serial connection or serial service, the HID (Human Interface Device) keyboard connection or keyboard service and the mono headset connection or mono headset service. The second Bluetooth® module is used for the stereo audio.

Pairing the ALVA BC640 Base Unit (without Braille Audio Feature Pack):

1. Activate the option “My Bluetooth Places”
2. A Wizard will automatically open in Windows Explorer.
3. Activate the “Find Bluetooth devices” option (we recommend you not to use the Bluetooth® setup wizard).

4. Locate the ALVA BC640 in the list with found devices, by using the arrow keys or click it once with the mouse. If the ALVA BC640 is not found, press the F5 key to search again.
5. Press ALT + B for the Bluetooth® menu and select "Pair device".
6. To pair the ALVA BC640 with your Bluetooth® dongle, you have to enter a PIN code. The default code is 0000. In case you changed this number in the Local Menu of the ALVA BC640, please enter your personal code and press Enter. The ALVA BC640 is now paired with your computer.
7. Select the ALVA BC640 from the list with found devices again and press Enter.
8. Select "Dev B on ALVA BC640" to connect the serial port. A dialog box with the serial port number will appear. Press Enter to close this dialog box. The ALVA BC640 will issue an audible signal and displays the message "BT connected" to indicate that the connection has been established and is functioning.

Pairing the ALVA BC640 with Braille Audio Feature Pack:

1. Activate the option "My Bluetooth Places"
2. A Wizard will automatically open in Windows Explorer.
3. Activate the "Find Bluetooth devices" option (we recommend you not to use the Bluetooth® setup wizard).
4. Locate the ALVA BC640 in the list with found devices, by using the arrow keys or click it once with the mouse. If the ALVA BC640 is not found, press the F5 key to search again.
5. Press ALT + B for the Bluetooth® menu and select "Pair device".
6. To pair the ALVA BC640 with your Bluetooth® dongle, you have to enter a PIN code. The default code is 0000. In case you changed this number in the Local Menu of the ALVA BC640, please enter your personal code and press Enter. The ALVA BC640 is now paired with your computer.
7. Select the ALVA BC640 from the list with found devices again and press Enter.
8. Select "Dev B on ALVA BC640" and press Enter to connect the serial port. A dialog box with the serial port number will appear. Press Enter to close this dialog box. The ALVA BC640 will issue an audible signal and displays the message "BT connected" to indicate that the connection has been established and is functioning.
9. Select "Headset on ALVA BC640" and press Enter to connect the mono audio. A dialog box confirming that the audio connection is established will appear. Press Enter to close this dialog box.
10. Select "Human Interface Device on ALVA BC640" and press enter to connect the Braille keyboard.
11. Press Backspace to return to the "Find Bluetooth Devices" screen.

12. To pair the second module of the Braille Audio Feature Pack for stereo audio, locate the ALVA BC640 Braille Audio Feature Pack in the list with found devices, using the arrow keys or clicking the mouse once. If the ALVA BC640 Braille Audio Feature Pack is not found, press the F5 key to search again.
13. Press ALT + B for the Bluetooth® menu and select “Pair device”
14. To pair the ALVA BC640 Braille Audio Feature Pack with the Bluetooth® dongle, you have to enter a PIN code. The default code is 0000. In case you changed this number in the Local Menu, please enter your personal code and press Enter. The ALVA BC640 with Braille Audio Feature Pack is now paired with your computer.
15. Select the ALVA BC640 Braille Audio Feature Pack from the list with found devices again and press Enter. A dialog box confirming that the stereo audio connection is established will appear. Press Enter to close this dialog box.

When the ALVA BC640 screenreader drivers are properly installed (see paragraph 3.2), you might have to reload your screenreader to establish the Bluetooth® connection.

4. The ALVA BC640 Local menu

This chapter describes the ALVA BC640 Local Menu. The ALVA BC640 Local Menu enables you to change several built-in settings of the Braille display. In case you are using an ALVA BC640 with Braille Audio Feature Pack, there will be two extra menus: the Bluetooth® connection menu and audio menu.

When the ALVA BC640 is switched on, the following message is displayed:

“ALVA BC640 ready” followed by the date and time.

This message will automatically disappear when the screenreader is loaded and a successful connection with the ALVA BC640 is established.

4.1. Activating the Local Menu

To enter the ALVA BC640 Local Menu press the keys Smartpad key 2 + Smartpad key 3 + arrow down key simultaneously. To exit the Local Menu press Smartpad key 1. The Smartpad keys are located on top of the ALVA BC640. The arrow keys and Enter key are located in the centre of the Smartpad. The four remaining keys are labelled from left to right Smartpad key 1 to 4. Refer to chapter 2 for a detailed description of the ALVA BC640 keys.

4.2. Changing settings in the Local Menu

The Local Menu allows you to configure the ALVA BC640 or review information about the current status of the ALVA BC640. Some menu items will only be active when a Braille Audio Feature Pack has been installed. These menu items will proceed with “FP”.

Local menu	Key Combination
Next menu item	Smartpad arrow up key
Previous menu item	Smartpad arrow down key
Select menu item	Smartpad Enter key
Next field in menu item	Smartpad arrow right key
Previous field in menu item	Smartpad arrow left key
Change value in menu field item	Smartpad arrow up / down keys
Accept changed value in menu field item and leave field	Smartpad Enter key
Exit menu without saving changes	Smartpad key 1

Please note that when the menu is activated, it will be underlined with dot 7 and 8 on the Braille display.

The Local menu has the following items:

Battery: 45%, charging

Battery presents the status of the battery and informs you if the battery is charging. In case the USB connection is used, the ALVA BC640 will start charging via USB automatically.

The status of the battery loading can be:

- Charging The battery is charging
- Not charging The battery is not charging

FP Notes

Notes. Press Smartpad Enter to activate the internal Notes function.

Date/Time: 01/01/2008 01:00:00

Date and time. This menu item shows the current date and time. To set the date and time and to change the date and time format, please refer to the set date, set time and date and time format menu options.

Bluetooth status: connected

Bluetooth® status displays if the ALVA BC640 Bluetooth® is connected or disconnected.

Set Bluetooth pairable

Set Bluetooth® pairable. Activating this menu item will disconnect the current Bluetooth® connection and make the ALVA BC640 visible to other Bluetooth® devices. Alternatively, switching the ALVA BC640 off and on will disconnect the current Bluetooth® connection and switches the Bluetooth® module into the pairable state too.

Bluetooth: on

Bluetooth®. The ALVA BC640 Bluetooth® module can be switched off. To turn Bluetooth® off, press Smartpad Enter key and select “off” with the Smartpad up and down keys. To switch Bluetooth® on again, select “on”.

Bluetooth pin: 0000

Bluetooth® PIN. The ALVA BC640 is protected by a Personal Identification Number (PIN) to prevent other people from setting up a wireless connection with your Braille

display. This setting allows you to enter your own PIN code. By default the PIN code is set to 0000. This number is required when setting up a wireless connection with your PC or other Bluetooth® device. You can change the digits in the PIN by using the Smartpad Enter key and the Smartpad arrow up and down, left and right keys.

Bluetooth name: ALVA BC640

Bluetooth® Name. Specify the Bluetooth® name for the ALVA BC640. This allows you to uniquely identify the ALVA BC640 in case you are setting up a wireless Bluetooth® connection with your PC or other Bluetooth® device. By default the Bluetooth® name is ALVA BC640. You can change the separate characters by scanning through the alphabet. To change the Bluetooth® name press the Smartpad Enter key. Use the Smartpad arrow up and down keys to scan forward and backward through the alphabet. Use the Smartpad arrow left and right keys to move to the previous or next character. The maximum length of the Bluetooth® name is 16 characters. In case a Braille Audio Feature Pack has been installed, the Bluetooth® name of the Feature Pack will have the extension “FP”. By default the name of the Braille Audio Feature Pack is ALVA BC640 FP.

Sleep: 30

Sleep. By default the sleep time function is set to 30 seconds. Press the Smartpad Enter key and set a time from 5 seconds to 255 seconds, indicating the time after which the Braille display puts itself into sleep mode if there is no activity. To disable the sleep function, select never. Consequently, the ALVA BC640 consumes virtually no energy when the Braille display is in sleep mode. Set the time by using the Smartpad arrow up and down keys. Increase the value by pressing the Smartpad arrow up key. Decrease the value by pressing the Smartpad arrow down key.

FP Braille configuration

Braille configuration. Press the Smartpad Enter key to open the Braille configuration submenu. The Braille configuration submenu contains three menu items: keyboard layout, Braille table and contracted (Grade 2) table. These three items are being used to translate the Braille input on the 8 key Braille keyboard into text. It is therefore important to select the same keyboard layout on the ALVA BC640 as used by Windows. Your Windows keyboard layout can be found in the *Windows Control Panel/ Regional and Language options/ Languages/ Details*. When a keyboard layout is selected, the ALVA BC640 will attempt to match the Braille table and contracted table if present. If there is no matching Braille table present, no changes will be made. Alternatively, you can manually select a different Braille table or Contracted table. The default Braille table is the NABCC (North American Braille Character Code) and the US English contracted table. The selected Braille table is used in the Local Menu. Please note that screenreaders have their own Braille tables. Therefore, the Braille that is being displayed on the Braille display cannot be controlled by the ALVA BC640. Selecting a different Braille table for displaying Braille should be done in the screenreader itself.

To enable or disable contracted (or Grade 2) Braille input, press Space + B1245 (Braille keys 1, 2, 4, 5). The ALVA BC640 will generate a series of beeps to indicate that contracted Braille is turned on or off. Pressing Smartpad 1 allows you to exit this submenu.

Dot pressure: 5

Dot pressure. Select the Braille reading comfort by setting the dot pressure. If you prefer hard and stiff Braille dots set the value to 5. If you prefer soft and weak Braille dots, set the value to 1. The default setting is 5.

Number status cells: 0

Number status cells. By default the number of status cells is 0. To enable status cells, select a value of 1 to 20 status cells. If the statuscells are set to 0, the statuscells will be disabled and all 40 Braille cells are used for displaying text and data. Status cells provide useful information about the text on the screen, Windows controls and screenreader mode.

Status cells: Left

Status cells. The default position of the status cells is on the left. Select the position of the status cells either at the left or right side on the Braille display.

Split point active: off

Split point active. When the ALVA BC640 split point **is enabled**, the Braille display will be divided into two segments. This enables you to read and interact with two devices simultaneously. For example, the Bluetooth® screenreader output of a mobile phone can be displayed on the first 10 Braille cells while the remaining 29 cells can be used to read the PC screenreader output via the USB connection. The remaining cell of the 40 cells will have all 8 dots raised to indicate where the first segment ends and the second segment starts. The first segment with a default length of 10 Braille cells will display the information send via the Bluetooth® connection. The second segment with a default length of 29 cells will display the information send via the USB connection. To enable or disable the split point feature select on or off using the Smartpad arrow up and down keys. Alternatively, you might enable or disable the split point feature by pressing the hotkeys Smartpad 2 + 3 + Enter. To be able to use the split point feature, please make sure your screenreaders supports this function. To find out if your screenreader supports the ALVA BC640 split point functionality, please refer to your screenreader documentation.

Split point: 10

Split point. Adjust the number of cells used for the split point function by selecting the number of cells with the Smartpad arrow up and down keys. The first segment can be set from 5 to 35 cells. By changing the length of the first segment, the length of the second segment will automatically be adjusted.

Repeat keys: on

Repeat keys. The repeat function of the thumb keys can be enabled or disabled. When holding down a thumb key with the repeat function switched on, the thumb keys will start repeating after half a second. This enables you to quickly scroll the display by keeping the thumb key pressed.

2nd cursor routing row: on

2nd cursor routing row. The ALVA BC640 is equipped with one row of cursor routing keys. However, by keeping the cursor routing key pressed the ALVA BC640 is able to emulate a second row of cursor routing keys. This will give you the possibility to perform a left and a right mouse click with one row of cursor routing keys. When enabled, the ALVA BC640 will differentiate between a short and a long cursor routing key press. A short press will be treated as a press of a key of the first row of cursor routing keys, while a long press will be seen as a press of a key of the second row of cursor routing keys. The length of a “long press” can be set by adjusting the delay time. The default time is half a second. The assignment of the first and second row of cursor key actions depends on the used screenreader. Refer to your screenreader documentation for the assignment of the first and second cursor routing key actions.

2nd cursor routing delay: 500 ms

2nd cursor routing delay. The length of the delay before the second cursor routing key is activated can be set by adjusting the delay time. By default the delay is set to 500 ms or half a second.

FP Default volume: 31

Default volume is the volume of the Braille Audio Feature Pack when the ALVA BC640 is switched on. The default volume is 31.

Sounds: off

Sounds can be switched on or off. When switched on, a short beep is being generated on each key press. If you prefer no sounds to be generated, select “off”.

For each group of keys the beeps can be switched on or off individually.

Cursor routing sound: on

Cursor routing sound can be switched on or off. When switched on, a short beep is generated on each cursor routing key press. If you prefer no sounds to be generated when a cursor routing key is pressed, select off.

eTouch key sound: on

eTouch key sound can be switched on or off. When switched on, a short beep is generated on each eTouch key press. If you prefer no sounds to be generated when an eTouch key is pressed, select off.

Thumb key sound: on

Thumb key sound can be switched on or off. When switched on, a short beep is generated on each thumbkey press. If you prefer no sounds to be generated when a Thumb key is pressed, select off.

Smartpad sound: on

Smartpad sound can be switched on or off. When switched on, a short beep is generated on each Smartpad key press. If you prefer no sounds to be generated when a Smartpad key is pressed, select off.

Set time: 01:01:00

Set time. This setting allows you to set the time. Press the Smartpad Enter key to enter the modification mode. Use the Smartpad arrow up and down keys to change the time. Use the Smartpad arrow left and right keys to select the hour, minute or seconds.

Set date: 01:01:2008

Set date. This setting allows you to set the date. Press the Smartpad Enter key to enter the modification mode. Use the Smartpad arrow up and down keys to change the date. Use Smartpad arrow left and right keys to select the day, month and year.

Date format: dd-mm-yyyy

Date format. This setting allows you to set the date format. By default this is set to day / month / year. Press the Smartpad Enter key to enter the modification mode. Use the Smartpad arrow up and down keys to change the date format. The date can be displayed in the formats dd-mm-yyyy, yyyy-mm-dd and mm-dd-yyyy.

Date separator: /

Date separator This setting allows you to set the date separator. By default the “/” symbol is used. Press the Smartpad Enter key to enter the modification mode. Use the Smartpad arrow up and down keys to change the date Separator symbol. The date separator can be changed in a “/”, “.”, “.” or “-”.

Time format: 24 hours

Time format. This setting allows you to set the time format. The default is the 24 hour format. Press the Smartpad Enter key to enter the modification mode. Use the Smartpad arrow up and down keys to change the time format. The time can be displayed in either the 24 hours or 12 hours format.

Time separator: :

Time separator. This setting allows you to set the time separator. The default is the “.” symbol. Press the Smartpad Enter key to enter the modification mode. Use the Smartpad arrow up and down keys to change the time Separator symbol. The time separator can be changed in “.” or “.”.

Language: **English**

Language. Select the language for the Local Menu.

Preferred interface: Bluetooth

Preferred interface. Select the preferred interface; Bluetooth® or USB. If the ALVA BC640 is connected to two mainstream devices using its Bluetooth® and USB connection, the ALVA BC640 will only “communicate” with the preferred interface. By default the preferred interface is Bluetooth®. In the default situation, when the ALVA BC640 is connected to two devices via Bluetooth® and USB, the ALVA BC640 will only “communicate” with the Bluetooth® connection. The USB connection will be ignored. The hotkey to change the preferred interface is Smartpad 2 + Smartpad 3 + Cursor Routing Key 1 for the Bluetooth® connection and Smartpad 2 + Smartpad 3 + Cursor Routing Key 40 for the USB connection.

BT MAC: xxxx:xx:xxxxxx

BT MAC refers to the ALVA BC640 Bluetooth® MAC address. A Bluetooth® MAC address uniquely identifies a Bluetooth® device.

Firmware: 1.10.64 BT: 1.11.0

Firmware version presents the version of the built in ALVA BC640 software and the software version for the Bluetooth® module. If a Braille Audio Feature Pack has been installed the second number after “BT” represent the firmware version for the Braille Audio Feature Pack. Please make sure that the firmware versions are equal to or higher than 1.14.0 BT: 1.16.0 1.4.0 when using the ALVA BC640 with Braille Audio Feature Pack. Please refer to chapter 8 for detailed instructions on how to update the firmware.

Hardware: 18

Hardware version presents the version of the ALVA BC640 hardware revision.

Firmware update

Firmware update. When activated with the Smartpad Enter key, the ALVA BC640 will be in “Firmware Update” mode. Please refer to chapter 8 for detailed instruction on how to update the firmware.

Restore Factory settings

Restore factory settings allows you to restore the factory settings. If you select this option the ALVA BC640 factory defaults will be restored. Alternatively, the ALVA BC640 factory defaults can be restored by deleting the SETTINGS.A6 file from the internal ALVA_BC6 drive. To delete the SETTINGS.A6 file:

- Connect your ALVA BC640 to the USB port of your PC using the supplied USB cable
- The ALVA_BC6 driver will automatically appear
- Locate and open the system folder
- Locate the SETTINGS.A6 file
- Delete the file

Note: The system folder is a hidden folder. To see hidden folders, you have to enable “show hidden files and folders” in Windows Explorer.

Test mode

Test mode Use the test submenu if you have the impression that certain keys or Braille cells are not responding properly. The ALVA BC640 has two built-in tests: a Braille cell test and a key test. When activating the key test, pressing a particular key provides a unique dot pattern on the array of Braille cells. The Braille cell test will cycle all dots continuously. Pressing Smartpad key 1 to exit the test mode.

4.3. The Feature Pack Menus

The ALVA BC640 with Braille Audio Feature Pack has two extra menus: the Bluetooth® connection menu and Audio menu.

4.3.1. The Braille Audio Feature Pack Bluetooth® Menu

To enter the ALVA BC640 Bluetooth® Menu press Smartpad key 2 + Smartpad key 3 + arrow left simultaneously. To exit the Bluetooth® Menu press Smartpad key 1. The Smartpad keys are located on the top of the unit below the Braille line. The arrow keys and Enter key are located in the center of the Smartpad. The four remaining keys are labelled from left to right: Smartpad key 1 to 4. Please refer to chapter 2 for a detailed description of the keys of the ALVA BC640 Base Unit and Braille Audio Feature Pack.

BT: Braille keyboard stereo mono

Bluetooth® Menu. The Bluetooth® Menu shows the status of the various Bluetooth® connections. If a connection is activated dot 7 and 8 will be raised to indicate the connected state. Clicking on a connection with the cursor routing keys will connect or disconnect the particular Bluetooth® connection. Note that the Windows drivers will try to reconnect the Braille connection. Therefore it is highly recommended to disconnect the Braille in the screenreader and not from the ALVA BC640. Secondly, dismantling the keyboard from the ALVA BC640 will prevent the ALVA BC640 driver to reconnect the keyboard until it has been connected from the PC again. In general, it is recommended to have the ALVA BC640 software taking care of the connections. Manually connecting and disconnecting any Windows Bluetooth® connection from the ALVA BC640, except the stereo connection, may result in unexpected behaviour.

The first time the ALVA BC640 is trying to reconnect the audio and keyboard, the Windows Bluetooth® software will generate a sound and show a popup window, asking for permission to have the ALVA BC640 automatically reconnect the audio and keyboard. To make sure the audio and keyboard automatically reconnect when a connection is lost or during a restart of Windows or the screenreader, please select the checkbox to always allow the connection and click the OK button. Please follow this procedure for the keyboard connection, mono connection and stereo connection. Should the popup not appear at the Bluetooth® notification sound, click on the Bluetooth® icon in the system tray next to the clock.

4.3.2. The Braille Audio Feature Pack Audio Menu

To enter the ALVA BC640 Audio Menu press Smartpad key 2 + Smartpad key 3 + the arrow right key simultaneously. To exit the Audio Menu press Smartpad key 1. The Smartpad keys are located on the top of the unit below the Braille line and cursor routing keys. The arrow keys and Enter key are located in the center of the Smartpad. The four remaining keys are labelled from left to right: Smartpad key 1 to 4. Please refer to chapter 2 for a detailed description of the keys of the ALVA BC640 Base Unit and Braille Audio Feature Pack.

Audio: BT USB Off Exit

Audio Menu. The Audio Menu shows the availability and status of the Bluetooth® and USB audio. And audio channel is active and available when shown in the menu and underlined with dot 7 and 8. To switch to another audio channel, click the cursor routing key on the audio channel. Clicking the cursor routing key on the menu item “off” wil disable the audio, whereas clicking on “Exit” will close the menu.

5. Power Sources

There are three ways to power the ALVA BC640:

- Using the USB power (Base Unit only);
- Using the internal battery;
- Using the external power adapter

In the subsequent sections, you will find more information about each power option.

5.1. *USB power*

If you connect your ALVA BC640 directly to a PC using the USB connection, the PC will provide power to the Braille display through the USB cable. To power the ALVA BC640 with Braille Audio Feature Pack, please use the power supply that came with the unit. Please note that only the ALVA BC640 Base Unit can be powered via the USB connection. The USB connection provides sufficient power to charge the base unit's internal battery and to operate the ALVA Braille Controller Base Unit.

5.2. *Internal battery*

The ALVA BC640 runs on an internal battery in case no external power is attached or no USB connection is used. The ALVA BC640 will approximately run 12 hours on a fully charged battery. Typically, the internal battery is used in case of a wireless Bluetooth® connection.

5.3. *External power adapter*

An external power adapter is provided with the ALVA BC640 package. As soon as the external power adapter is connected to the Braille Controller, the battery will be charged and the ALVA BC640 will draw its power from this adapter. Please note that the ALVA BC640 with Braille Audio Feature Pack cannot be charged via the USB connection. Charging an empty battery will take approximately 4 hours. The external power can be used in combination with wireless Bluetooth® technology or the USB connection.

Warning! Do not use another external power adapter than the one supplied in the ALVA BC640 package.

6. The ALVA BC6 Bluetooth® settings application

This chapter describes the use of the ALVA BC6 Bluetooth® settings application. This application can be found in the *Start Menu/ Programs/ Optelec/ ALVA_BC6* folder. When having paired the ALVA BC640 with Braille Audio Feature Pack with the PC and having installed the screenreader drivers and Bluetooth® software, this application enables configuration of the connection settings of the Braille keyboard and mono and stereo Bluetooth® audio.

The ALVA BC6 Bluetooth® settings application enables configuration of the automatic connection of the keyboard and audio when activating the screenreader. In case the connection gets lost, the application automatically instructs the ALVA BC640 driver to reconnect the keyboard and audio.

In order to enable this functionality, you will have to create the keyboard and audio connections once, using “My Bluetooth® Places”. The first time the ALVA BC640 driver creates the audio and keyboard connection, the Windows Bluetooth® will generate a sound and a popup window asking for permission to have the ALVA BC640 automatically reconnect the audio and keyboard.

To make sure the audio and keyboard automatically reconnect when a connection is lost or during a restart of Windows or the screenreader, please select the checkbox to always allow the connection and click the OK button. Please follow this procedure for the keyboard connection, mono connection and stereo connection. Should the popup not appear at the Bluetooth® notification sound, click on the Bluetooth® icon in the system tray next to the clock.

6.1. *Configuring Audio*

There are two settings:

- Automatically connect Bluetooth® Braille keyboard;
- Manually connect Bluetooth® Braille keyboard.

If the first setting has been selected, the ALVA BC640 will automatically connect the Braille keyboard when activating the screenreader. In case the connection gets lost, the keyboard will attempt to reconnect automatically.

6.2. *Using the Braille keyboard*

The Bluetooth® audio settings are only available in Windows XP. The audio settings are greyed out in Windows Vista. The security restrictions of Windows Vista prevent the ALVA BC640 software from configuring the audio connections. It is therefore recommended to use the stereo connection in Windows Vista, as the stereo connection is more reliable and provides a better audio quality.

There are three settings:

- Automatically connect mono Bluetooth® audio;
- Automatically connect stereo Bluetooth® audio;
- Manually connect Bluetooth® audio

In order to have the sound correctly routed to the ALVA BC640 in Windows XP, please select the audio playback devices for the default audio device, the Bluetooth® mono playback device and the Bluetooth® stereo playback device. Please note that these options are not available in Windows Vista.

The name of the default audio device depends on the brand and model of the sound card manufacturer. The mono Bluetooth® device is in most situations labelled as Bluetooth-audio. To select the stereo Bluetooth® audio device first, you will need to create a Bluetooth® stereo connection:

1. Open "My Bluetooth Locations";
2. Go to "Find Bluetooth Device";
3. Search for ALVA BC640 FP and make sure it is paired;
4. Press Enter on ALVA BC640 FP to create the stereo connection;
5. (Re)start the ALVA BC6 Bluetooth® settings application;
6. Select high quality Bluetooth® audio to set the stereo Bluetooth playback to audio device.

To save the settings to all users, make sure the "Save to default user" is checked and press Enter on the "OK" button. To save the settings to the current user, make sure the checkbox is not selected.

7. Using the Braille keyboard

This chapter gives an overview of the Braille key combinations to emulate keyboard function keys. There are two types of function keys; chord commands to emulate keyboard function keys and prefix chord commands to emulate function key combinations.

7.1 Chord commands

Chord commands are Braille key commands in combination with the space bar. These chord commands are used to emulate keyboard function keys such as Home, End, Backspace, Delete and Escape. Please refer to the figure below for an overview of the chord commands:

Chord command	Keyboard key
Backspace	Space + b12
Escape	Space + b15
Delete	Space + b45
Page up	Space + b13
Page down	Space + b46
Home	Space + b123
End	Space + b456
Tab	Space + b56
Shift-Tab	Space + b23
Shift-Space	Space + b7

7.2 Prefix chord commands

There are three prefix chord commands to emulate the function keys F1 to F12 (Space + b124 for activating F, followed by the number 1 to 12), Insert (Space + b24) and Caps Lock (Space + b14) combinations.

To press one of the function keys F1 to F12, press Space + b124 (Space + F) followed by a character. The function key number is represented by the characters “a” to “l” whereby the character “a” represents F1 and “l” represents F12. To press F1 using the ALVA BC640 Braille keyboard, first type space + b124 followed by the character “a” (Braille key 1). Press the Space + b124 combination again to exit this function.

Similarly, the Insert and Caps Lock can be typed by pressing Space + b24 for Insert and Space + b14 for Caps Lock. To type an Insert + J command, press Space + b24 followed by a “j” (Braille keys 2, 4, 5). Typing the Insert or Caps Lock command twice will switch the status of this command on or off, depending on the current status.

7.3 *Contracted Braille on/off*

The chord + g command or Space + b1245 will turn contracted Braille (Grade 2) on or off. The contracted Braille table can be best set in the *Local Menu/ Braille configuration/ Contracted Braille table*.

8. Notes

The Notes application enables writing and saving notes on the internal memory of the ALVA BC640 with Braille Audio Feature Pack. Notes contains a few elementary text adding and editing functions, such as copy, cut, paste and a block function.

To use Notes no PC connection is required. In case the ALVA BC640 is connected to a PC, you can continue working with the PC after closing Notes. The ALVA BC640 Braille keyboard has to be used to enter Notes. It is not possible to enter text in Notes via the PC keyboard.

Text entered in Notes will be saved in a proprietary Braille format. The typed Braille characters will be saved in the same format as they have been entered. This means the characters are independent of any Braille code. Please keep in mind that Notes is not compatible with Microsoft Word or other similar text application. Therefore, it will not be possible to open Word documents in Notes or Notes documents in Microsoft Word. Text that has been entered in Notes can only be opened, read and edited in Notes.

8.1 Starting Notes

To start the Notes application:

1. Open the Local Menu by pressing Smartpad key 2 + Smartpad key 3 + the arrow down key simultaneously. Use the arrow down key to navigate through Notes and press the Smartpad Enter key to start Notes.
2. Press Smartpad key 2 + Smartpad key 3 + the arrow up key simultaneously to start Notes directly.

In both situations Notes will open a new page.

8.2 Using Notes

This paragraph will describe the various edition options of Notes in more detail.

8.2.1 Adding text

Notes can be entered using the ALVA BC640 Braille keyboard. An Enter will be displayed as an empty cell followed by a cell with all 8 dots raised followed by another empty cell.

8.2.2 Editing text

Use the following chord commands to edit text.

Chord command	Keyboard Key
Backspace	Space + b12
Delete	Space + b45

8.2.3 Navigating in Notes

Please use the following commands to navigate through the text:

Chord command	Keyboard Key
Previous character	Space + b3
Next character	Space + b6
Previous line	Space + b1
Next line	Space + b4
Previous word	Space + b2
Next word	Space + b5

Alternatively, the cursor routing keys can be used to move the cursor and navigate through the text.

8.2.4 Defining a block

Please follow the procedure below to define a text block:

1. Place the cursor at the beginning of the defined block;
2. Press Ctrl+J to mark the beginning of the block;
3. Move the cursor to the end of the block;
4. Press Ctrl+K to mark the end of the block.

The defined block can be copied or cut and pasted to a preferred position in the text. The short keys for cut, copy and paste are, Ctrl+X, Ctrl+C and Ctrl+V. Alternatively, the commands can be found in the Notes menu.

8.2.5 Reading in Notes

The following commands have been defined for reading text in Notes:

Command	ALVA BC640 Key
Previous word	Thumb key 1
Start of line	Thumb key 2
Start of next line	Thumb key 3
Next word	Thumb key 4

Pressing thumb key 2 while the cursor is within the line will move the cursor to the beginning of the current line. Pressing thumb key 2 again will move the cursor to the beginning of the previous line.

In case the entered text in Notes does not contain any new line commands (Enters) thumb key 2 and 4 will move the respectively the beginning and end of the text.

8.3 *The Notes Menu*

The Notes menu can be activated by pressing the Alt-key. Click on the cursor routing key to activate a menu item. To exit the menu, click the Alt-key again. Notes menu options can be activated with short keys. Please refer to paragraph 8.4 for a list with short key combinations.

8.3.1 Notes main menu

To enter the main menu, press the Alt-key. The Notes main menu consists of the following three submenus:

- File
- Edit
- Help

Press the Alt-key again to exit the main menu.

8.3.1.1 File menu

To enter the File menu, click on File in the main menu, using the cursor routing keys. The File menu contains the following options:

- New
- Open
- Save
- Save As
- Exit

New Note

To create a new note, click on New in the File menu. In case the current note contains unsaved text, Notes will ask you to save the changes.

Open a Note

To open an existing note, click on Open in the File menu, using the cursor routing keys. A list with existing notes will be displayed. Navigate through the list using the Smartpad arrow up and down keys. Press the Smartpad Enter key to open the selected note.

Save a Note

To save a note, click on Save in the File menu. If the note has no name yet, you will be prompted to enter a new name.

Save As

Clicking on Save As will prompt you to give your current note a name. If the name already exists you will be asked to overwrite the existing file. Clicking on "No" will prompt you to enter a different file name.

Exit Notes

Clicking on Exit will close Notes. If the current note is not saved, you will be prompted to save the file. While clicking “Yes” you will be asked to enter a file name. Clicking “No” will close Notes without saving the file.

8.3.1.2 Edit menu

The Edit menu contains the following options:

- Undo
- Cut
- Copy
- Paste

Undo

To undo the last change, click on “undo” in the Edit menu.

Cut

By clicking “cut”, the selected block or character will be cut. The selected block or character can be pasted on a different location in the text. The shortcut for “cut” is Ctrl+X.

Copy

By clicking “copy”, the selected block or character will be copied. The selected block or character can then be pasted on a different position in the text. The short key combination for “copy” is Ctrl+C.

Paste

To paste a “cut” or “copied” block, click “paste”. Alternatively, use the short key combination Ctrl+V.

8.3.1.3 Help

The Help menu displays the about screen. Press Space+b15 to exit this screen.

8.4 Notes Short Key Combinations

Below is an overview of the most frequently used short key combinations:

Command	Keyboard Key
Previous line	Space + b1
Previous word	Space + b2
Previous character	Space + b3
Next line	Space + b4
Next word	Space + b5
Next character	Space + b6
Backspace	Space + b12
Delete	Space + b45
Escape	Space + b15
Copy	Space + b14
Mark start of block	Space + b245

Mark end of block	Space + b13
New note	Ctrl + b1345
Open note	Ctrl + b135
Save note	Ctrl + b234
Paste	Ctrl + b1236
Open / Close Notes Menu	Alt

9. Updating the Firmware

The ALVA BC640 firmware can be updated using the supplied firmware update software. The ALVA BC640 Firmware Updater is located in the *Start Menu/ Programs/ Optelec/ ALVA BC640 Updater*. Please read this chapter carefully before attempting to update your ALVA BC640 firmware.

Caution:

- *An incorrect update of the ALVA BC640 firmware will seriously damage the ALVA BC640, which may result in having to return the unit to the factory for repair.*
- *An incorrect applied firmware update is not covered by warranty.*
- *Please check if the battery is charged and if the ALVA BC640 is connected to the power supply before applying the firmware update.*
- *Please remember to never switch off the ALVA BC640 while the firmware update is in progress.*


To start the firmware update procedure, please follow the following steps:

1. Switch the ALVA BC640 on;
2. Connect the ALVA BC640 to the PC using the USB cable that came with the unit;
3. Activate the ALVA BC640 Local Menu by pressing Smartpad key 2 + Smartpad key 3 + arrow down simultaneously;
4. Start the ALVA BC640 Firmware Updater, which is located in the *Start Menu/ Programs/ Optelec/ ALVA BC640 Firmware Updater* and follow the on screen instructions carefully;
5. When the ALVA BC640 Firmware Updater asks you to put the ALVA BC640 in firmware update mode, locate the "Firmware Update" in the Local Menu, by pressing Smartpad key 2 + Smartpad key 3 + arrow down simultaneously. Navigate through the Local Menu, using the arrow up key and press the Smartpad Enter key to select "Firmware Update". If this is the first time you update the firmware, a hardware wizard will appear asking to install the necessary USB drivers. Please follow the instructions to install the USB drivers.
6. Please wait until the firmware update is finished and the message "Bluetooth base module updated" is being displayed on the ALVA BC640. This might take several minutes.


10. Mounting the Braille Audio Feature Pack

The ALVA BC640 Braille Audio Feature Pack clicks onto the back of the ALVA BC640 after which it will be secured with two 2.5 mm hex screws (see picture below).


1. Switch off the ALVA BC640



2. Click the Feature Pack to the ALVA BC640



3. Tighten the 2 Feature Pack screws with the 2.5 mm hex key



Regulatory notice

This equipment has been developed and tested in order to comply with the following regulations:



and



for a Class B digital device, pursuant to part 15 rules.

These limits are designed to provide reasonable protection against harmful interference in a residual installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Re-orient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the distributor or an experienced technician for help.

Medical electrical equipment needs special precautions regarding EMC (Electro Magnetic Compatibility) and needs to be installed and put into service according to the EMC guidelines.

WARNING: The use of accessories, transducers and cables other than those specified and sold by the manufacturer of the ALVA BC640 as replacement parts for components, may result in increased emissions or decreased immunity of the ALVA BC640.

By design the ALVA BC640 will shut down in case of an ESD 8KV test pulse. The ALVA BC640 will need to be switched on again to resume.

CE Caution:

- ◆ That SIP/SOPs need to be connected properly and that any unused SIP/SOPs shall not be accessible for the users to contact to other nonrecommand devices in the user area .
- ◆ The unit is for exclusive interconnection with IEC/EN 60601 certified equipment outside of user environment and IEC/EN 60601-1 certified equipment inside the user environment.
- ◆ This device comply with EN 60601-1-2, to minimize the interference from other equipment.
- ◆ Equipment connected to the analog or digital interfaces must comply with the respective IEC/EN standards (e.g. IEC/EN 60950-1 for data processing equipment and IEC/EN 60601-1 for medical equipment).
- ◆ To reduce the risk of electric shock, make sure power cord is unplugged from wall socket. To fully disengage the power to the unit, please disconnect the power cord from the ac outlet. Do not remove cover(or back). No user

serviceable parts inside. Refer servicing to qualified service personnel. The AC outlet shall be readily available and accessible.

- ◆ Furthermore all configurations shall comply with the current version of the standard for SYSTEMS IEC/EN 60601-1-1. Everybody who connects additional equipment to the signal input part or signal output part configure a medical system, and is therefore responsible that the system complies with current version of the requirements of the system standard IEC/EN 60601-1-1. If in doubt, consult the technical service department or your local representative."

Explanation Of Graphical Symbols:



This symbol warns user that uninsulated voltage within the unit may have sufficient magnitude to cause electric shock. Therefore, it is dangerous to make any kind of contact with any part inside this unit.



This symbol alerts the user that important literature concerning the operation and maintenance of this unit has been included. Therefore, it should be read carefully in order to avoid any problems.



Stand-by



Direct Current


Safety Logo:



ALVA BC640 complies with the 93/42/EEC, EN60601-1, EN60601-2 of related European standards.

Disposing of your old product:

- ◆ Within the European Union



EU-wide legislation, as implemented in each Member State, requires that waste electrical and electronic products carrying the mark (left) must be disposed of separately from normal household waste. This includes monitors and electrical accessories, such as signal cables or power cords. When you need to dispose of your display products, please follow the guidance of your local authority, or ask the shop where you purchased the product, or if applicable, follow any agreements

made between yourself.

The mark on electrical and electronic products only applies to the current European Union Member States.

- ◆ Outside the European Union

If you wish to dispose of used electrical and electronic products outside the European Union, please contact your local authority so as to comply with the correct disposal method.

- ◆ Please follow the national guidelines for unit disposal.

Classification:

Power by Class I power adapter: 18UT-05-3000

- ◆ No Applied Part.
- ◆ No protection against the ingress of water: IPX0

- ◆ Mode of operation: Continuous Operation
- ◆ The equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide: Not AP or APG Category:

Product Specifications:

<i>Specifications</i>	
Power Adapter	Input Voltage Range: 100 – 240 VAC
	Input Frequency Range: 50 – 60 Hertz
	Fuse: Internal and not user replaceable
	Power Output: 5V 1.5A
	Max Input AC Current: 90-264VAC 0.3A
Environment	Operating Temperature: +0°C ~ +40°C
Conditions	Storage/Transportation Temperature: -20°C ~ +60°C
	Operation Relative Humidity: 20% ~80% Non-Condensing
	Storage/Transportation Relative Humidity: 10% ~95% Non-Condensing
	Atmospheric pressure: 500~1060 hPa

FCC :

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced technician for help.

You are cautioned that changes or modifications not expressly approved by the party responsible for compliance could void your authority to operate the equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference and (2) this device must accept any interference received, including interference that may cause undesired operation

FCC RF Radiation Exposure Statement:

This Transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. To maintain compliance with FCC RF exposure compliance requirements, avoid direct contact to the transmitting antenna during transmitting.

Any changes or modifications (including the antennas) made to this device that are not expressly approved by the manufacturer may void the user's authority to operate the equipment.

Warranty conditions

Optelec warrants the ALVA BC640, effective from the date of delivery, to be free of any defects in material and workmanship.

Warranty is not transferable and does not apply to groups, multi users or agencies.

The ALVA BC640 has been designed for the individual purchaser to be used in home and office environments. Optelec reserves the right to repair or replace any purchased ALVA BC640 with a similar or better product.

In no event shall Optelec or its distributors be liable for indirect or consequential damages. The original user's remedies are limited to replacement of ALVA BC640 modules. This warranty is valid only when serviced in the country of original purchase and with intact seals. For additional warranty claims or service during or after warranty period, please refer to your distributor.

Optelec takes no responsibility for use of this device other than described in this manual.

Optelec offices worldwide

Optelec U.S. Ltd

3030 Enterprise Court, Suite C
Vista, CA 92081
U.S.A.
T. +1 (800) 826-4200
www.optelec.com

Optelec Canada

1832, Boulevard Marie-Victorin
Longueuil, QC J4G 1Y9
Canada
T : +1 (800) 665-3005
www.optelec.com

Optelec Limited

2, Millfield House, Woodshots Meadow
Croxley Green Business Park
Watford, WD18 8YX
United Kingdom
T: +44 (0)1923 23 13 13
www.optelec.co.uk

Optelec Nederland BV

Breslau 4
2993 LT Barendrecht
The Netherlands
T: +31 (0)88 6783 555
www.optelec.nl

Optelec NV

Baron Ruzettelaan 29
B-8310 Bruges
Belgium
T: +32 (0)50 35 7555
www.optelec.be

Tieman GmbH

Fritzlarer Strasse 25
34613 Schwalmstadt
Germany
T: +49 (0)6691 96170
www.tieman.de

Optelec Corporate Headquarters: Breslau 4 - 2993 LT Barendrecht – The Netherlands – T: +31 (0)88 678 3444 – www.optelec.com