

Electronic Reading Aid




User's Guide

Open box and unpack contents

- 1 x Wireless Electronic Reading Aid (WERA)
- 1 x carry bag
- 1 x user's guide

A quick look at the WERA



1. Magnify
2. Mode
3. Freeze & D-freeze
4. De-magnify
5. Power on/off
6. Charging indicator: red-charging, off-charging finish
7. DC-in

How to use the WERA?

1. Insert the small end of the power adaptor into black end of the Y-shaped power/video cable from the Receiver.
2. Plug the power adaptor into a wall socket.
3. Insert the yellow RCA plug on the video cable into one of your television's yellow video input connectors. Refer to your television owner's manual for more information locating these connectors.
4. Turn on the television. If necessary, switch your television to the appropriate input channel.
5. Turn on the WERA by sliding the on/off button to ON position.
6. Put the WERA on newspapers, magazine etc and then you can read their big images from TV wirelessly.

7. Key definitions:

- 1) M: Switch reading mode among Full color / Black on White / White on Black
- 2) + : Magnify
- 3) - : De-magnify
- 4) * : Freeze & D-freeze

How to charge the Wireless Electronic Reading Aid?

When the battery electricity is low, the WERA will not be able to work normally, then you have to recharge its built-in lithium battery.

1. Switch off the WERA by sliding on/off button to off.
2. Plug adaptor into DC-in hole of the WERA.
3. Plug the power adaptor into a wall socket.
4. Then charging indicator will be in red; when the charging finished, the indicator will be off.

Safety:

- The WERA is not waterproof, so keep it dry.
- Use the WERA only at -5°C-50°C.
- The white LEDs which illuminate the WERA target field are very bright. Do not stare directly into these LEDs as it may damage your eyes.
- Use only the power adaptor that came with the WERA. Other adaptors may damage the unit.
- When battery charging is full, disconnect the adaptor with WERA timely.

Warranty Card

Your Reading Aid has a one year guarantee against any manufacturing defect!

Date of purchase:		Product No.	
Seller's info	Name:		
	Address:		
	Tel:		
Defect description			
Cause			
User's info	Name:		
	Address:		
	Tel:		



NOTICE

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Warning:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment