
# eScope Instruction Manual



# Computer System Requirement

Microsoft Windows XP, Vista or Windows 7
Pentium® 4 or AMD Athlon™ processor or Higher
100MB free hard disk space
512MB RAM or more
16Bit color Display
USB interface
CD or DVD Drive


#### TV Mode

TV or AV Monitor

#### Accesories

Adjustable Stand User Guide CD ROM TV Cable AC Adaptor USB Cable

### **Product overview**






# Copy Amcap Tool

 Insert enclosed CD into CD-ROM driver, Copy the "Amcap .exe" to your desktop or other location of your Computer. Then you can use the tool Capture photo and video.





#### PC Connection

Plug the small end of the USB cable into the eScope, plug the large end into your working computer.



When the TV receive box light flashing, turn on your eScope main unit.

After synchronize, TV box & main unit LED light up continuity.

\* First time connect USB Dongle to PC,PC will find new device and install the software of device automatically. Then prompt the device is ready to use.

### Viewing

1. Open "AMCap" icon in your PC.



2. If your computer has another webcam, you can select eScope from



3. In the pull down menu > option, you can select "Video Capture Filter" or "Video Capture Pin" which to set properties that you want.

# Focusing

Put eScope on the subject and turn Focus adjust knob slightly to find the sharpest focus position.




## Capture image

 a) After foucsing, you can capture image by Snap shot button

A Windows(PhotoDisplay) will appear, you click "SavePhoto" save the captured image to destination path.

 b) Click Photo>Start photo, A Windows (PhotoDisplay)will appear, you click "SavePhoto" save the captured image to destination path.

Save the photo is BMP format.



#### capture video

1. Click "Capture > start capture" on the pull down menu.



2. Confirm or change the video destination folder, click OK.



- 3. Confirm or change the video file name.
- 4. Capture start, until you stop it by pull down menu "Capture > stop capture".



#### TV Mode

- 1.Use AC adaptor power supply to TV receive box.
- Connect TV Cable between TV receive box and TV, Connect TV Cable yellow plug to TV yellow video input interface.
- 3.Set TV in AV mode Auto system.
- 4. Power on the main Unit; After synchronize, TV box & main unit LED light up continuity.
- 5. The real-time image appear on TV screen, you can see the object what you want.
- 6.In TV mode, you only see real-time image, can't save photo or Video.

# **Specifications:**

PC Mode:		
Image Sensor	0.3Mega Pixels	
Video Capture Resolution	640x480	
Photo Capture Resolution	640x480	
Prevrew Resolution	640x480,320x240,160x120	
Color	24bit RGB	
Magnification Rate	Microscope Lens 10X-20X,120X	
Focus Range	Manual Focus From 10mm to 500mm	
Main Device Power Supply Built-in Li-on 3.7V 400mAH		
Dongle Power Supply	5V DC from USB Port	
Battery Continues Working 2Hour		
time		
Transport distance	10M	
Load Frequency	2.4G	
Flicker Frequency	50Hz/60Hz	
Frame Rate	15fps@VGA	

LED	Main:2pcs; Red-working indication, Blue-Charge Indication
	TV box:1pc; Red-working indication
Button	Main: 1pc For Snap shot & Pair, TV box: 1pc For Pair
Photo Format	BMP; JPEG
Video Format	AVI format
White Balance	Auto
Exposure	Auto
Light Source	8 LED White Lights (High /Low )
Interface	USB 1.1
Operation System	Windows XP/Vista/Windows 7
Size & Weight	Main:130mm(L) x37mm(D)
	TV box: 80mm(L) x50mm(w)x12.6mm(H)
AV Mode:	
Adaptor	AC 100-240V Input, 5V DC Output
Video Interface	TV(Auto system) or AV monitor

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not causeharmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

#### NOTE:

Changes or Modifications not expressly approved by the party responsible could void the user's authority to operate this device.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- -- Reorient or relocate the receiving antenna.
- -- Increase the separation between the equipment and receiver.
- -- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- -- Consult the dealer or an experienced radio/TV technician for help.



