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RENPHO <https://renpho.com/pages/help-and-manual>

(EN) RENPHO promotes environmental conservation by minimizing paper use. Access the full electronic User Manual via the QR Code or URL.

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QUICK START GUIDE

| Button / Display | Description |
|------------------|---|
| | Reset Button Press to reset the scale. Note: Resetting the scale will erase all stored user data and settings. |
| | Left / Volume Key Tap to cycle from volume off / volume low / volume middle / volume high or scroll content left. |
| | Right / Unit Conversion Key Tap to change the unit or scroll content right. |

MorphoScan Nova

Model: R-MSC04

| Button / Display | Description |
|------------------|--|
| | BMI |
| | Body Fat Percentage |
| | Skeletal Muscle Percentage |
| | Body Water Percentage |
| | Protein Percentage |
| | Bone Percentage |
| | Weight |
| | Units Switched via the Unit Conversion Key. |
| | Overloaded Maximum Capacity: 400 lbs / 180 kg. |
| | Low Battery Charge the scale. |

USING YOUR PRODUCT

Note: Before first use, press the Reset Button or connect the scale to a power source using the provided charging cable and charger (not included) to initialize the scale.

- Place the scale on a flat, hard surface and step on the scale's electrodes with bare, dry feet. Ensure there is space between your knees and thighs.

- Grasp the handle with your thumbs on the Thumb Electrodes and the other fingers on the Four-finger Electrodes.

- Lift and hold the handle to the proper level with your arms naturally hanging, slightly bent and away from your body and abdomen. Stand upright, keep still and wait for your final measurement(s) to be shown on the scale's display.

- After measuring, carefully put the handle back to the scale in place. The scale will transmit the data to the app. Open the app to sync and view your detailed measurements, track your progress, and access additional features and insights.

* Create your correct user profile in the RENPHO Health App before taking your first measurement. This ensures your measurements are precisely recorded and associated with your profile.

FCC Regulatory Compliance

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ISED Regulatory Compliance

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- 1: This device may not cause interference.
- 2: This device must accept any interference, including interference that may cause undesired operation of the device.

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