User Manual


Smart Body Weight Scale

Please read this User Manual carefully before using the device


IOS QUICK START GUIDE


1. Search "Vitafit" via Apple Store.

2. Convenient login by Apple or Facebook account.


3. Choose suitable Mode. Note: Weight Only Mode specialized for pregnant or having a pacemaker, etc.

4. Turn on Bluetooth on your phone.

5. Pairing with Vitafit Scale.


Notes:

1. For body fat and heart rate measurement, make sure your feet bare and dry when measuring.
2. For excellent experience, recommend to reassemble batteries for the first beginning use.
3. UI will be different depended on scale type: Smart Weight, Smart Body Fat or Smart Body Fat & Heart.



IOS QUICK START GUIDE



6. Turn on VitaFit App and step on for measuring.

7. Measurement data.

8. Set up for sync with Apple Health.


9. View measurement data in Apple Health after sync.


Notes:

- For body fat and heart rate measurement, make sure your feet bare and dry when measuring.
- For excellent experience, recommend to reassemble batteries for the first beginning use.
- UI will be different depended on scale type: Smart Weight, Smart Body Fat or Smart Body Fat & Heart.


ANDROID QUICK START GUIDE


1. Search "Vitafit" via Google Play Store.


2. Convenient login by Google or Facebook account.


3. Choose suitable Mode. Note: Weight Only Mode specialized for pregnant or having a pacemaker, etc.

4. Turn on Bluetooth and Location on your phone, **OR ELSE THE APP DON'T GET RESULTS ACCORDING TO ANDROID REGULATIONS.**

5. Pairing with Vitafit Scale.


Notes:

1. For body fat and heart rate measurement, make sure your feet bare and dry when measuring.
2. For excellent experience, recommend to reassemble batteries for the first beginning use.
3. UI will be different depended on scale type: Smart Weight, Smart Body Fat or Smart Body Fat & Heart.



ANDROID QUICK START GUIDE



6. Turn on VitaFit App and step on for measuring.

7. Measurement data.

8. Set up for sync with Google Fit.

9. View measurement data in Google Fit after sync.

Notes:

1. For body fat and heart rate measurement, make sure your feet bare and dry when measuring.
2. For excellent experience, recommend to reassemble batteries for the first beginning use.
3. UI will be different depended on scale type: Smart Weight, Smart Body Fat or Smart Body Fat & Heart.

CONTENTS

01-06		English
07-12		Deutsch
13-18		Français
19-24		Español
25-30		Italiano

Warning

1. This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
2. Children shall not play with the appliance.
3. Cleaning and user maintenance shall not be made by children without supervision.
4. Model of battery that appliance used is 3×AAA.
5. Non-rechargeable batteries are not to be recharged.
6. Different types of batteries or new and used batteries are not to be mixed.
7. Batteries are to be inserted with the correct polarity.
8. Exhausted batteries are to be removed from the appliance and safely disposed of.
9. If the appliance is to be stored unused for a long period, the batteries are removed.
10. The supply terminals are not to be short-circuited.
11. Take the scale, turn it over, open the battery compartment, put 3×AAA batteries according to the right polarities, then return the battery compartment cover.

VitaFit Smart Body Weight Scale

Thank you for choosing VitaFit Smart body weight scale. This scale is your personal health assistant that could track and record your weight. We sincerely hope that you enjoy our VitaFit scale.

Safety Cautions


1. This scale is not intended for commercial use.
2. This scale should not be used to diagnose or treat any medical condition. Consult your primary care physician before making changes to your diet, exercise plans, or physical activities.
3. If you are pregnant, have a pacemaker, or other internal device, you should not heavily rely on the data provided by this product. Please note that all data should be used as a reference and is not substitute for medical advice from a licensed healthcare professional.
4. Check the scale before use. Do not use the scale if it is damaged in any way.
5. The glass platform is slippery when wet. Make sure the platform and your feet are dry before stepping on for measurement.

Use and Care

1. To weigh safely and accurately, stand with your feet shoulder-width apart while on the scale.
2. DO NOT strike or drop the scale.
3. DO NOT stand on the edge of the scale or jump on it.
4. DO NOT stand on the scale platform for a prolonged period of time or when not in operation.
5. DO NOT store or operate near heat sources or in humid environments.
6. Handle the scale with care, as it is a precision instrument.
7. Clean the scale with a clean and dry cloth. Avoid getting water or any other liquid inside the scale.
8. DO NOT use chemical or abrasive cleaning agents to clean the scale
9. Keep the scale in a cool and dry environment.
10. Always keep the scale in a horizontal position, do not store vertically.
11. If the scale fails to turn on, check whether the batteries are installed correctly. Replace the batteries if necessary.
12. If an error appears on the display or if the scale does not turn off automatically, remove the batteries for over 30 seconds, then reinstall them.
13. DO NOT overload the scale (maximum weight capacity: 182kg/400lb/28st).
14. DO NOT mix new and used batteries or different battery brands to power the scale.

First Time Setup

IMPORTANT! BEFORE USE: Remove the isolator tab from the battery compartment.


Operation

Weighing without VitaFit App

1. Remove the isolator tab from the battery compartment.
2. For best accuracy, place the scale on a hard and flat floor. Carpeted or uneven floors can negatively effect the accuracy .
3. The scale operates with step-on technology but needs to be calibrated automatically under following occasions.
 - A. when first time use, moving the scale or battery exchange, the scale will display scale itself weight, and appear "CAL" after around 1second, and then display "0.0" or "0.00" after around 2seconds , now it's ready to step on to weigh normally.

B. If step on the scale before calibration, the scale will display an wrong weight and appear "CAL" after around 1 second, and then display "0.0" or "0.00" after around 2seconds.

After calibration, please step again on to weigh and get correct weight.

4. Under other occasions, please step on to weigh directly.
5. When weighing, please remain standing evenly on the scale without shifting feet or body to get stable and right weight.
6. Step off the scale, the weight readings remain for several seconds. Then the scale will turn off automatically.
7. When Weighing without App, switch the unit by press the unit button on the backside of scale.

Weighing with VitaFit App


The VitaFit Bluetooth Scale can connect with your mobile phone or Ipad via Bluetooth. The Body Weight and BMI data will be showed at the App named "VitaFit", which is easy to download from IOS App Store or Google Play.

Download App and Register ID

Scan the QR code or Download "VitaFit" via IOS App Store or Google Play.

Register Your own ID, Then add personal data and confirm OK.

Note: When weighing with VitaFit App, please switch the unit via App, the scale will keep same unit with App. Let's make more clear, the unit button on the backside of scale will be nonfunctional when the scale connected to VitaFit App.

Pairing the VitaFit Scale with your smartphone

Turn on Bluetooth® on your smartphone.

NOTE: For Android users, Location and Bluetooth® must be turned on for Android 6.0 and later.

1. Tap + on the measuring page to discover the scale.

2. Step onto the scale, then add the device named "VitaFit BodyWeight".

3. Keep standing on the scale for 5-10 seconds. Pairing is successful when your final weight measurement appears on the VitaFit app.

Viewing Reports

1. At the bottom you will find a bar when you can click compare data to display.


2. Click Weight or BMI bar, you can view history data for the previous three times.

3. To view comparison data for a specific time, select starting date and ending date to generate the report, you can save or share your progress with friends via Facebook and Instagram by Click the "share" button on this page.

4. Click Trend to view your progress history. You can select your progress for the week, month, or year.

5. To view data for a specific day, click calendar logo the top right corner of the screen.

6. To view reference data from a specific time frame, click Data Select, then select date.

Deleting Reference Data

Adding/Deleting Users and Other Setting


1. To delete a reference data, click Select and select date, then click the "X".


1. On the Measurement page, tap My Account at the bottom of the screen.

2. To add new users, click Member and select "Add a user" and then enter their information.

3. To delete a user, first click Edit and delete, and then confirm.

4. To delete the account, first click SET logo and click "Delete Account".


5. other setting: You can set goal weights, change weight units, select sounds, select languages, and set weighing reminder on the current page.

Connecting with Other Fitness Apps

Google Fit


Authorization Steps on the VitaFit App

1. On the page of My Account, click "Health Connect".
2. Please Install Health Connect if the device system below Android 14, and then start authorize.
- If the system is Android 14 and above, please skip the install step, just directly click "Authorize Immediately".
3. click all allow button, then the authorization is successful.


Connect steps on Google Fit

1. Open the Google Fit, and click "Get started".
2. Click "Set up".
3. Click "Done" to finish connected.
4. After successfully connected, your measurement data will sync with Google Fit automatically.




Samsung Health

1. Firstly, please make same steps of authorization on the VitaFit app as above Google Fit mentioned.
2. After VitaFit app connected with Health Connect, please start Settings on Samsung Health as follow steps.
 - 2.1. Open the Samsung Health, click "Settings"
 - 2.2. Click "Health Connect"
 - 2.3. Click "App permissions"



- 2.4. Click "VitaFit" to allow access
- 2.5. Click "Allow all" button to finish connecting
- 2.6. After successfully connected, your measurement data will sync with Samsung Health automatically

Apple Health



1. On the page of My Account, Select Apple Health from the menu.
2. Click Authorize and the App will authorize immediately.
3. Open the Health App on your iPhone, then tap Sources. Be sure to turn on your Bluetooth®.
4. Select Vitality from the list.
5. Turn on all categories to allow the Vitality app to work.

Note: In IOS system, couldn't be allowed to Sync with Google Fit.


Fitbit

1. On the page of My Account, Select Fitbit from the menu.
2. Turn on Fitbit, then sign into your Fitbit account.

3. Tap Allow to access the Fitbit app.
4. After successfully connecting, your data will sync to the Fitbit app. You will be able to view your progress both on the app as well as online.

Frequently Asked Questions

1. Why doesn't the scale turn on?
 - Check if the batteries are properly installed.
2. Why does the data flash on the screen and then disappear?
 - The batteries are low and need to be replaced.
3. Why doesn't the scale connect with the VitaFit app?
 - Make sure your phone's operating system is running on iOS 9.0+ or Android 4.3+.
 - Download and open the newest version of the VitaFit app.
 - Turn on Bluetooth® on your phone (for Android users, make sure your Location is turned on as well).

Package Contents

Serial Number	Item	Quantity
1	Smart Body Weight Scale	1
2	1.5V AAA Battery	3
3	User Manual	1

Specifications

Capacity: Min 5 kg/11lb/0.11lb

Max 182kg/400lb/28st

Weight unit: kg/lb/st

Graduation: 0.05kg/0.1lb/0.1lb

Dimensions: 280x280x24mm or 300x300x24mm(Depended on the model)

Display: LCD or LED (Depended on the model)

Platform: Tempered glass

Batteries: 3xAAA

Warranty Information

VitaFit guarantees all products to be of the highest quality in material and craftsmanship, and provides all products with 3 years warranty effective from the date of purchase.

For more information please visit www.51vitaFit.com, www.vitaFit.com

Email: service@51vitaFit.com

FCC Warning

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE 1: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

NOTE 2: Any changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

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Mail: service@51vitafit.com

“This device complies with Industry Canada's license-exempt RSSs. Operation is subject to the following two conditions:

- (1) This device may not cause interference; and
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.”

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- (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

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Maximum operating
temperature :45℃

Manufacturer company: Shenzhen Shine Industrial Co.,Ltd.

Address : 2-3/F, Bldg 5,1st Industrial Zone, Changzhen Community, Yutang town, Guangming District, Shenzhen, Guangdong, China

Importer company: Shenzhen Shine Industrial Co.,Ltd.

Address : 2-3/F, Bldg 5,1st Industrial Zone, Changzhen Community, Yutang town, Guangming District, Shenzhen, Guangdong, China