



# iRemember User Guide

Congratulations,  
you're committed to  
better health.




- 1. Status button
- 2. USB cover
- 3. Speaker



- 4. Pods (x7)
- 5. Sound button
- 6. Timer button


1

Open USB cover (2) to  
connect charger.




2

Once charged, open cap  
and pods by turning  
counter-clockwise.



3

Stack pods and place  
upcoming pod on top.




4

**Set Light & Sound**

Press left button (5) to  
select:

Sound On  
Flight Mode  
Sound Off



5


**Set timer**

Press right side of  
button (6) until  
selection is heard.  
To Cancel timer, select  
"Clear"

**Hours**

· 24  
· 12  
· 6\*  
· Clear


\* Start with morning dose



6


For 'Last opened'  
notification press status  
button (1):

Ok  Missed 



# Tips



Sort 7 days of pills and place doses in pods.



Download our app and get the most out of iRemember.

Search ‘iRemember’ in the Google Play App Store.

 Remembers location of last connection

 Activity log  
 Open/Closed and Status button

 Battery level indicator

# FAQ

## 1. Do I need a mobile phone to use iRemember?

No, iRemember will work by itself for daily reminders. For additional features, pair with mobile phone app.

## 2. Does iRemember work with multiple medicines?

Yes, multiple pills and multiple doses can be organized with iRemember.

## 3. May I use liquids with iRemember?

No, liquids are not permitted.

# Certifications

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help. This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. The device has been evaluated to meet general RF exposure requirement. The device can be used with minimum distance 0mm between the radiator and your body.

IC WARNING  
This radio transmitter has been approved by Industry Canada.  
This Class B digital apparatus complies with Canadian standard ICES-003. This device complies with Industry Canada License exempt RSS standard(s). Operation is subject to the following two conditions.  
This device may not cause interference and (2) this device must accept any interference, including interference that may cause undesired operation of the device.  
This equipment complies with ISED radiation exposure limits set forth for an uncontrolled environment.  
The device has been evaluated to meet general RF exposure requirement. The device can be used with minimum distance 0mm between the radiator and your body.

Industrie Canada (IC)  
Le présent appareil est conforme aux CNR d’Industrie Canada applicables aux appareils radio exempts de licence. L’exploitation est autorisée aux deux conditions suivantes: (1) l’appareil ne doit pas produire de brouillage, et (2) l’utilisateur de l’appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d’en compromettre le fonctionnement.

AVISRELATIFS À L’EXPOSITION RF  
Le dispositif a été évalué et homologué comme répondant aux exigences d’exposition RF. L’appareil peut être utilisé avec une distance minimale de 0 mm entre le radiateur et votre corps.

# Contact us!

We are happy to help you.

 [hi@wealthtaxi.com](mailto:hi@wealthtaxi.com)  
 [@WealthTaxi](https://twitter.com/WealthTaxi)

6 Month Manufacturer Warranty.

FC IC CE

