



Size:248\*159\*64mm  
250克铜版纸（双铜） 1200灰板

产品使用说明书

超声波洁牙仪

·感谢您购买我们的产品。

·使用前请务必仔细阅读此说明书并妥善保管，以备将来参考。

·本产品仅用于家庭个人牙齿清洁，不可用于任何医疗用途。

产品介绍

感谢选择本产品作为您的口腔健康助理。请按下图所示产品清单检查相关部件是否齐备。

平头洁牙头 X2

尖头洁牙头 X2

充电器

底座

刮舌器

探针

口镜

小扳手

保护罩

充电器底座

刮舌器

探针

口镜

小扳手

产品使用

警告

开启洁牙仪时,切勿将金属清洁头对准眼睛等人体敏感部位以及电源、灯具等物。产品开启时,请勿用皮肤触碰金属清洁头以免受伤。

注意:首次使用前请将产品充满电,因为出厂时为了安全考虑,主机并未充电。充电时,产品不能使用。

拆装清洁头

1. 安装清洁头

将清洁头对准底座上的金属触点,轻轻按下,直到听到“咔哒”声,表示安装成功。

2. 拆卸清洁头

将清洁头对准底座上的金属触点,轻轻向上提起,直到听到“咔哒”声,表示拆卸成功。

切换模式

1. 长按电源键 3 秒,进入模式切换菜单。

2. 使用方向键选择所需模式,并按下确认键。

3. 屏幕将显示所选模式的名称。

清洁牙齿

1. 清洁牙龈上方牙龈边缘的牙垢和牙石。

2. 清洁牙龈下方牙龈边缘的牙垢和牙石。

3. 清洁牙齿表面,特别是牙缝和牙冠部分。

4. 清洁牙齿表面,特别是牙缝和牙冠部分。

5. 清洁牙齿表面,特别是牙缝和牙冠部分。

6. 清洁牙齿表面,特别是牙缝和牙冠部分。

产品充电

将底座充电器侧面的 USB 接口通过 USB 充电线连接至电脑或充电器。

充电时,充电指示灯将亮起。

充满电后,充电指示灯将熄灭。

产品参数

产品名称: 超声波洁牙仪

额定功率: 3W

防水等级: IPX6

产品规格: 3档调节

额定电压: 5V

无线充电电流: 1A

电池容量: 800mAh

产品尺寸: 208.5\*28\*28mm

重要注意事项

危险

1. 如果电源插头或电源线损坏,请勿使用。

2. 请勿将本产品放入水中或任何液体中。

3. 请勿将本产品放入火中或任何热源附近。

4. 请勿将本产品放入微波炉中。

5. 请勿将本产品放入洗碗机中。

6. 请勿将本产品放入洗衣机中。

7. 请勿将本产品放入烘干机中。

8. 请勿将本产品放入烘干机中。

9. 请勿将本产品放入烘干机中。

10. 请勿将本产品放入烘干机中。

11. 请勿将本产品放入烘干机中。

12. 请勿将本产品放入烘干机中。

警告

1. 请勿将本产品放入水中或任何液体中。

2. 请勿将本产品放入火中或任何热源附近。

3. 请勿将本产品放入微波炉中。

4. 请勿将本产品放入洗碗机中。

5. 请勿将本产品放入洗衣机中。

6. 请勿将本产品放入烘干机中。

7. 请勿将本产品放入烘干机中。

8. 请勿将本产品放入烘干机中。

9. 请勿将本产品放入烘干机中。

10. 请勿将本产品放入烘干机中。

11. 请勿将本产品放入烘干机中。

12. 请勿将本产品放入烘干机中。

注意

1. 请勿将本产品放入水中或任何液体中。

2. 请勿将本产品放入火中或任何热源附近。

3. 请勿将本产品放入微波炉中。

4. 请勿将本产品放入洗碗机中。

5. 请勿将本产品放入洗衣机中。

6. 请勿将本产品放入烘干机中。

7. 请勿将本产品放入烘干机中。

8. 请勿将本产品放入烘干机中。

9. 请勿将本产品放入烘干机中。

10. 请勿将本产品放入烘干机中。

11. 请勿将本产品放入烘干机中。

12. 请勿将本产品放入烘干机中。

故障排除

故障现象

故障原因

解决方法

产品合格证

品名: \_\_\_\_\_ 检验结论: \_\_\_\_\_

生产日期: \_\_\_\_\_ 检验员: \_\_\_\_\_

产品保修卡

用户名字: \_\_\_\_\_

购买时间: \_\_\_\_\_

联系方式: \_\_\_\_\_

用户地址: \_\_\_\_\_

故障: \_\_\_\_\_

Instruction Manual

Ultrasonic therapeutic apparatus

· Thank you for purchasing our products.

· Be sure to read this manual carefully before use and keep it safe for future reference.

· This product can only be used for personal dental cleaning at home, not available for any medical purpose.

Product description

Thank you for choosing this product as your oral health care assistant. Please check the product list shown below to check whether the relevant parts are available.

Flat head cleaning head x 2

Tip cleaning head x 2

Charging Cable

Protective cover

Scaler hand

Charging base

Tongue scraper

Probe

Mouth mirror

Small wrench

Protective cover

Scaler hand

Charging base

Tongue scraper

Probe

Mouth mirror

Small wrench

Product use

Warning

When turning on the scaler, do not point the metal cleaning head at sensitive parts of the human body such as eyes, as well as power supplies, lamps and other equipment. When the product is turned on, please touch the metal cleaning head with your skin to avoid injury.

Note: Please fully charge the product before using it for the first time, because the main unit is not fully charged when it leaves the factory for safety reasons. While charging, the product cannot be used.

Disassembling the cleaning head

1. Install the cleaning head

Align the cleaning head with the interface on the top of the transducer, use a small wrench to gently rotate and tighten it clockwise.

2. Disassemble the cleaning head

Align the cleaning head with the interface on the top of the transducer, use a small wrench to gently rotate and loosen it counter-clockwise to replace or remove.

Switch mode

1. Long press and hold the Power Mode switch button to enter the mode selection menu.

2. Use the direction keys to select the desired mode, and then press the confirm key to enter the selected mode.

3. The screen will display the selected mode name.

Clean the teeth

1. When cleaning the calcification, select the Scaler or the Probe. Do not use the Mouth Mirror.

2. When scraping the surface of calcification, use the Scaler or the Probe. Do not use the Mouth Mirror.

3. When using the Mouth Mirror, use the Scaler or the Probe. Do not use the Mouth Mirror.

4. When using the Mouth Mirror, use the Scaler or the Probe. Do not use the Mouth Mirror.

5. When using the Mouth Mirror, use the Scaler or the Probe. Do not use the Mouth Mirror.

6. When using the Mouth Mirror, use the Scaler or the Probe. Do not use the Mouth Mirror.

Product charging

Connect the USB port on the side of the base charger to the power adapter via a USB charging cable, then place the scaler handle into the base to charge.

When fully charged, the charging indicator is always on.

When fully charged, the charging indicator is always on.

When fully charged, the charging indicator is always on.

When fully charged, the charging indicator is always on.

When fully charged, the charging indicator is always on.

When fully charged, the charging indicator is always on.

Important Notes

Danger

1. Do not use if the power plug or power cord is damaged.

2. When turning on the scaler, do not point the metal head at sensitive parts of the human body such as eyes, and damaged equipment such as power supplies and lamps.

3. Do not touch the skin with the metal cleaning head for a long time to avoid injury.

4. The product contains a built-in rechargeable battery, please do not heat in fire or charge, use or place it at high temperature.

5. The product is not suitable for people under the age of 18, please keep this product away from children.

6. Others who do not have full capacity for civil conduct or lack relevant experience and knowledge must be under the supervision and guidance of the supervisor use, clean and maintain the product.

Warning

1. Do not immerse the product in water or any liquid. If it falls into the water when connected to the power supply, please unplug the power plug first and then take out the product.

2. Do not press the metal cleaning head hard on the gums or teeth to avoid damage.

3. If you need to charge, please confirm that the power plug and the interface are completely dry, and then connect the power.

4. Do not modify, disassemble or repair this product by yourself.

5. Before connecting the power supply, please check whether the local power supply voltage is consistent with the voltage marked on the product.

6. Do not use the product in a bath or bath.

7. Do not wrap the power cord around the product when it is powered on.

8. Ultrasonic scaling is contraindicated in patients with cardiac pacemakers to avoid symptoms such as dizziness and arrhythmia caused by the interference of electromagnetic radiation.

Note

When using this product for the first time, there may be soreness, itching or slight bleeding. This is because the gums are not used to the efficient cleaning method. It is recommended to use the gentle mode of the product to clean first. This phenomenon will generally stop within 1-2 weeks. If the above phenomenon improves significantly over 2 weeks, stop using this part and consult a dentist in time.

Troubleshooting Summary

Troubleshooting Summary

Trouble causes

Cause Solutions

Product certification

Product name: \_\_\_\_\_ Inspection conclusion: \_\_\_\_\_

Production date: \_\_\_\_\_ Inspector: \_\_\_\_\_

Username: \_\_\_\_\_

Purchase time: \_\_\_\_\_

Contact information: \_\_\_\_\_

User address: \_\_\_\_\_

Fault: \_\_\_\_\_

FCC Warning Statement: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC Radiation Exposure Statement

The antennas used for this transmitter must be installed to provide a separation distance of at least 20 cm from all persons and must not be co-located for operating in conjunction with any other antenna or transmitter.

Size:100\*80mm  
128克的铜版纸