

Product Manual



SPECIFICATION

| Item | Description |
|------------------------|---|
| Wi-Fi Type | Wireless LAN 802.11 b/g/n |
| OS Support | Android, iOS, Windows |
| Image Sensor | 1/2.9 inch Color CMOS (2.19M pixels) |
| Sensor Pixel | Effective pixel : 1985×1105Pixel |
| Image Capture Size | VGA (640 X 480) |
| Image Frame Rate | VGA (640X480) - 30Frame Max |
| LED | Chip LED |
| Camera Button Function | Power, Image Capture |
| Measurement Mode | Moisture, Roughness, Pore, Oil, Sebum, Dead skin, Sensitivity, Pigment, Wrinkle |
| Battery | Type: Li-Polymer 960mmAH Normal Voltage: 3.7V |
| DC IN Adaptor | 5V DC / 2A |
| Dimension | 50.4(Ø)X87.3(H) mm |
| Weight | About 120g |
| Cradle | Stand |
| Operating Temperature | -10 ~ 50°C |
| Storage Temperature | -30 ~ 50°C |
| Operating Humidity | 5 ~ 85% |

PRODUCT COMPONENTS

- mi.moa device
- C-Type charger with cable
- Cradle

DIMENSION & KEY



OK- OK key (Capture key)
> Capture / Retry key

P - Power key
> Power On / Off (Press and hold for more than 2 seconds)

YOUR ONLY ONE SKINCARE MANAGER

Skin Measuring Instrument - mi.moa



mi.moa

mi.moa

The right standard for the emerging beauty market
mi.moa skin measuring instrument



Suggests personalized skincare solutions
based on accurate analysis



Checks various categories of your skin with a
professional camera developed by an expert



Records and compares your skin data in real-time

Needs



When you are too busy to visit skincare professionals



When you find it hard to afford frequent skincare services



When you do not know which home care products
to use because there are just too many cosmetic
products in the market




If you want to find the right skin care product
that works best for your skin,

It is essential to understand
YOUR SKIN CONDITION FIRSTLY.



Features

Reliable analysis results for various skin categories




If you understand your skin, you no longer need to worry about any skin troubles



Optical Microscope for Professional Inspection

Looking into all three layers of your skin



Three layers of skin



Outer skin



Stratum basal



Sebum production in pores



Personalized functions for continuous skincare



Push notifications sending skincare information

A community for users sharing beauty information


Chat consultations with skin experts

Ranking service among mi.moa users

User-friendly UI



Reliable Data



BIG DATA



To produce personalized skincare solutions based on specific skin conditions, environmental factors, lifestyles, etc.

The more tests you have, the more reliable your data will be. We will understand and approach your skin condition based on 100,000 test data.

User Guide




Open the mobile application and select categories that you want to measure.



Gently place the device on your skin.



mi.moa analyzes your skin condition by looking into all three layers of your skin.




Based on your skin condition, mi.moa suggests a personalized skincare method.

How to download the app

<http://isee.kcscopes.com/mimoa/auth/download.html>


Access the above link on your Android device and click the download image to download the APK.
If you download it, you must log in after registering as a member to use it normally.

Login



- This is the first screen that appears when the app is launched.
- Please click Sign Up to proceed with membership registration.
- After that, you can log in by entering your registered ID or E-mail and Password.

Sign Up



mi.moa

당신에 대해 알려주세요.
별표(*)표시는 필수입력 항목입니다.

* 이름을 띄어쓰기 없이 적어주세요.

홍길동

* 당신의 성별을 선택해주세요.

남성 여성

* 당신의 생년월일을 입력해주세요.


20020106

* 주소를 선택해주세요.

서울특별시 강남구

Already registered? [Login](#)

Next button highlighted with an orange box.



mi.moa

프로필 사진

Upload

당신의 피부 고민을 선택해주세요. (복수 선택)

Pore Acne

수부지 Sensitivity

Oil Wrinkle

당신의 두피 고민을 선택해주세요. (복수 선택)

Oily Dry

Inflammatory 탈모

각질 가려움증



mi.moa

안전한 비밀번호입니다.

* 비밀번호 확인
비밀번호가 일치합니다.

* 닉네임 입력하신 닉네임은 사용가능합니다.
hong 중복확인


* 이메일 입력하신 이메일은 사용가능합니다.
honggd@gmail.com 중복확인

* 휴대폰번호
1234 인증하기

* 인증번호
인증번호 입력 확인

Already registered? [Login](#)

Next button.



mi.moa

어떠한 직군에서 일을 하시나요?
(주변 환경에 따른 더욱 효과적인 맞춤형 피부 솔루션을 위한 질문입니다.)

IT·인터넷

mi.moa 인증번호로 인증하기 (선택)
기기 인증을 해아 mi.moa 측정이 가능합니다.
기기번호와 인증번호는 재를 보증을 확인하여 주세요.

기기번호
z1234567

기기인증번호
기기인증번호 인증하기

전체동의 ☒


[필수] moa 이용약관 및 개인정보처리방침 내용을 확인하였으며, 이에 동의합니다. ☒

[선택] 광고성 정보 수신 전체 동의 ☒
마케팅 정보 수집에 동의하시면, 이벤트 및 할인혜택 안내를 누구보다 빠르게 받아보실 수 있어요!

동의하고 가입하기


- Items marked with an asterisk are required fields.
- On screen 3, you need to double-check ID, nickname, and email. Cell phone number verification is not performed. Please enter a password of 6 or more digits.
- Please enter 6 or more random digits for the device number on the 4th screen. If the device number is not entered, subsequent device access is restricted.

Home




- This is the home screen you see after you log in normally.

Settings – LED on / off



- You can run the setting screen by clicking the rightmost setting button at the bottom.
- Click Settings > LED to turn LED on / off.
However, the LED on/off changes only when mi.moa and the device are connected.
- If it is not connected to mi.moa, the wifi setting popup is displayed.



FCC Information to User

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution

Modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

FCC Compliance Information : This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation

Specific Absorption Rate (SAR) information:

This Device meets the government's requirements for exposure to radio waves. The guidelines are based on standards that were developed by independent scientific organizations through periodic and thorough evaluation of scientific studies. The standards include a substantial safety margin designed to assure the safety of all persons regardless of age or health.

FCC RF Exposure Information and Statement The SAR limit of USA (FCC) is 1.6 W/kg averaged over one gram of tissue. Device types: **Skin tester (FCC ID: 2A5D5MMS1)** has also been tested against this SAR limit. The highest SAR value reported under this standard during product certification for properly worn on the body is **0.527 W/kg**. This device was tested for typical body-worn operations with the back of the handset kept **0 cm** from the body. To maintain compliance with FCC RF exposure requirements, use accessories that maintain a **0 cm** separation distance between the user's body and the back of the handset. The use of belt clips, holsters and similar accessories should not contain metallic components in its assembly. The use of accessories that do not satisfy these requirements may not comply with FCC RF exposure requirements, and should be avoided.

Body-worn Operation

This device was tested for typical body-worn operations. To comply with RF exposure requirements, a minimum separation distance of **0 cm** must be maintained between the user's body and the handset, including the antenna. Third-party belt-clips, holsters, and similar accessories used by this device should not contain any metallic components. Body-worn accessories that do not meet these requirements may not comply with RF exposure requirements and should be avoided. Use only the supplied or an approved antenna.