

Aumacom Sonic Electric Toothbrush Operating Manual



Product Specification:

1. Product Name: Sonic Electric Toothbrush Operating Manual.
2. Use environment: This product is used for indoor bathroom.
3. Product type: 5VDC, 1.7A.
4. Product size:
Handle: 30.45*30.99*183.39 mm;
Charging base: 49.72 * 59.97 * 33.7mm;
Brush head: 18.3 * 18.3 * 85mm.
5. Product weight:
charging holder: 100g ;
Handle: 116g;
Toothbrush head: 4g .
6. Operating ring temperature: 0~40℃.
7. Storage temperature: -10 ~ +60℃.
8. Waterproof ability: Handle IPX7, base IPX7.
9. Working voltage: 3.7V.
10. Working current: 600max.
11. Working time: the continuous working hours of the battery shall be divided into ≥ 60 minutes.
12. Charging time: ≤ 10 hours,
battery specification:
Li-ion battery 3.7V 600mAh.
13. Static current: $\leq 30\mu\text{A}$: $\leq 30\mu\text{A}$.

14. Function: 30S reminder, 2 minutes automatic shutdown function.
15. The functional indicators are white; Charging indicator light: red and green.
16. Operation dry sound: distance 30cm \leq 60dBA.
17. Pulling force: The pulling force of single plant should be \geq 1.8kg.
18. Packaging requirements: Carton packing + blister packing (1X host, 1X charging base, 2X brush head (including color ring), 4 color ring).
19. Safety requirements: Meet the safety requirements of CE/ROHS/REACH/WEEE/EMC.

FCC Warning

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. During the operation of device a distance of 15 cm surrounding the device and 20 cm above the top surface of the device must be respected.