

正面

   <p>Thanks for choosing Comfytemp</p> <p>SKU: K5102</p>	<h3>Specifications</h3> <table border="1"> <tr><td>Item Name</td><td>Smart Pill Organizer</td></tr> <tr><td>Weight</td><td>0.33lb</td></tr> <tr><td>Product Size</td><td>4*3.5*1.14in</td></tr> <tr><td>Material</td><td>Food-grade materials PP, ABS</td></tr> <tr><td>Color</td><td>White</td></tr> <tr><td>Compartments</td><td>4</td></tr> <tr><td>Input</td><td>5V/1A</td></tr> <tr><td>Battery Capacity</td><td>1000mAh</td></tr> <tr><td>Package Includes</td><td>1*Smart Pill Organizer 1*USB-C Charging Cable 1*User Guide</td></tr> </table>	Item Name	Smart Pill Organizer	Weight	0.33lb	Product Size	4*3.5*1.14in	Material	Food-grade materials PP, ABS	Color	White	Compartments	4	Input	5V/1A	Battery Capacity	1000mAh	Package Includes	1*Smart Pill Organizer 1*USB-C Charging Cable 1*User Guide	<h3>Product Diagram</h3> 	<h3>Product Features</h3> <ol style="list-style-type: none"> Display of number of pills taken; Medication reminder; Storage of medication records; Memory function: the alarm clock set before the shutdown will be restored by default when the kit is turned off and turned on again; Reminding method: sound effect+light reminder+vibration; <h3>How to Use</h3> <ol style="list-style-type: none"> Short press ON/OFF to start the machine, long press the ON/OFF button for 3s to shut down the machine, press the ON/OFF button three times to enter the DUT mode, and then press the ON/OFF button three times to exit the DUT mode; Display of the number of pills taken: within 30 m after the time point of taking the medicine, only the number of pills taken in the set medicine bin will be displayed after opening the medicine box, and other medicine bins will not be displayed. During this period, the number of pills taken in this period will not be displayed if the medicine bin is opened once and then closed again; End the display after closing the outer cover of the kit; When the time for taking medicine is up, the medicine box will emit a drip sound and vibration lasting for 1 minute, and the green light will last for 1 minute. The above behaviors will stop when the user opens the medicine box or presses the button of the medicine box; If there is no action to open the kit, repeat the above reminders after 10m, a total of 3 reminders; Open the medicine box within 30 minutes of the set time reminder to record that the medicine has been taken; Storage of medication record: open the cover for the first time within the medication time, and trigger a medication record; Open the cover at other times and trigger the record of opening the cover once; Only records when the phone is connected to the kit will be stored; Bluetooth connect: find the Bluetooth K5102 Pill Dispenser on your mobile phone and connect it; APP control: Download Comfytemp App from the Google Play or App Store and register an account(Supports iOS 12 or Android 5.0). Adjust Bluetooth to "on" status in the "Settings" menu of your mobile phone. Connect the pill organizer in the app and enjoy the app control. 	<h3>Charging the Pill Box</h3> <p>Connect the pill box with a Type-C cable included for charging.</p> <h3>Daily Maintenance</h3> <ul style="list-style-type: none"> When the surface of the pill dispenser is dirty, wipe it with a soft cloth dipped in a little water. During the cleaning process, please prevent water or other liquids from entering the body to avoid short circuit; Do not wipe this product with volatile liquids such as alcohol; If do not use it for a long time, please cover the product to avoid dust. <h3>Product Scrap Disposition</h3> <p>Please support the recycling of rechargeable batteries (lithium-ion batteries)</p> <table border="1"> <tr> <td style="text-align: center;">  Danger </td> <td style="text-align: center;"> When the surface of the pill dispenser is dirty, wipe it with a soft cloth. The battery of this product is a special rechargeable battery which cannot be used in other products. Please do not charge the battery that has been taken out. </td> </tr> </table>	 Danger	When the surface of the pill dispenser is dirty, wipe it with a soft cloth. The battery of this product is a special rechargeable battery which cannot be used in other products. Please do not charge the battery that has been taken out.
Item Name	Smart Pill Organizer																							
Weight	0.33lb																							
Product Size	4*3.5*1.14in																							
Material	Food-grade materials PP, ABS																							
Color	White																							
Compartments	4																							
Input	5V/1A																							
Battery Capacity	1000mAh																							
Package Includes	1*Smart Pill Organizer 1*USB-C Charging Cable 1*User Guide																							
 Danger	When the surface of the pill dispenser is dirty, wipe it with a soft cloth. The battery of this product is a special rechargeable battery which cannot be used in other products. Please do not charge the battery that has been taken out.																							

反面

 Danger <ul style="list-style-type: none"> Do not throw into the fire, make it overheated Do not pierce with sharp objects, gravitational impact, or disassemble or modify Do not connect battery anode and cathode with metal Do not store it with necklaces, hairpins, etc Do not charge, use and place in high temperature places such as near the stove or under the scorching sun to prevent heat generation, fire, rupture, etc Please keep the removed battery out of reach of infants and children to avoid swallow by mistake When the battery fluid leaks, please do not touch it with bare hands, and handle it as follows: When battery fluid enters the eyes, there is a danger of blindness. Do not rub your eyes directly with your hands, rinse eyes with clean water immediately and consult a doctor When battery fluid gets on your body or clothing, it may cause skin inflammation or injury. Please rinse with clean water. If you still feel unwell, consult a doctor immediately 	 Danger <ul style="list-style-type: none"> Only the product can be disassembled when it is discarded. Unauthorized disassembly may cause fire, electric shock, injury and other accidents. If repairs are required (except for battery replacement), please consult the store where you purchased the product <p>The lithium-ion battery used in this product is a valuable resource that can be recycled. Please do not randomly throw. Please dispose them as per local environmental requirements.</p> <p>When disposing of the appliance, the battery must be taken out; When removing the battery, the appliance must be powered off; Battery Removal Steps:</p> <ol style="list-style-type: none"> Use professional tools to pry open the shell on the switch side; Use a screwdriver to unscrew the fixing screws of the internal bracket; Unscrew the screws on the circuit board; Use scissors to cut the battery connection leads; Remove the battery 	<h3>Precautions</h3> <ol style="list-style-type: none"> Inside with a lithium battery, please keep it away from the high-temperature environment. Please keep this product out of reach of children. As a personal product, please do not mix it with others. Please store in a well-ventilated place, keep away from the environment with strong sunlight and corrosive gas and liquid erosion. Please place the product in a safe location, avoid placing heavy objects on the product to prevent damage. 	<h3>Warranty</h3>     <p>If you are not satisfied with our products, please do not hesitate to contact us by sending an email to help@comfytemp.com</p> <p>We offer all customers 5 years warranty, 365 days free replacement, 90 days free return and 24 hours quick customer service.</p>	 www.comfytemp.com <p>Product Name: Smart Pill Organizer Brand: Comfytemp Model: K5102</p> <p>Company Name: UKCA SERVICE LTD Address: 21 Eillesmere Avenue,Worsley, Manchester,UK,M28 0AL E-mail: ukca.service@gmail.com</p> <p>Company Name: TopToll GmbH Address: Einsteinstr.2,Neuss 41464, Germany USt-ID:DE339651736 Email: top tollde@gmail.com</p> <p>Manufacturer: Shenzhen Yicaiwangluo Technology Co.,Ltd Company Name: Xincashenheidianzhangwyuxiangongsi Address: Yousong Second Industrial Zone, Fukang community, Longhua, 416Rm, 4F, 3Bldg, 1980 Technology&Cultural Industry Park Shenzhen GuangDongSheng 518000 CN</p> <p>FCC Statement This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment. This equipment has been tested and found to comply with the limits for a Class B digital device pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on. The user is encouraged to try to correct the interference by one or more of the following methods: •Reorient or relocate the receiving antenna. •Increase the separation between the equipment and receiver. •Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. •Consult the dealer or an experienced radio/TV technician for help. RF Exposure Information The device has been evaluated to meet general RF exposure requirement. This equipment should be installed and operated with minimum distance 0mm between the radiator & your body.</p>    
---	---	--	---	---