


User Manual for Smart Hygiene Technology Indicator

This user manual will introduce how to use urine detector and indicator. Please read before using and keep this user manual.

1. Introduction

Urine Indicator



Item	Description
④ Reset	If indicator is not working normally, please press the reset button
⑤ Cooling hole	Cooling hole
⑥ Stand	Stand
⑦ Speaker	Speaker
⑧ Power Button	Power on/off, adjust brightness, enter pairing mode, mute
⑨ Indicator	White: stand by; Green: dry. Red: Wet and need to change
⑩ Charge Jack	Please use micro-USB cable for charging

2. How to Use

2.1 Urine Indicator

2.1.1 Power button

- Press the button for 2 seconds to power on the device. The light will turn on.
- When, keep pressing for 3 seconds more, it will enter pairing mode for pairing urine detector.

2.1.2 Low battery

- The internal rechargeable battery can only work for 3-5 hours. When the battery exhausted, the device will shut off automatically.
- Plug in USB cable with the charger. It will charge the battery automatically. The cable plugged-in all the time for in house use is recommended.

2.1.3 Brightness


- Press and hold the power button will dim the brightness while the device is in the working mode.
- To brighten the brightness after dim, just press and hold the power button again.

2.1.4 Mute (including wetness alarm, and sleep posture alarm)

- Press the power button once to mute, but it will alarm again if it happens again

2.1.5 Color of urine indicator

- The color will change according to the amount of the urine.



4. Product Specification

Item	Urine Indicator	
Size	R 70 * H 70 mm	
Weight (battery included)	Around 80 g	
Battery	Lithium Polymer Rechargeable	
Transmission	BLE 4.0: Bluetooth Low Energy	Wi-Fi
Distance	15 meters	10-20 meters
Frequency	2.4 GHz	2402MHZ-2480MHZ
Operating Temperature	10 °C – 40 °C	10 °C – 40 °C
Storage Temperature	0 °C – 50 °C	0 °C – 50 °C

FCC WARNING

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

To maintain compliance with FCC's RF Exposure guidelines, This equipment should be installed and operated with minimum distance between 20cm the radiator your body: Use only the supplied antenna.