

Bluetooth body fat scale

Instructions

Body Fat Scale (USB Charging)

Read before use. Bluetooth body fat scale supports IOS and Android versions of two systems

▲ Precautions

- Packaging of weighing scales one month after sale (except for human factors such as humidity, breakage and self-unloading)
- When the numbers are zeroed, the measurement can be started
- Keep away from water sources and fire sources, and do not use chemical cleaners to avoid product failure or damage.
- Please keep it out of the reach of children

▲APP connection instructions

- Download the [OKOK] application on your mobile phone
- The mobile phone turns on Bluetooth, then enters the APP, searches for devices in my device, stands on the scale surface and keeps the screen on, and can link and bind

▲ Operating instructions

- Stand on the scale with both feet smoothly
- Let stand for about 5 seconds and wait for the numbers to stabilize

- measurement data is for reference only

- the weighing range is 0.2-180kg

▲Charging instructions (please ignore this article for battery models)

- Please charge for 1-2 hours for the first time
- The upper limit of charging time is 2 hours
- Please use the 5V1A mobile phone plug. 2A fast charging head is not recommended. When charging, the battery symbol will jump, and stop jumping when it is full.

▲ Weigh items above 0.2kg

(Please press the button on the back first to switch the required units)

- Place the weighing scale on a hard and flat ground
- Apply pressure to the weighing scale to make its screen light up
- After waiting for the number to return to zero, you can display the number by putting the item, and the number of the light item will automatically return to zero after flashing 2 times.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.