LedSonic®

Manual | Warranty

DRAFT

FRA

LedSonic®

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10. Étiquetage

 ESTHER MICHÈLE LedSonic®	
Model No.: QM-EF18 12V  500mA IP22 EN62471:2008 Optical radiation: 0.7W (RED) / 0.4W (BLUE) Wavelength: 660nm (RED) / 415nm (BLUE) Ultrasonic frequencies: 1MHz / 5MHz This device complies with part 15 of the FCC Rules FCC ID: 2A-BNVEF18 Risk group is exempt group according to EN 62471: 2008 Non-continuous operation with red/blue LED bead: 9 minutes on, 8 minutes off Continuous operation with ultrasonic head: 15 minutes operation	
 0120	    

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LedSonic®

4 In 1 Skin Care Therapy

Red, Blue Light Therapy and Ultrasonic in 1MHz/ 5MHz

Instructions for Use

Special Attention To Safety!

Do not focusing on the light!

When using the LEDs attachments, please do not stare directly into the light, otherwise, it can cause eye uncomfortable or injury.

When using the LEDs attachments, do not use on your eyes or eyelids!

When using the LEDs attachments, do not use over the thyroid gland (neck area).

When using the ultrasonic attachments, keep away from your eyes at least 1cm away at all times!

Children should be supervised to ensure that they do not play with the appliance.

Please read all warnings and safety precautions before use, and strictly follow them when using LedSonic®.

Keep away from the children.

ENG

ENG

Table of Contents

1. Introduction
2. Benefits
3. Before You Start
 - 3.1 LedSonic® device description
 - 3.2 Package contents
4. Important Safety Cautions
 - 4.1 Warning
 - 4.2 Cautions
 - 4.2.1 Avoid to use LedSonic® if
 - 4.2.2 Areas not to treat
 - 4.3 Possible side effects
5. Preparing for use
 - 5.1 Charging
 - 5.2 Preparing treatment areas
 - 5.3 Patch test
6. Using LedSonic®
 - 6.1 LedSonic® setup
 - 6.1.1 How to change the attachment
 - 6.1.2 Operate the system
 - A) Ultrasonic therapy mode
 - B) Light therapy mode
 - 6.2 Treatment tips
 - A) Ultrasonic therapy mode
 - B) Light therapy mode
7. Cleaning, maintenance and storage
 - 7.1 Cleaning
 - 7.2 Maintenance and storage
 - 7.3 To protect the environment
8. Troubleshooting
9. Specification
10. Labeling
11. Service and Warranty
12. Copyright Reserved

1. Introduction

Congratulations on your purchase of LedSonic®! Welcome to experience the **4 in 1 therapy** effects of LedSonic®!

How does Ultrasonic works?

The application of ultrasonic frequency is commonly used by salon and clinic professionals but now it is reach of the home user by using LedSonic®. It offers safe and effective facial toning treatment in the comfort of your own home. It can plump up and smoothen the fine lines and wrinkles, skin tightening, shrink enlarged pores, relieve puffy eyes, improve the dark eye circles, deep penetrate the skin care products, stimulation of cellular metabolism, improvement and reduction of scars, freckles, age spots and redness. The process is safe, painless and gentle for skin care treatment.

What is the difference between 1MHz and 5MHz?

Generally, lower ultrasonic frequencies will penetrate deeper below the skin. Thus, the 1MHz mode is designed for tissues and muscles are thicker (i.e. cheeks, chinetc.). The 5MHz ultrasonic mode is appropriate for treatment around the most delicate and thinner area (i.e. under eyes circle, temple, forehead and both sides of the nose)

How does the ultrasonic treatment offers the improvement of skin?

When ultrasound absorbed by muscle, tendons and bone, the harder surfaces like bone will absorb and reflect ultrasound energy. This in turn dramatically increases needed oxygen to the cells and aids in the removal of waste products. This antibacterial cleansing action also helps to improve the acne and blemishes. The skin then only left feeling instantly energized, firm, tighten up, refreshed and noticeably softer – even after just one application. (See Fig.1) The intensive rejuvenating cyclic sound waves penetrate below the skin reaching the dermis. Generally, 1MHz frequency mode will penetrate deeper below the skin. It is designed for face (except under eye circle and around the eyes areas) and body where tissues and muscles are much thinner. The 5MHz frequency mode is appropriate for treatment around the eyes which is the most delicate area. With using the skin care products, the ultrasound can penetrates the ingredients deep into the skin. It can maximize the effects of your skin care products!

ENG

ENG


Fig. 1

How does the light therapy by LedSonic® works?

LED photobiomodulation was originally discovered by NASA about 40 years ago to stimulate plant growth. Specific wavelengths stimulate specific responses in the body. Red light therapy was found to encourage cell repair and healing. Blue light therapy was also found to be effective in treating acne vulgaris. The FDA approved blue light therapy for acne treatment in August 2002. LedSonic® uses a super-luminous LED array to deliver concentrated beams of red light for fractional phototherapy effect. The light therapy by LedSonic® is very gentle and effectively treatment can be performed at home.

How does the Red light therapy works?

Specific wavelengths of red light (660nm) help boost overall cellular energy. These approaches to skin renewal rely on high-powered coherent red LEDs to create heat energy which can result in phototherapy effect to tissues. Cellular activity and energy levels fall as we age. With the delivery of the 660nm light wavelength, the dermal cells receive a rejuvenating, anti-aging boost from the light energy. This safe, painless, natural healing process converts light energy within skin cells, very similar to the way photosynthesis takes sunlight and turns it into food energy in plants.

What does the Red light therapy by LedSonic® benefits to the skin?

LedSonic® uses a super-luminous LED array to deliver concentrated beams of red light for fractional phototherapy effect.

The phototherapy enhances blood circulation, helping bring nutrients and cellular growth factors to the treated areas. It stimulates fibroblast and collagen production to strengthen the dull and wrinkled skin. These cellular and metabolic effects can help to reduce the localized hyperpigmentation, improve neocollagenesis to fight against skin aging and improve the skin tone, increase in tissue oxygenation for the clarification of dull skin.

Red light therapy is safe for all skin types and ages, there is no down time or recovery time.

It is one of the most effective and non-invasive ways to help rejuvenate and maintain the skin.

How does the Blue light therapy works?

Specific wavelengths of blue light (415nm) penetrate the skin and which effectively destroys acne-causing bacteria which is approved by FDA. It is showing the reductions in inflammatory acne lesions and rendering phototherapy with blue light an effective safe treatment for acne.

What does the Blue light therapy by LedSonic® benefits to the skin?

LedSonic® uses a super-luminous LED array to deliver concentrated beams of blue light for fractional phototherapy effect.

Once delivered, this light energy is absorbed by the targeted bacteria that produces the inflammation associated with acne known as Propionibacterium Acnes, or "P.acnes". P.acnes releases porphyrins, which are naturally occurring molecules in the body. When porphyrins absorb the LedSonic® 415nm, free radical damage is produced which destroys the bacteria.

Without P.acnes around to cause inflammation, acne clears up. (See Fig.2)

Note: The results of the above treatment can be vary by individual and skin type. Continued daily application will provide more long-term and lasting results.


Fig. 2

ENG

ENG

2. Benefits

LedSonic® provides the following benefits :

Experience for your own treatment schedule

With the ultrasonic technology and phototherapy, you are now able to create the most suitable treatment schedule for yourself. The condition of skin is vary by weather, physical period change and stress.....etc. You are the only one who know yourself the most. Now, you can be your best therapist!

Advanced and Innovative 4 in1 technology for home use

Ultrasonic frequencies in 1MHz/ 5MHz, Red light and Blue light therapy are used in the professional beauty salons and clinics for skin rejuvenation. But now, you can also use the LedSonic® which brings this innovative technology in a single device to the comfort of your home.

LedSonic® provides you a total skin care solution:

Mode	Treatment
Ultrasonic - 1 MHz (2 Level)	For cheeks, chin, the lower jaw (See Fig.3 "U1")
5 MHz (2 Level)	For under eyes and around eyes, forehead, temple, neck (See Fig.3 "U2")
Red light	Skin rejuvenation, tighten up and ingredients penetration (See Fig.3)
Blue light	Acne Cleaning (See Fig.3)

U2-5MHz

U1-1MHz

Red light


Blue light

Fig. 3

Unique pattern design to maximize the effectiveness

With the unique pattern design for the surface of the ultrasonic attachment, it allows to hold effectively the ingredients of your skin care products in the contact with the skin. (See Fig.4)

High intensity energy

High intensity 60mW/cm² illuminated surface produce 6 joules/cm² dose in 120 seconds.

Design for convenient and flexibility

With the light weight, portable and rechargeable design, you can use LedSonic® at home anywhere and anytime!

Gentle treatment

During the process you will only feel a little heat. It does not contain any harmful ultraviolet (UV) light which can damage skin and does not result in dryness or discomfort.


ENG

ENG

3. Before You Start

3.1. LedSonic® device description

LedSonic® consists of the applicator, attachments for ultrasonic/ red LEDs/ blue LEDs, charger base and an adaptor:

General Description (See Fig. 5)

1. Main Body
2. Attachment connector base
3. Touch key for "Power On"
4. Touch key for "Select" the ultrasonic mode or light on
5. Ultrasonic mode indicator
6. Charge contact point
7. Charger Base
8. Red LEDs attachment
9. Blue LEDs attachment
10. 1MHz/ 5MHz Ultrasonic attachment
11. Power adaptor/Vents

3.2 Package contents

- Applicator
- Attachment for Ultrasonic/ red LEDs/ blue LEDs
- Charger Base
- Adaptor
- This instructions manual
- Soft Cloth

4. Important Safety Cautions

 Please read this instructions manual carefully before you use the appliance and save it for future reference. If you are unsure if LedSonic® is safe for you to use, please consult your doctor or dermatologist!

4.1. **WARNING:**

As with most electrical appliances, electrical parts are electrically live even when the switch is off. To reduce the risk of injury or death by electric shock:

- Do not reach for the unit if it has fallen into water or other liquid. Unplug unit immediately.
- Do not place or drop into water or other liquid.
- The adaptor and charger base are not waterproof, keep away from the water
- Keep the appliance and the adaptor dry at all times
- Do not attempt to open or repair the appliance by yourself.
- If the appliance is broken, do not touch any inner part to avoid electric shock.
- Never insert any objects into the appliance.
- Do not use the appliance if it is not working properly or if it appears damaged.
- Do not use near the shortwave or microwave equipment
- Make sure the adaptor do not get wet.
- Keep away from flammable aerosol products being used or where oxygen is being administered.
- Only use the parts supplied with the appliance.
- Only charge the appliance with the adaptor supplied.
- Do not use the appliance or the adaptor if it is damaged.
- Do not charge when exposed to direct sunlight or heat
- Keep this product out of reach of children. The use of this appliance by persons with reduced physical, sensory or metal capabilities or lack of experience and knowledge can give cause to hazards. Unless they have been given supervision or instruction concerning the use of appliance by a person responsible for their safety.
- The adaptor contains a transformer. Do not cut off the adaptor to replace it with another plug as this causes a hazardous situation.
- Turn off the appliance if it becomes unusually warm to touch..
- Switch off the appliance before cleaning
- Stop immediately if you experience any adverse effects or develop prolonged redness. You may be experiencing a photosensitive reaction.
- LedSonic® is designed for occasional home use (as opposed to exhaustive commercial or professional use). Avoid such exhaustive use.
- LedSonic® is only designed for personal use, **DO NOT** share with others.

ENG

ENG

4.2 **Cautions**

 Never use the appliance with the LEDs attachments around and on the eyes to avoid the risk of eye damage.

 Read all instructions before using LedSonic®, use this appliance only at settings suitable for your skin type. LedSonic® is not designed for everyone. Certain conditions may prevent the use of LedSonic®.

DO NOT USE LedSonic® if any of the following currently apply to you.

 Please note that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

 This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.

 * **Do not use** a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.

 **Caution:** This unit has been thoroughly tested and inspected to assure proper performance and operation!

 * **Caution:** this machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

-  Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance $d = 3,3$ m away from the equipment.

Note: This list is not exhaustive. If you are not sure whether you can use the appliance, please consult to a doctor.

-  This equipment has been tested and found to comply with Part 18 of the FCC Rules. Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

-  **CAUTION:** Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the product.

Note: This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna. Increase the separation between the equipment and receiver. Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. Consult the dealer or an experienced radio/TV technician for help."

4.2.1 Avoid to use LedSonic® if

- Do not use if you are pregnant, breast feeding or nursing
- If you have a pacemaker or internal defibrillator
- If you have superficial metal or other implants in the desired treatment area
- Do not use on dry or fragile skin caused by the use of chemical peels, glycolic peels, Alpha Hydroxy Acids (AHAs), Beta Hydroxy Acids (BHAs), Retin-A, Accutane and topical isotretinoin and azelaic acid.
- If you have a history of any type of cancer, or pre-malignant moles.
- If you have varicose veins in the treatment area
- If you have pre-cancerous lesions or multiple atypical moles in the areas to be treated.
- If you have dark brown or black spots, such as moles, large freckles, birthmarks or warts on the area to be treated

ENG

ENG

- If you have eczema, psoriasis, lesions, open wounds or active infections, such as cold sores in the area to be treated. Wait for the affected area to heal before using LedSonic®.
- If you have an impaired immune system due to immunosuppressive diseases such as AIDS and HIV, or use of immunosuppressive medications
- If your skin is sensitive to light and easily develops a rash or an allergic reaction
- If you have had surgery in the areas to be treated
- If you epilepsy with light sensitivity
- If you have diabetes, or other systemic or metabolic diseases or congestive heart disease
- If you have any other condition which in your physician's opinion would make it unsafe for you to be treated
- If you have Epilepsy
- If you are on gold therapy
- If you have received radiation therapy or chemotherapy treatments within the past 3 months

Note: This list is not exhaustive. If you are not sure whether you can use the appliance, please consult to a doctor.

4.2.2 Areas not to treat

- Do not use on eyes and should keep away from eyes for at least 1cm
- Do not use on the middle of neck (throat), it is for ultrasonic mode only.
- Do not use on nipples, areola, or genitals
- Do not use on dark brown or black spots such as moles, birth marks or freckles
- Do not use on an area of recent surgery, deep peel, laser resurfacing, scars, or skin that has been damaged with burns or scalds
- Remove the applicator from the skin if it is unusual warm

4.3 Possible side effects

When used according to the instructions, side effects and complications associated with use of LedSonic® are rare. However, every energy based systems may involve some degree of risk. Therefore, it is important that you understand and accept the risks and complications that can occur. LedSonic® is a gentle home use appliance, most users only sense the moderate heat applied by the device during use with no further discomfort.

- Skin may become slightly irritated after using LedSonic®. This is a typical and normal phenomenon. If the irritation does not disappear within 3 day, please consult a physician or a dermatologist.
- Skin may become slightly red. This should disappear a few hours after treatment.

5. Preparing for use

5.1 Charging

- Fully charge the batteries before you use LedSonic® for the first time and when the batteries are empty. Fully charging the batteries takes up to 5 hours
- Low battery alert: Charging LedSonic® when hearing the continuous "D..D.." sound during use to indicate that the battery is low.
- Fully charged batteries provide operation time at least approx. 1 hr 20 min. for ultrasonic mode and approx. 1hr. 40min. for LEDs light therapy mode.
- Fully charge the appliance every 3 to 4 months, even if you do not use LedSonic® for a longer time.

Directions for charging the LedSonic®:

1. Switch off the appliance
2. Insert the small plug from the adaptor into the back of the charger base and put the adaptor in the wall socket. (See Fig. 6)


Fig. 6

Notes: The LEDs attachments into the charger base are only for storage purpose.

ENG

ENG

3. Charging indicator light (See Fig. 7)

- The "Select" key is blinking in blue while charging
- The "Select" key is light on in blue when fully charged


Fig. 7

Notes: It is normal when the adaptor and the appliance feel warm during charging.

Warning: NEVER cover the appliance and adaptor during charging

4. Remove the adaptor from the wall socket and pull the small plug out of the appliance when fully or after charged.

5.2 Preparing treatment areas

- Start by cleaning your face with non-oil based cleanser and free of any powder on the skin.
- Use the provided soft cloth to clean the surface of the attachment.

5.3 Patch test

- Familiarize yourself with the features of the LedSonic® device
- Turn on the LedSonic® with the ultrasonic mode and test on a small patch of skin and wait 24 hours.
- Light Test : Turn on the LedSonic® with the Red LEDs attachment (and test again in Blue LEDs attachment) and test on a small patch of skin with 2 minute treatment on the same area. Wait 24 hours
- The patch test is to ensure there are no adverse reactions. If irritation occurs do not continue treatment and please consult the doctor or dermatologist.

6. Using LedSonic®

6.1 LedSonic® setup

6.1.1 How to change the attachment?

Press the buttons on the left and right sides and place the selected attachment onto the attachment connector base. (See Fig. 8)

6.1.2 Operate the system

(Please check with section 3.1

"LedSonic® device description")

A) Ultrasonic Therapy Mode

- Step 1: Ultrasonic attachment is well equipped on the attachment connector base
- Step 2: Long touch "ON/OFF" for power on.

Note: The buttons are touch sensors, please use a dry finger while touching.

- Step 3: Touch "SELECT" one time for selecting difference frequency mode (1MHz/5MHz) and energy level (Normal/ Strong) for difference purposes. (See. Fig. 9 for the sequences of the default setting when press "SELECT")

Fig. 9

Note: For details, please see section 7.2 "Treatment Tips"

- Step 4: To change the frequency mode or energy level, touch "SELECT" again. The light indicator (See Fig. 6) is blinking in difference status in each selection :

Frequency Mode	Level	Indicator Status
1 MHz	Normal	"1MHz" Blinking slowly
	Strong	"1MHz" Blinking quickly
5 MHz	Normal	"5MHz" Blinking slowly
	Strong	"5MHz" Blinking quickly


Fig. 8

ENG

ENG

- Step 5: Long touch "ON/OFF" for power off when finished treatment.
- Step 6: Use a damped soft cloth to clean the surface of the attachment.
- Step 7: The appliance will stop function after 15 mins for reminding the treatment time, you have to press "SELECT" again to continue the treatment.

Note: As it is recommended to do 15 mins for each treatment, the appliance will be power off after 15 mins automatically.

B) Light Therapy Mode

- Step 1: Red/ Blue LEDs attachment is well equipped on the attachment connector base.
- Step 2: Long touch "ON/OFF" for power on.
- Step 3: Touch "SELECT" for light on.
- Step 4: Long touch "ON/OFF" for power off when finished treatment.
- Step 5: Use a damped soft cloth to clean the surface of the attachment.
- Step 6: LedSonic® will sound "b" for each one minute to remind the treatment time.
- Step 7: LedSonic® will stop function after 9 mins for reminding the treatment time, you have to press "SELECT" again to continue the treatment.

Note: The appliance will be power off after 9 mins automatically and please wait 8 mins for cool down to power on again.

6.2 Treatment Tips

A) Ultrasonic Therapy Mode

- Apply the skin care products on the ultrasonic attachment surface (See Fig. 10) and suitable effective frequencies area (See. Fig. 3)
- Apply on a serum/ lotion/ aqueous skin care product. The energy level is recommended to use in difference structure of the skin care product to maximize the penetration to the skin :

Fig. 10

Fig. 11

B) Light Therapy Mode

LedSonic® uses a super-luminous LED array to deliver concentrated beams of light for fractional phototherapy effect. This safe, painless, natural healing process converts light energy within skin cells, very similar to the way photosynthesis takes sunlight and turns it into food energy in plants. To maximize the effect, we highly recommend you:

- For acne treatment: Use the blue LEDs light therapy for 3 minutes over affected areas twice daily. To maximize the results, follow with the red LEDs light therapy use for 3 minutes over affected areas twice daily.
- Skin rejuvenation and firming: Use the red LEDs light therapy, holding the appliance for 3 minute on difference of each area of the face (see section 4.2.2 "Areas not to treat"), then move to another area of the face. The appliance will sound "B" to remind you in every minute.
- Use it for 3 times a week.
- If the area that require more rejuvenation, it is suggested to use for 3 minutes over there.

7. Cleaning, maintenance and storage

7.1 Cleaning

⚠ CAUTION: Before cleaning your LedSonic®, make sure that the power switch is OFF and the power cord is disconnected from the charger base.

- Never clean LedSonic® or any of its parts under the tap or in the dishwasher.
- Never use abrasive cleaning agents or aggressive liquids such as petrol or acetone to clean the appliance.
- Use the soft cloth included with your LedSonic® appliance to clean the following parts:

- Maximize the effects of your skin care products – ultrasonic wave is effective for the ingredients penetration, and specific frequencies for difference area:
A) Using with products like essence serum, lotion: Repeat each step a total of 5 times glide LedSonic® very slowly and slightly. (See Fig. 11)
B) Using with products like facial mask: Repeat each step a total of 10 times glide LedSonic® very slowly and slightly. (See Fig. 11)

ENG

ENG

- The surfaces of the main body and charger base
- The surfaces of the attachments

7.2 Maintenance and storage

- Store LedSonic® in a dry place at a temperature between 0°C and 60°C
- If you want to replace a part or want to purchase an additional part, please go to your local provider or contact us at www.esthermichele.com

7.3 To protect the environment

- At the end of its life, **do not** discard LedSonic® with the normal household waste. Dispose of it at an official recycling collection point.

8. Troubleshooting

Refer to this troubleshooting guide if you experience any problems with LedSonic®, as this section addresses the most common problems you could encounter with this appliance. If you still have problems, please contact go to your local provider or contact us at www.esthermichele.com


Problem	Possible reason	Solution
LedSonic® does not work	The rechargeable batteries are empty The appliance is defective Improperly touch the "ON / OFF" key	Charge the appliance (see section 5.1 "Charging") Contact the nearest EM dealer Use a dry finger to touch on the key for few seconds.
LedSonic® suddenly switches off	The rechargeable batteries are empty	Charge the appliance
The appliance temperature rises during operation and recharge	This is normal	No action is required
Battery charging indicator light ("SELECT" key blue light) does not turn on when connecting the PLUG into the CHARGER BASE SOCKET	The POWER ADAPTOR is not plugged in the electrical outlet	Insert the PLUG into the CHARGER BASE SOCKET. Connect another wall socket to check if the electrical outlet is live If it is still does not charge, contact your nearest EM dealer

9. Specification

Ultrasonic Frequencies		1MHz/ 5MHz
LEDs Light wavelength	RED	660nm
	BLUE	415nm
LEDs Light Energy	RED	50mW/ cm ²
	BLUE	30mW/ cm ²
Treatment area		Approx. 11.5cm ²
Lithium-ion battery		11.1V / 750mAh
Charging time		Approx. 6 hours
Operating time	Ultrasonic Therapy	1 hr. 20 mins
	LEDs Light Therapy	1 hr. 40 mins
Rated voltage		100V-240V
Rated Frequency		50Hz-60Hz
Adaptor spec.		AC/DC adapter (S08-006-0120-00400): Input: 100-240V, 50/60Hz, 0.2A; Output: 12Vd.c., 400mA
Operating conditions		Temp. +15°C to +30°C
		Relative humidity: 25% to 75%
		Barometric pressure: 800-1060hPa
Storage and transport conditions		Temp. 0°C to +60°C
		Relative humidity: 10% to 90%
		Barometric pressure: 800-1060hPa

10. Labeling

ENG

ENG

Symbol	Description
	Warning/ Cautions
	Keep away near the eyes
	Class II Equipment
	Keep away from the water
	CE Mark
	FCC Mark
	WEEE Waste Electrical and Electronic Equipment
	Refer to instruction manual
	Type BF equipment

Safety Standards

- Medical CE: IEC 62471/EN 62471,
IEC 60601-2-57:2011/EN 60601-2-57:2011,
IEC 60601-1:2005/EN 60601-1:2006,
IEC 60601-1-2:2007/EN 60601-1-2:2007,
IEC 60601-2-5:2000/EN 60601-2-5:2000,
IEC 60601-1-11:2010/EN 60601-1-11:2010,
- FDA: ISO10993-5:2009, ISO10993-10:2010
- EMC: EN55011:2009+A1:2010, EN61000-6-1:2007,
EN61000-3-2:2006+A1:2009+A2:2009, EN61000-3-3:2008
- FCC CFR47 PART15 SUBPART B:2012
- Manufacturer System: ISO13485:2003 and 93/42/EEC

11. Service and Warranty

We warrant this product against any defects that are due to the faulty material or workmanship for a 2 year period from the original date of consumer purchase.

If the product should become defective within the warranty period, we will repair any such defect without charge provided there is proof of purchase. This does not mean an extension of warranty period.

The warranty shall apply in all countries in which our product was sold via an authorized dealer.

This warranty does not cover the damage is due to by accident or misuse, abuse, alteration to the product or use inconsistent with the technical and /or safety instructions required.

This warranty shall not apply if the product has been dismantled or repaired by a person not authorized by us.

12. Copyright reserved

Copyright © 2012 Esther Michèle SAS
All rights reserved.

Esther Michèle reserves the right to make changes to its products or specifications to improve performance, reliability, or manufacturability. No part of this document may be produced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Esther Michèle SAS. Esther Michèle SAS has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this documents. System design and specifications may change without prior notification.

Contact information for customer service:

Esther Michèle SAS

Website: www.esthermichele.com

E-mail: info@esthermichele.com

ENG

ENG

Certificate of Purchase/Certificat d'achat

ESTHER MICHÈLE

Mode Name/Nom de Mode: LedSonic®

Date of Purchase/
Date d'achat: _____

Dealer stamp/
Cachet du revendeur