EN READ ALL INSTRUCTIONS BEFORE USING THE KOLIBREE TOOTHBRUSH

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
EN BEFORE USING THIS PRODUCT, PLEASE READ THESE INSTRUCTIONS CAREFULLY AND UNDERSTAND THE SAFETY INSTRUCTIONS.

1. This product is not serviceable. To avoid electric shock, do not attempt to open the outer product covering.
2. This product is for indoor use only, do not expose to liquid, vapor, or rain
3. CAUTION-Risk of Electric Shock, disconnect power to the receptacle before installing or removing the unit. When removing receptacle cover screw, cover may fall across plug pins or receptacle may become displaced.
4. This power unit is intended to be correctly orientated in a vertical or floor mount position.

EN CORRECT DISPOSAL OF THIS PRODUCT

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection system or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

Created by Kolibree in France, Made in the Philippines.
©2016 All rights reserved. iPhone, iPod and iPad are registered trademarks of Apple Inc, registered in the US and other countries. Android is a trademark of Google Inc.

EN READ ALL INSTRUCTIONS BEFORE USING THE KOLIBREE TOOTHBRUSH

DANGER– To reduce the risk of electrocution:

1. Do not use while bathing.
2. Do not place or store product where it can fall or be pulled into a tub or sink.
3. Do not place in or drop into water or other liquid.
4. Do not reach for a product that has fallen into water. Unplug immediately.

WARNING– To reduce the risk of burns, electrocution, or injury to persons.

1. Close supervision is necessary when this product is used by, on, or near children or invalids.
2. Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
3. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to a service center for examination and repair.
4. Keep the cord away from heated surfaces.
5. Never place it on a soft surface, such as a bed or couch, where the air openings may be blocked.
6. Never use while sleeping or drowsy.
7. Never drop or insert any object into any opening or hose.
8. Do not use outdoors.

WARNING– For the purposes of recharging the battery, only use the detachable supply unit provided with this appliance. Model No KLTB002 supply unit.

1. ⚠️: Use only with model No. KLTB002 supply unit.
2. The battery must be removed from the appliance before it is scrapped.
3. The appliance must be disconnected from the supply mains when removing the battery.
4. The battery is to be disposed of safely.

IMPORTANT SAFEGUARDS

Toothbrushes can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Cleaning and user maintenance shall not be made by children without supervision. The appliance is only to be used with the power supply unit provided with the appliance. The external exible cable or cord of this adaptor cannot be replaced; if the cord is damaged, the adaptor shall be scrapped. The socket-outlet should be installed near the equipment and be easily accessible.

WARRANTY

You can view the warranty of your purchased product by visiting our website at: www.kolibree.com/warranties.

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
"This device complies with FCC radiation exposure limits set forth for general population (uncontrolled exposure). This device must not be collocated or operating in conjunction with any other antenna or transmitter."

The device must not be co-located or operating in conjunction with any other antenna or transmitter.

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
FCC ID: 2ACPT-KOLIBREE16

EN FCC STATEMENT

WARNING: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. Consult the dealer or an experienced radio/TV technician for help. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures: Reorient or relocate the receiving antenna. Increase the separation between the equipment and receiver. Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. Consult the dealer or an experienced radio/TV technician for help. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

IC STATEMENT

This device complies with Industry Canada licence- exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device. Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada. To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p.) is not more than that necessary for successful communication.

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- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Created by Kolibree in France, Made in the Philippines,
©2016 All rights reserved, iPhone, iPad and iPad are registered trademarks of Apple Inc.,
registered in the US and other countries. Android is a trademark of Google Inc.