MATP0816 V3

Additional information and multilingual versions: scan the QR Code

How to use the automatic fall detector Vibby Oak S from Telecom Design.

The Vibby Oak S is configured in wrist wearing or pendant mode* by your installer. You should not change the wearing mode without prior notice to your installer.

To ensure your effective protection, it needs to be worn permanently. In case of high sensitivity and increased risk of allergy, wear the Vibby Oak S in pendant mode.

FUNCTIONS

The Vibby Oak S is an information device in the event of a heavy fall or when pressing the SOS button. Its function is to transmit an alert to a call centre via a telecare device installed in the house.

Manual alarm:

In the event of an accident or if the user feel in danger, you can manually raise an emergency call to your telecare provider by pressing the SOS button. It is located in the middle of the Vibby Oak S and identified by tactile raised dots. The device flashes and vibrates for a few seconds. It is then no longer possible to cancel the alarm.

Automatic alarm:

Following a heavy fall** the Vibby Oak S checks if the wearer is unable to stand up by analyzing movement for 30 seconds. If the fall is confirmed, the alarm is sent to the telecare provider via your telecare device. Your Vibby Oak S alerts you that it is sending an alarm by vibrating intermittently and flashing for 30 seconds. To manually cancel the alarm, you must cover the Vibby Oak S with your palm hand and then remove it as soon as the vibrations stop.

The Vibby Oak S does not detect all falls. In the event of a fall, if you are able to do so, you should always press the SOS button to ensure the alert is sent.

* Based on the offer.

** A heavy fall is a sudden fall occurring after a period of normal activity in an upright position. It is followed by a strong impact with the floor characterized by a period of immobility on the floor.

TECHNICAL SPECIFICATIONS

Model: VB-WTC-10
Weight: 35g
Waterproof: IP67

Shock resistance: IK04
Cleaning: a damp tissue***
Plastics processing: hypoallergenic

| | Radio frequency |
|------------|-------------------|
| SRD | 902 – 928 MHz |
| BLE | 2400 – 2483.5 MHz |

Radio range: 300m in free field
Operating temperature: +41°F / +104°F (+5°C / +40°C)
Battery life: Up to 3 years. The battery must be changed by your installer (CR2477 / LiMnO2)

Storage recommendations:

- **Temperature:** +41°F / +104°F (+5°C / +40°C)
- **Humidity:** 80% maximum
- **Maximum duration:** 6 months

***No household cleaning product

IMPORTANT: The device's technology does not allow for the analysis and interpretation of all situations and the detection of all falls. Some falls, such as soft falls, controlled falls against a wall or chair, or some heavy falls, are not detected by the Vibby Oak S. Exceptionally, certain fall patterns and/or movements of the bracelet or pendant may be analyzed as a heavy fall even though they do not meet the heavy fall scenario. If a false alarm is triggered, quickly clarify the situation. For more information on automatic fall detection or the operating of the product, please consult: <https://telecomdesign.com/downloads/>

REGULATIONS AND CERTIFICATIONS

All our official documents and reports (radio, hypoallergenic, REACH, RoHS, etc.) are available on <https://telecomdesign.com/downloads/>

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

NOTE: The grantee is not responsible for any changes or modifications not expressly approved by the party responsible for compliance. Such modifications could void the user's authority to operate the equipment.

Declaration of conformity

We hereby certify that the Vibby Oak S designed by Telecom Design, complies with the main requirements and other associated regulations presented in Directives 2011/65/EU (RoHS). A copy of the declaration of compliance is available for consultation at <https://telecomdesign.com/downloads/>

Address of the Manufacturer

The Vibby Oak S is a product manufactured by Telecom Design Zone Actipolis II, 2bis rue Nully de Harcourt, 33610 Canéjan, France

Warning: Battery

- Do not ingest battery, Chemical Burn Hazard
- This product contains a coin battery. If the coin battery is swallowed, it can cause severe internal burns in just 2 hours and can lead to death.
- Keep new and used batteries away from children. If the battery compartment does not close securely, stop using the product and keep it away from children. If you think batteries might have been swallowed or placed inside any part of the body, seek immediate medical attention.
- Risk of explosion if the battery is replaced by an incorrect type.
- Disposal of a battery into fire or a hot oven, or mechanically crushing or cutting of a battery, that can result in an explosion.
- Leaving a battery in an extremely high temperature surrounding environment can result in an explosion or the leakage of flammable liquids or gas.

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