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①智能变色刷毛
(Intelligent color bristles)

②刷头
(Brush head)

③开关
(Switch)

④电池盖
(Battery cover)

⑤换刷头
(Replace brush head)

⑥换电池
(Replace battery)

顺时针
刷头安装
Clockwise install

逆时针
刷头取下
Counterclockwise
take down

按住电池盖一侧向外推
Press the side of the battery cover to push out

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保修说明（Warranty Description）：
1、本产品保修期为一年，在保修期内，凡属于正常使用情况下由于产品本身质量问题引起的故障，本公司将负责给予免费维修。
This product is guaranteed to keep in good repair for one year. During this period, any belongs to under normal usage circumstance cause because of the product quality' s problem of breakdown, our company will be responsible for giving free maintain.
2、非本公司产品质量原因引起的维修，如使用不当、保管不妥、擅自拆机等原因造成的损失，公司维修需要另外收取费用。
The damaged caused by other reasons such as don' t follow the request of the Manual or persevere inappropriate, the maintain needs extra charge
3、保修卡及购机发票一经涂改，保修即时失效
Once the Warranty Card or purchasing invoice has been changed, the warranty immediately expired
4、保存保修卡及购机发票作为产品的保修凭证，请用户妥善保管，遗失不补
This card and the purchasing invoice are both considered as Protection to fix warrantees, so please reserve them carefully. Lose don' t repair.
5、服务热线（Customer service hotline）：400-689-9282
6、本服务承诺最终解释权归深圳嘀嗒互动所有
The service commitment to the final interpretation of Shenzhen Tita Interactive Technology Co., Ltd. all

Welcome to Thomas!
This is a sonic toothbrush that connects to an interesting brushing mobile device application via Bluetooth,which can draw the attention of the child while developing a good habit of brushing your teeth on time.
Before operating this toothbrush, please read these instructions and save this manual for future reference.

Important

1.Periodically check the battery cover is securely fastened. A damaged or non-functioning unit should no longer be used. Do not modify or repair the product. This may cause fire, electric shock or injury.
2.If you don' t do not intend to use this product for a long period, the battery should be removed.
3.Cleaning and maintenance should not be performed by children.
4.Use this product only for its intended use as described in this manual. Do not use attachments which are not recommended by the manufacturer.
5.Replace brush head every 3 months or sooner if brush head becomes loose due to damage.

Warning

1.Do not submerge the unit in water or liquid.
2.When replacing batteries pay attention that water does not come in contact with the surrounding area.
3.This toothbrush is a child' s care device and is not intended for use on more than one children.

Replacement

1.The toothbrush adopts intelligent hierarchical color bristles,more than the usetime,brush the top color will have significant change,please change in time.
2.Replace Thomas brush heads every 3 months to achieve optimal results.

Using your toothbrush

1.Scan the QR code below or download the “i brush” application from the App Store or Google Play

2.Place the brush head firmly onto the toothbrush shaft.
3.Install new battery paying close attention to the battery orientation.

Sweet tips

1.The toothbrush has been cleaned before leaving the factory but it is recommended you clean again before use.
2.Product contains batteries and/or recyclable electric waste. For environment protection do not dispose with household waste, rather take it to the electric waste recycling collection points provided in your country.
3.Your device must support Bluetooth 4.0, otherwise you cannot connect to the App. Proposed system version: Later than Android 4.0 or IOS 7.0
4.App functions maybe changed due to software upgrades. Please go to www.titakid.com for latest information.

Cleaning recommendations

Turn off the main unit and remove the brush head before cleaning them with water
Brush head:
1.Carefully wash the inside and the bottom areas of the brush head.
2.Shake well to remove water after washing.
3.Do not pull the bristles of the brush.
4.Do not wash with water hotter than 80°C.
Main unit:
1.Keep the cover of battery away from water when washing.
2.Do not immerse the main unit under water.
3.Do not wash with water hotter than 40°C.

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FCC Warning

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions

(1) this device may not cause harmful interference, and

(2) this device must accept any interference received, including interference that may cause undesired operation. Changes or modifications not expressly approved by the party responsible for compliance could void the user’s authority to operate the equipment.

Note 1: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

—Reorient or relocate the receiving antenna.

—Increase the separation between the equipment and receiver.

—Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

—Consult the dealer or an experienced radio/TV technician for help.

Note 2: 1.Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user’s authority to operate the equipment.