Piezo Instructions

Set the PiezoRxD Clock

It is important your PiezoRx is set to the correct time. Your PiezoRx uses the time function to automatically reset your daily data to '0' at midnight every day.

Remember to adjust your time if Daylight Savings Time *Daylight Savings Time affects your region.*

In your step display (the shoe will be in the top left of the screen) hold

+/- until the screen stops flashing

Press **+/-** to set the hour (AM & PM)

Press **MODE** to display min

Press **+/-** to adjust min

Press to finish

If you don't wear a watch... no problem.

With the PiezoRx you always have the time!

If you are travelling in another time zone, the PiezoRx uses your home timezone.

Clear any Data & Remove the Protective Film

Your Piezo may have gathered some "data" during shipping.

To clear all of the data and memory, press & hold the large button in the bottom right corner for 15 seconds.

When you receive your Piezo, you can remove the protective film on the LCD.

Position your PiezoRx

Slide one corner of the PiezoRx clip onto your waist band or belt, and position it so that it is secure and level. Be sure to attach the loss prevention strap by clipping it to a belt loop or to your clothing.

*The Piezo can also be worn on the bottom band of a bra.

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Icons & Main Displays

Steps, Your Eventual Goal, 8,000 a day

Check your steps throughout the day... when you can, make choices that will increase your number!

* Averaging 8,000 steps a day? Great! Now it's time to start tracking your MVPA...

Press for Next Display

MVPA

Moderate/Vigorous Time

Your Eventual Goal?

150 min/wk

Your total minutes of physical activity at a moderate intensity (or higher) for the past 7 days from today.

Minutes of MVPA is strongly associated with your health. Learn more...

* Averaging 150 minutes a week? Great! Now it's time to start tracking your BOUTS...

BOUTS

Your Eventual Goal?

3 most days, 15 /wk

Each star represents a 10 minute BOUT of continuous activity time, at or above a moderate level (100 Steps/Minute)

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Additional activity Time Displays

Moderate + Vigorous Time Today

(Hrs:Min:Sec)

Your total minutes of physical activity at or above moderate intensity (100 steps per minute) for TODAY.

Moderate Activity Time Today

(Hrs:Min:Sec)

Your total minutes of physical activity at or above a moderate intensity (>100 steps per minute) & below vigorous intensity (<120 steps per minute) for TODAY

Vigorous Activity Time Today

(Hrs:Min:Sec)

Your total minutes of physical activity at or above vigorous intensity (>120 steps per minute) for TODAY.

Total Activity Time Today

(Hrs:Min:Sec)

Your total minutes of all physical activity at or above 60 steps per minute for TODAY.

(your awake hours) - (T) = your sedentary time

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Viewing your 33 days of display memory

You can access the memory of all your data *except MVPA-7 (which is an average of the last 7 days)*.

Simply press the **M** button on any screen to access up to the last 33 days of data!

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4 M = 4 Days Ago Total

... etc up to 33 Days Ago

Press the large square button to return to your last display.

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Adjusting the SPM thresholds

The SPM (Steps per Minute) thresholds is the pace you need to step to earn moderate/vigorous physical activity time. Pre-set SPM are based on research & apply to the average individual. Moderate SPM = 100

Vigorous SPM = 120.

The SPM thresholds are sometimes changed based on height, health status and fitness levels that are lower or higher than average.

Adjusting your SPM thresholds

On your or display, hold +/- x4 sec, then press again to adjust.

Changing your 3V battery

When the battery icon shows under the shoe, place a coin in the slot along the bottom & twist to pop off the cover.

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Synchronizing your PiezoRxD

Activate your LogYourSteps account and enter the PiezoRxD digital ID on your account page.

Insure your computer has the PiezoRxD Data link software installed.
(Download from the LogYourSteps home page)

Place your PiezoRxD within 6 inches of your Apple desk top or laptop computer.

Press and hold the Bluetooth button on the far right of your PiezoRxD for 3 seconds

Follow the prompts on the screen to sync your data.

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Follow the prompts on the screen to sync your data.

StepsCount Inc.

PO Box 430
Deep River, Ontario
K0J1P0
866.342.2328

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions:

(1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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This device complies with FCC Part 15 of the FCC Rules.
Operation is subject to the following two conditions:

(1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Caution:

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.