User Manual

Warning

Due to the unaccustomed operation some gum bleeding may initially occur, this is normal. Should this continue for more than 2 week, please stop using the appliance and contact your dentist.

Do not permit liquids of any kind to enter the charging station;

After cleaning, ensure the charger is completely dry;

This appliance contains rechargeable battery, please do not put it into fire or placed in heated surface;

This appliance cannot be knocked or dropped, if the appliance is damaged, stop using it.;

Do not repair or disassemble the product by youself;

Stay away from strong electromagnetic field when charging.

Danger

Check before you use it, do not use if the appliance is damaged;

Before you connect the power supply, please check the local voltage is consistent with the rated voltage;

If the appliance is damaged in any way (brush head,toothbrush handle,etc), stop using it;

Do not use the charger outdoor or near heated surface;

Do not use in the bath or the shower;

Notices

Do not share the brush head with others:

Do not use the power cord around on the adapter;

Do not put the metal objects on the charger, such as coins or clip, etc;

Nnplug the power supply to avoid fire or danger ,if you do not use for a long time.

Product description

The power shaft protection cap

The brush head protection cover

Brush head

handle

Power & mode switch button

The bluetooth indicator

Battery level indicator

Charger

First charging the appliance

Allow the hand element to charge up for 16 hours before using it for the first time, until the harging indicator light turns to green.

Using the appliance

step 1

Align the front of the brush head with the front of the handle, ensure the gap between two parts.

step 2

Wet the toothbrush bristles and apply a small amount of toothpaste

step 3

Place the toothbrush bristles against the teeth at 45° angle towards the gumline

with the mode button select the desired brushing mode

Notice

Press the power on/off button for one time to turn on the appliance;

Press the button for 1second to stop the appliance;

Press once to switch mode while brushing your teeth

Cleaning brush head

Rinse the brush head under running water after each use;

To keep the toothbrush head clean, put the brush head cap on it after cleaning

Charging

Place the handle on the charger;

The flashing light of the battery level indicator shows that the toothbrush is charging. the power light will be on after full charge.

Data synchronizations

Open the APP, drop down the homepage and press "update the data". The bluetooth indicator light turns on when data are transmitted

Disposal

the battery must be removed from the appliance before it is scrapped, and hand it in at an official collection point for recycling.

If you can't remove the battery please send the product back to Soocare service center. We will help you to remove the battery from the appliance and deal it with in a environment-harmless way.

Technical data

model name: SOOCARE X3

rated voltage: 3.7V power consumption: 2w protection rating: IPX7 CMIIT ID: XXXXXX

standard: GB4706.1-2005,GB4706.59-2008

origin: China

model name: charging station

protection rating: IPX7

power supply: 100-240V-50/60Hz

power consumption: 3w

charging time: approx 20 hours

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- -- Reorient or relocate the receiving antenna.
- -- Increase the separation between the equipment and receiver.
- -- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- -- Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.