

Graphene intelligent physiotherapy bra

Use manual for NM- WX001

Thank you for the use of SHENZHEN T&M HEALTH TECHNOLOGY CO., LTD. produced physical therapy bra (full name: graphene intelligent physiotherapy bra)! Before use, please read this manual carefully and correctly grasp the method of using the intelligent physiotherapy bra, in order to make full use of its excellent performance, I believe it will give you health and enjoy the incomparable figure brought. Please keep the instruction manual for reference. If you encounter problems in the use of the process, please contact us or the company's after-sales service outlets or the company, we will be happy to serve you, thank you for your cooperation!

Product introduction:

Shenzhen's beautiful graphene intelligent physiotherapy bra, different from the traditional massage bra, specially with the market pain points developed "intelligent professional physiotherapy functional bra, combined with the mobile phone APP control, which mainly solves the problem of women's breasts: such as: breast hyperplasia, shaping, menstrual period pain, lactation pain, blood circulation and eliminate the masses. The prevention and improvement of breast disease in women is a great boon!

Product features:

The graphene unique far IR heating massage.

The key points of days Creek point simulation breast massage.

A variety of hot compress and massage mode with the body needs.

The exclusive built-in graphene and massage module process without trace.

The super circuit management, security worries.

Product specification parameter list:

Time to market: 2016

Brand: Song Yue Fei Ni

Name: physiotherapy bra

Executive standard: FZ/T 73012-2008

Material: Lycra, cotton, breathable cotton, nylon, spandex imports etc.;

Built in: graphene heating film and massage module

Main color: color, black, red, purple

The thickness index: moderate

Style: simple

Applicable season: Four Seasons

Origin: Shenzhen

Shelf life: 12 months

Product illustration:



Usage method:

- 1, Download APP for user registration, the mainstream download platform, such as 360, the application of treasure, Tencent, apple, etc..
- 2, wearing a bra will be connected to the power supply, and then open the phone APP for equipment connection.

3, for the first time, set the Bluetooth connection is successful, the next use, Bluetooth automatic connection.

4, according to the APP provided by the physical therapy program, choose their own options.

Cleaning and maintenance:

Good things need maintenance, please use neutral detergent cleaning (or bath) hand washing, can be used when necessary auxiliary toothbrush cleaning, washing can be placed in a ventilated place, proper maintenance can keep good comfort and adjustment effect. Graphene physiotherapy bra implantation of sophisticated technology, with high efficiency infrared light and nano silver adhesion to fibers, will avoid a bra to throw in the washing machine, but also avoid the bra on exposure in the sun, it will make micro molecular loss and diffusion and weaken the adjustment function, can not be dry after washing. Put in the shade to dry wear.

Matters needing attention:

1, it is recommended to use no more than 15 minutes a day.

2, breast-feeding women should not use vibration mode, should be pure hot compress.

3, half an hour after a meal should not be used.

4, in the use of the process if there is discomfort, please stop using.

After-sale service:

After sales service location, see warranty card.

1, the company provides a year warranty, the machine in the normal use of the process will be able to failure, consumers with warranty cards and invoices to our local maintenance points to get free repair services, except for human failure.

2, such as the emergence of the following circumstances within the warranty period, the company will provide you with free artificial service (free inspection fees, maintenance fees), only to charge replacement parts costs;

A, did not follow the instructions on the requirements of the operation and cause trouble; B, in case of natural disasters or

C, due to fall, collision caused by fault; D, self loading and repair and move into a fault;

E, due to injury accident, misuse, abuse, fault; cause F, commercial discoloration, abrasion, fouling overheating damage etc..

3, the warranty card is lost without replacement, please safekeeping.

4, free maintenance; please show the warranty card and the purchase invoice (such as the card after the modification or purchase did not fill, are invalid).

5, the user needs to repair, please contact the company at the local maintenance point, our professional maintenance staff will provide you with the enthusiastic service.

FCC Statement

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions: 1) this device may not cause harmful interference, and 2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

FCC Radiation Exposure Statement

This modular complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.