

Remedy Usage Guide

- ▶ Table of Contents
- ▶ Precautions for management
- ▶ Be aware!

01_ Check the function of the Ready program

02_ Check the components

03_ View Product Details

04_ Check how to use the cradle

- ▶ Use Ready

01_ Please check before using.

02_ Check the care mode before using

03_ Remedy Application Guide

- ▶ Preparation for Ready Care

01_Setting Voice Guide

02_ Setting a Language

03_Measuring Skin

- ▶ Know the care mode

01-High frequency

02_ Iontophoresis

03_EMS

04_LED

- ▶ Charging the battery
- ▶ Product Management
- ▶ Precautions for safe use
- ▶ Check before reporting a fault
- ▶ Product Warranty
- ▶ Consumer Dispute Resolution Standards
- ▶ Responsibility for paid service
- ▶ Appendix
- ▶ Product specification
- ▶ Precautions for management

Please make sure to check for safe use.

It is a face-only product, so it is not recommended to use it on other parts of the body, especially avoid thyroid areas.

After plastic surgery or skin-related procedures and treatment, it is recommended to use it while the skin has recovered after about 2 months.

However, if you are seriously suffering from a skin disease or are being treated, it is recommended to consult a specialist before using it.

Women who are pregnant or lactating may have skin problems due to hormonal effects, so we do not recommend using this device.

Regardless of whether the product is used or not, the skin condition can be temporarily changed due to environmental factors such as seasonal changes or physical factors such as drinking and stress. If the skin condition has changed, it is recommended to stop using it and return to normal or consult a specialist before using it.

In the case of pregnant women, long exposure to radio waves can affect the fetus, so it is recommended to avoid using them.

Depending on the skin condition of each individual, there may be skin irritation such as redness or tingling of the skin when using the product.

If irritation persists, stop using immediately and consult a dermatologist before using it.

If you have skin diseases or frequent dermatitis on your face, consult a specialist before using it.

It is not recommended to be used by people wearing pacemakers or transplanting devices.

A person who cannot clearly express his/her intention, including a child, is recommended to use it under the guidance of his/her guardian

I'm proud of you.

Don't use it for people with photosensitive diseases such as sun allergy.

Do not use the head of the product in one place, but move it according to the guide in this manual.

It is not recommended to use anything other than skin care.

Before using the product, remove the accessories worn on the area where the product touches.

Do not use this product for anything other than what is in the instruction manual.

Do not look directly at the LED light during operation. It can cause vision damage.

► Be aware!

01_ Check the function of the Ready program

This device is a product with high frequency/EMS composite technology for elastic care and iontophoresis/LED composite technology that promotes absorption of cosmetic ingredients.

High frequency

Technology to improve skin elasticity by generating deep heat inside the skin tissue

EMS

Technology that stimulates skin tissues and muscles with wavelengths similar to biocurrent to help improve elasticity

Iontophoresis

A technology that induces absorption of active ingredients by using electrical properties caused by potential differences in the skin as a function of effectively transmitting active ingredients in cosmetics

LED

Technology that boosts skin care through light therapy with different energy transmitted depending on the wavelength of light

Red (Play), Green (Whitening+Calming), Blue (Sterilization), Yellow (Playing+Whitening+Calming)

02_ Check the components

Please check the components before using it.

the main body

USB

Cable (Type C)

rests

instruction

03_ View Product Details

Clean up the Remedy Operation Button

Power button : Used to power on and off or change the mode.

LEVEL button: Used to adjust treatment strength.

LED/NEXT Button : Used to change LED color or to proceed with LED off function, next step.

Press for about 2 seconds to change the language in the order of Korean, English, and Chinese.

04_ Check how to use the cradle

You can charge the product by placing it on the stand.

A. Hold the product on the holder so that the operating button of the device is in front of you.

B. Place the device by aligning the charging wire hole of the holder with the charging port of the device.

C. Please charge the battery by connecting the USB cable and the charging adapter to the charging port on the back of the product.

► Use Ready

01_ Please check before using

precautions

Do not look directly at the LED light from a close distance during operation.

When wearing glasses and accessories, please remove them before using the device.

After washing your face clean, apply basic products in the order of skin-serum-cream (lotion) that suits your existing product or skin type.

When using EMS, it is good to apply vitamin serum when managing moisture gel and iontophoresis before using remedy.

Please remember that.

Remedy automatically provides suitable daily care by measuring the moisture content of your skin and the condition of your skin barrier and reflecting your menstrual cycle. Please take care of it according to the set care program, listen to the notification sound that notifies you for 3 seconds per part, and take care of the device and the face.

02_ Check the care mode before using

Customized care mode

Remedy offers a daily care program that fits your skin condition. Customized programs are set according to your skin barrier, skin moisture, and skin conditions by menstrual cycle. Customized care mode is recommended according to the instructions because the time of use and recommended treatment are different depending on the program.

Single care mode

This mode allows you to select tightening, lifting, and ion boosting. Tightening is a high frequency, lifting is an EMS, and ion boosting is an ion topophoresis function. If you select a mode that requires centralization, treatment proceeds. It takes less than 3 minutes and you can change the strength and color with the LEVEL button and LED button. Treatment proceeds by selecting the mode you want to focus on your skin.

App Choice Care Mode

This mode allows you to select the care you want.

First, turn on the Remedy and select the treatment that you are worried about or want to improve from the app and send it to me.

Then press the power button to select App Choice Care to proceed with the treatment.

Depending on the program, the usage time, intensity, and LED will be different, so we recommend that you use it according to the instructions.

03_ Remedy Application Guide

Please download the application for effective use of Geepiskin Remedy.

You can download it by searching Remedy on Android and App Store.

► Preparation for Ready Care

01_ Skin measurement

Press the power button to activate the ready to enter skin measurement mode.

Please measure it before washing your face. If you measure it after washing your face, apply the basic product and measure it 30 minutes later.

Place the measuring sensor of the device on the cheek for 10 seconds to measure the skin according to the voice guidance.

Depending on the skin condition measurement results, the customized care mode of the device is set.

Before starting the care, please wash your face thoroughly, including facial wastes and makeup.

Apply enough cream and lotion on your face and start the care according to the program's instructions.

02_Treatment

Apply enough serum, cream, and other products on your face.

Please proceed with the treatment according to the application's usage guide.

Keep that in mind

All care is guided by voice, and if care needs to be interrupted during use, press and hold the device's power button for at least 2 seconds to turn it off.

To restart, press power to start custom care from scratch.

► Care Mode Effects Precautions

01-Tightening

Care effect

It is a care that helps regenerate collagen and elastin by generating deep heat with high frequency care. It is recommended to be used in combination with anti-aging products as an effective solution for skin problems such as reduced elasticity and wrinkles.

02_ Lifting

Care effect

It is a treatment that helps skin elasticity by supplying microcurrents similar to biological currents of 40 to 60uA that flow through the human body through EMS care.

As aging progresses, you can stimulate the slowly decreasing biological current in the human body to activate the skin.

03_ Ion Boosting

Care effect

It is a technology that helps absorb cosmetics and solves skin problems through skin activity with Iontophoresis Care. It can be used in combination with various types of skincare products such as wrinkles, nutrition, whitening, and moisture to increase the absorption rate of active ingredients in cosmetics, which can be expected to be better than applying products in general.

04_LED

Care effect

Four basic optical therapy treatments that help with management and boosting.

Red - Playback

Green - Whitening + Soothing

Blue – Sterilization

Yellow - Playback + Whitening + Soothing

05_ Precautions for Care

- It is a treatment that generates deep heat, not the surface of the skin, and the skin heat may rise after use, but it is a natural phenomenon that occurs after care.
- Do not use treatment continuously for more than 3 minutes per day.
- Don't stay in one area for more than 3 seconds
- Please move your face according to the notification button sound of the device.
- Applying additional moisturizing cream to the care area after care helps to improve the care effect.
- Don't look directly at the light from the device during maintenance.

► Charging the battery

Charging

Please charge the battery before using the product.

- Connect the USB cable to the charging port of the device.
- Connect to the charging adapter using a USB cable.
- Charging begins when the charging adapter is connected to a power outlet.

Keep in mind when charging

- Be sure to use the charging cable provided.
- Disconnecting power after fully charging the battery saves energy.
- The product cannot be used while charging.
- You can charge the product by placing it on the stand. To use a cradle, refer to 'Check how to use the cradle'.
- Charging time and operating time may vary depending on battery conditions or operating conditions.

To check the battery status

Blue : If you are using the product

Red : When the product is being charged

Green : When the product is buffered

Charging time and operating time may vary depending on battery conditions or operating conditions.

► Manage Your Products

Manage the main body

Please remove any foreign substances such as cosmetics after the treatment is finished. Wipe the main body and head of the device lightly with a soft cloth. If you wipe it with a wet tissue, dry it in a cool place without moisture.

Strong chemical products such as alcohol, benzene, and acetone are prohibited. Be careful not to get water or cosmetics in the gaps in the product.

Please be careful for your safety.

The following "Careful for Safety" is to prevent unexpected risks or damage in advance, so please make sure to check it. Caution for safety is divided into two categories: Warning and Caution.

WARNING: Violation of instructions can lead to serious injury or death.

Caution: Violation of instructions may cause minor injury or damage to the product.

► Precautions for safe use

Associated with adapters when charging

Ensure that the charging adapter is firmly plugged into the outlet to prevent shaking and incomplete contact can cause fire. Make sure to use a 5V USB charging adapter with KC safety certification.

Stop charging immediately if there is smoke or a strange smell during charging. Please turn off

the power and remove the charging adapter and contact the customer service center.

If there is smoke during charging, ignoring it and continuing to charge it will cause fire and electric shock.

Do not damage, process, bend, twist, or heat the charging adapter and USB cable. Damage to the cable can cause fire and electric shock.

Warning

With the charging adapter connected to an outlet, do not insert conductors such as chopsticks into the outlet groove on the other side. Also, do not touch the charging adapter immediately after removing it from the outlet. It causes electric shock.

When you remove the charging adapter from the outlet, hold the charging adapter body and pull it out. It can cause damage or short circuit when you hold and pull the cable.

attention

Be sure to use the charging cable provided. Other charging cables may cause failure and cannot be compensated for repair.

Place on blankets, electric pad, or carpet and do not use or charge for long periods of time. Increased internal temperature can cause failure and fire.

Battery related

Internal batteries cannot be removed or replaced at your discretion. Also, do not dispose of batteries arbitrarily, but follow the disposal procedures in the area.

Don't shock the product, put it in water, store it in a car for a long time, or leave it unattended. There is a risk of explosion if the battery is not properly maintained.

attention

Do not store or use Remedy in hot or humid places (car rooms, jjimjilbangs, etc.) above 60 °C. There is a risk of explosion or fire. If the battery (charging terminal) leaks or stinks, be sure to contact the customer service center. There is a risk of explosion or fire.

Do not charge the battery using other power supplies or circuits. There is a risk of battery damage or fire.

If the battery is completely discharged and not used for a long time, the battery life may be reduced.

This product has a battery built into the system. When you use it for the first time after purchasing, charge the battery completely.

The battery is consumable and will gradually reduce the battery's usable time when used for a long time.

Internal batteries cannot be removed or replaced at your discretion, and free service is not available if the battery is removed at your discretion.

If you need to remove or replace the battery, please contact G Power Co., Ltd. The removal and replacement of the battery is covered by a paid service.

► To check before reporting a fault

Troubleshooting

If you experience the following symptoms while using the product, please check it again. It may not be a fault.

Symptom	Causes and Solutions	Symptom	Causes and Solutions	Symptom	Causes and Solutions
The product is not working properly.	<ul style="list-style-type: none">- Please check if you pressed the power button correctly.- Please charge the product according to the instruction in the instruction manual.- Make sure that the battery status indicator lights up when charging the product.- Make sure that the skin barrier sensor, skin moisture sensor, and treatment electrode head of the product are in close contact with	No power on.	<ul style="list-style-type: none">- Please charge the battery when the battery runs out.- Connect the charging adapter to the power supply.- Please check if the charger is working properly.- Please check if the voltage indicated on the product corresponds to the main voltage of the region.	I can't hear the audio	<ul style="list-style-type: none">- There is a speaker at the bottom of the front of the product. Hold the device and use it so that it doesn't cover it.

	the skin. Accurate measurements are made only when the skin moisture sensor and treatment electrode head are in correct contact with the skin.				
--	--	--	--	--	--

► a product warranty

Consumers have a valuable right to receive free service for one year for failures that occur while using the product.

This product can only be serviced free of charge within the warranty period in Korea.

Model name		Quality Assurance Period: 1 year
Date of purchase date		Internal Battery Quality Assurance Period: 6 months

※ When you request a repair, you need to provide this warranty with the date of purchase to receive the service, so please keep it well.

- Quality Assurance for this product is covered under warranty as stated in the warranty.
- The warranty period is calculated from the date of purchase, so please write down the date of purchase.
- If the date of purchase is not confirmed, the warranty period will be applied from 3 months after the date of manufacture.
- If a household product is used for business purposes or as a means of business, or in an abnormal environment, the warranty period is reduced by 50 percent.
- This warranty will not be reissued.
- In the case of a separate contract supply (procurement, Host, Single Goods, etc.), the warranty is applied in accordance with the main contract.

► Consumer Dispute Resolution Standards

Types of Consumer Damage		Compensation details	
		Within the warranty period	After Quality Assurance Period
	Important repairs required within 10 days of purchase	Replacement or refund of purchase price	paid repairs
	If the important components are that need be repaired within one month of purchase.	Product exchange or free repair	
	Replacement requires critical repair within one month	refund of purchase price	
	Operator is not possible.		
	In the event of a defect	a free repair	
	If the same defect has been repaired but the failure has recurred (3rd time)	Product exchange	
	In the case of a total of 4 repairs due to a failure in several areas, but the failure recurred (5th time)		
	In the event of damage during transportation or installation after purchase of the product		
	In case repair is not possible	Replacement or refund of purchase price	Refund by adding 10% of purchase price to the fixed depreciation amount
	If the product requested by the consumer is lost by the business operator or is not delivered for more than one month,		Refund by adding 10% to the fixed depreciation amount

► Liability for paid services

In the following cases, the consumer will have to pay for the warranty period when applying for A/S service. service.

- If a foreign object such as a toy or drink enters the inside of the product and causes a failure
- In the event of product damage or functional failure due to impact, such as dropping the product,
- Damage or deformation of product appearance due to customer carelessness or organic solvents such as thinner, benzene, etc
- If the product fails using non-genuine parts (please purchase genuine parts from GPOWER Co., Ltd.)
- Replace consumable parts that wear normally for use
- Product fault due to incorrect voltage selection
- If the customer arbitrarily dismantles the product and the accessories are lost or damaged
- Failure due to failure to comply with 'Care for Safety' within the product manual (Keep 'Care for Safety' will allow you to use the product for a long time.
- In the event of a failure due to natural disasters such as lightning, fire, wind and flood, salt damage, gas, earthquake, etc
- In the event of a failure due to other customer's fault
- For external reasons, such as a faulty power outlet itself or a failure to connect the power plug, other than a defect in the product itself
- Failure due to abnormal use of power supply or poor connection equipment;
- If you request a service with your own product, but you receive the service due to a third-party or other product defect
- In the event of a failure due to repair or alteration of the product by a person other than the official service point of G Power Co., Ltd
- If the service is requested due to a furniture defect, not a product defect
- Repairing defects caused by abnormal installation and alteration is a paid service, so please check the installation manual and installation checklist before installation.

If you have any questions, please contact GPOWER Co., Ltd. (T. 1588-3029, Kakao Talk @ gpskinHomeTherapy, info@GPOWER .kr). Contact the Customer Service Office (T. 1588-3029) to receive any suggestions or complaints.

The revised User Manual does not contain the FCC required warning statements.

a. Rule Part 15.19(a)(3): This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

b. Rule Part 15.21: The users manual or instruction manual for an intentional or unintentional radiator shall caution the user that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.


Note: This equipment has been tested and found to comply with the limits for a Class B digital

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help

► Product specification

General matters	
Device body	Product : GPSKIN REMED-
	Manufacturer : (주)지파워
	Input : 5 V, 750mA
Size (Length x Height x Width)	Device : 약 41.26 * 189.31 * 56.97 (mm)
	cradle : 약 81.9 mm x 81.9 mm x 68.43mm
Temperature used	5 °C to 35 °C
Humidity of use	0 % to 90 %
Recommended number of uses	Use 1-2 times a day (based on custom care mode) <ul style="list-style-type: none"> • It can be fully charged and used for approximately 60 minutes on a customized care basis. (Based on internal experimental results) • Hours of use may vary depending on the warranty period of the battery, ambient temperature, or conditions of use.
Full charge time	Within 3 hours Based on 5 V, 2.0 A charging adapter <ul style="list-style-type: none"> • Charging time may vary depending on charging adapter specification, remaining battery capacity, ambient temperature or charging conditions

KC certification		
Product	GP1I_H2101_V0R	
Product name (product name)	Specific small-power wireless equipment (Wireless equipment for wireless data communication systems)	
the trade name/manufacturer	GPOWER INC.,	
Country of Manufacture	Republic of Korea	
Identification code	R-R-G10-GP1I_H2101_V0R	

- Product images in this manual are intended to help you understand and may differ from the actual product. The contents indicated in this manual are subject to change without notice to improve the quality of the product.
- Refer to the label affixed on the outside or under the product for safety indications including the name of the product and the power supply
- This product is designed to open the inside of G Power Co., Ltd. It is dangerous for consumers to arbitrarily open or repair the cover of the product, so please refer it to the GPOWER Customer Center.
- Service parts are not included internally.

GPOWER Customer Center Representative

Phone Number T. 1588-3029

Kakao Talk @ gpskinHome Therapy

- If you know the exact product model name, fault condition, phone number, and address, you can get service faster.
- Please contact the above contact number when you receive any suggestions or complaints.

☐ Customer consultation time

- Weekdays 10:00-17:00 Closed on weekends and holidays