
For FAQs and more information, please visit:
www.eufylife.com

T9150 51005003825 V01

Smart Scale P3

Owner's Manual

English	01
Deutsch	12
Español	25
Français	38
Polska	51
Türkçe	64

Table of Contents

01	Important Safety Instructions	09	Maintenance
05	What's in the Box	09	Troubleshooting
06	Inserting the Batteries		When Measuring
06	Selecting the Measuring Unit	10	Specifications
07	Installing the EufyLife App and Pairing	11	Customer Service
08	Starting a Measurement		
08	Weight-Only Measurement		
08	Data Service		
	Third Party Services		
	Data Export		

Thank you for purchasing Smart Scale P3. Carefully read all the instructions below before using this device, and keep this manual for future reference.

Important Safety Instructions

- This device is powered by 4 AAA batteries. Do not attempt to recharge non-rechargeable batteries as this can cause the batteries to leak or explode.
- Do not mix and use different types of batteries together, e.g. lithium-ion (Li-ion) and nickel-metal hydride (NiMH) batteries.
- Do not mix new and used batteries.
- If the device is not used for an extended period of time, remove the batteries before storing the device in a cool and dry place.
- This device can be used by children aged from 8 years and above, and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of this appliance in a safe way and understand the hazards involved. Children should not play with this appliance. Cleaning and user maintenance should not be carried out by children without supervision.
- Used batteries should be placed in a sealed

plastic bag and disposed of safely according to local environmental regulations.

- Read all instructions before using the device. Keep the manual for future reference.
- Only use the device as described in this user manual. Any misuse may result in electric shocks, burns, fires, or cause other unexpected hazards.
- Be aware that this device is designed for individuals to measure their weight and to monitor their body composition. Any information provided by this device is in no way meant to treat, cure, or prevent diseases or illnesses from happening.
- DO NOT use this device if you are pregnant or if you suspect you may be pregnant. The effects of this device on a fetus are unknown and untested.
- DO NOT use this device if you use a wearable or implanted electronic device such as a defibrillator or pacemaker.
- This device should not be used by individuals who are acutely or chronically ill, suffering from diseases or taking medications that can affect

water levels. The accuracy of the readings for these individuals has not been verified. Specific medical advice should be obtained from a physician.

- The device is equipped with a data transmission function that emits electromagnetic energy. Portable Radio Frequency (RF) communications equipment (including peripherals such as antenna cables and external antennas) should not be used within 30 cm (12 inches) of the device, otherwise, performance may be affected.
- Do not modify this device.
- Do not expose this device to extreme temperatures, direct sunlight, moisture, or corrosive environments.
- To prevent slipping, do not step on the device after bathing or showering.
- This device should not be used when it is adjacent to or stacked on top of another device because it could cause the device to operate improperly. If such use is necessary, this device and the other device should be observed to verify that they are operating normally.

- If the batteries start to leak, put on your gloves and remove the leaking batteries from the battery compartment immediately. Place them in a clear plastic bag and recycle. If you don't wear gloves, you should wash your hands thoroughly with soap afterwards. In case of accidental contact with eyes, seek medical attention immediately.

Usage Tips:

- To ensure user profile accuracy, follow the instructions below when you start measurements:
- Place the scale on a flat, hard surface. Soft surfaces such as carpeting will affect the performance of the scale.
- Step on the platform with bare feet. Stand still until the measurement is complete.
- For the most consistent measurements, weigh yourself at the same time of day and location.
- Make sure your user profile (including birth date, gender, height, etc.) is accurate. Entering inaccurate information will affect the accuracy of the scale's measurements.
- Avoid weighing yourself immediately after

strenuous exercise, sitting in the sauna, taking a bath, or drinking and eating.

- The condition of the skin on the bottom of your feet can affect the measurement. The natural effects of aging or activity can harden the skin. Take your reading with clean, slightly damp feet for optimal accuracy.
- Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or overhydration due to factors such as alcohol consumption, menstruation, illness, intense exercise, etc.
- Considering the vast individual differences among children under the age of 16, the measurement of body fat or other body composition measurements may be less accurate. Their weight-only measurements will remain accurate.
- This scale is not intended for use by professional athletes.

Trademark Notice

- Apple, App Store and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries.
- Google Play, Google Home, Google Fit and the Google Play logo are trademarks of Google Inc.
- Fitbit is a registered trademark of Fitbit, Inc. in the U.S. and other countries.

Notice

This symbol on the device or its packaging indicates:
Do not dispose of electrical appliances as unsorted municipal waste; use separate collection facilities. Contact your local authority for information regarding the collection systems available. If electrical appliances are disposed of in landfills or dumps, hazardous substances can leak into the groundwater and get into the food chain, which are damaging to health and well-being. When replacing old appliances with new ones, the retailer is legally obliged to take back your old appliance for disposal at least free of charge.

FCC Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures: (1) Reorient or relocate the receiving antenna. (2) Increase the separation between the equipment and receiver. (3) Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. (4) Consult the dealer or an experienced radio / TV technician for help.

The following importer is the responsible party (for FCC matters only).

Company Name: POWER MOBILE LIFE, LLC
Address: 400 108th Ave NE Ste 400, Bellevue, WA 98004-5541
Telephone: +1 (800) 994 3056

ISED Compliance Statement

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1. This device may not cause interference.
2. This device must accept any interference, including interference that may cause undesired operation of the device.

Declaration of Conformity

This product complies with the radio interference requirements of the European Community. Hereby, Anker Innovations Limited declares that the radio equipment type T9150 is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following website: <https://support.eufy.com/s/article/Recommend?type=Download>

GB Declaration of Conformity

Hereby, Anker Innovations Limited declares that the product type T9150 is in compliance with Radio Equipment Regulations 2017 and The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012. The full text of the GB declaration of conformity is available at the following website: <https://support.eufy.com/s/article/Recommend?type=Download>

The following importer is the responsible party (contract for EU and UK matters)

Anker Technology (UK) Limited, GNR8, 49 Clarendon Road, Watford, Hertfordshire, WD17 1HP, United Kingdom

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Anker Innovations Limited

Room 1318-19, Hollywood Plaza, 610 Nathan Road, Mongkok, Kowloon, Hong Kong

Australia WARRANTY

This document details the terms and conditions of the warranty provided with this product. Please disregard any other representations contained on this product's packaging or provided with this product, which constitute or appear to constitute a warranty.

This warranty is offered by Anker Innovations Limited, Room 1318-19, Hollywood Plaza, 610 Nathan Road, Mongkok, Kowloon, Hong Kong, customer service support number (US) +1 (800) 988 7973 Mon-Fri 9:00am - 5:00pm (PT), (UK) +44 (0) 1604 936200 Mon-Fri 6:00 - 11:00 (GMT) and customer service email address support@anker.com.


Subject to the exclusions and limitations set out below, we warrant to you that if you purchase a product directly from us, or, from any of our authorized resellers, that it will not contain defects which are a result of faulty manufacturer workmanship or materials. This warranty subsists for 18 months from the date of purchase.

Subject to any rights you have under the Australian Consumer Law, we will not be liable for defects which are a result of damage or loss caused by factors beyond our control, any product which has not been maintained according to any care or maintenance instructions provided with the product, any alterations to a product which are not performed by us, or damage or defects caused to the product due to unusual or non-recommended use. We will not be liable for any special, indirect, consequential or economic loss or damage or loss of profits whatsoever suffered by you or any other person resulting from any act or omission by us. Our total liability in respect of this Warranty against Defects is limited to the replacement of the goods.

Any claim made under this warranty must be made by returning the defective goods to us, at your expense. Please provide details of the defect and evidence of purchase to support@anker.com and we will provide instructions on how to obtain a replacement. If you do not have Internet access please telephone us for further instructions.


This Warranty against Defects is provided in addition to other rights and remedies that you may have at law. Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

What's in the Box

Inserting the Batteries


- 1 Turn the scale over and place it on a flat surface.
- 2 Slide up to open the battery compartment door and insert 4 AAA batteries. Make sure the positive and negative ends are facing the correct polarity directions marked in the battery compartment.
- 3 Slide down to close the battery compartment door. Make sure the battery compartment door is fixed in place.

- **Low Battery**
 - will appear on the TFT display when battery is low. Replace all batteries at the same time. Do not mix new and used batteries.
 - When not measuring, the TFT display will turn off automatically after being idle for 10 seconds.

Selecting the Measuring Unit

- 1 Press the UNIT button at the bottom of the scale to cycle through and select the measuring unit: lb (pound) / kg (kilogram) / st (stone).


- To restore factory settings, press and hold the UNIT button for 10 seconds.


- 2 The measuring unit can be selected and saved within the EufyLife app. Open the EufyLife app and select the weight unit in the **Device Settings** menu.

Installing the EufyLife App and Pairing

- 1 Download the EufyLife app from the App Store (iOS devices) or Google Play (Android devices).

- 2 Enable Bluetooth and Wi-Fi on your smartphone.
- 3 Open the EufyLife app and create a user account. Please enter accurate information in your user profile.
- 4 Select **Add Device** on the **Me** page. Follow the in-app instructions to add the Smart Scale P3.

*Only with 2.4GHz Wi-Fi connection

- 5 Open the EufyLife app and follow the on-screen instructions to pair the scale with your device.

Once paired, the app and scale will automatically connect the next time the app is opened.


	Status Icon	Indicates
	Flashing	Pairing
	Solid	Connected
	Off	Not connected
	Flashing slowly	Ready for connecting
	Flashing quickly	Connecting
	Solid	Connected

- This scale is compatible with devices running on iOS 11.0 and above versions or Android 6.0 and above versions.
- To complete the setup, you need to allow the EufyLife app to use your location.
- Do not pair the scale from your smartphone's Bluetooth menu. Use the EufyLife app to pair the scale with your smartphone.

Starting a Measurement

- 1 Place the scale on a hard and flat surface.
- 2 When measuring for the first time, please make sure you have selected your user profile.
- 3 Step barefoot onto the scale.
 - An analysis of your weight and your body's composition will begin automatically.
 - After measuring, your weight and other available information will be shown on the scale's screen display. You can find your detailed body composition analysis within the EufyLife app.

- Step barefoot onto the scale. Stand still until the measurement is complete.
- The scale can automatically recognize users based on their previous stats.
- If the current user is not recognized, simply tap your feet to rotate between users.

Weight-Only Measurement

This scale can also be used as a conventional weight-reading scale.

- 1 Open the EufyLife app and tap **Device Settings**, then select **Simple Mode**.
- 2 Step on the scale to measure your current weight.

Data Service

Third Party Services

EufyLife App supports three third-party data service software: Apple Health, Google Fit, and Fitbit. You can open EufyLife App and find **Third Party Services** on your **Me** page to synchronize data.


Note: Only data of registered accounts can be synchronized. Data of family members will not be synchronized.

Data Export

EufyLife APP supports exporting of your data. You can open EufyLife App and find **Privacy and Data - Export All Data** on your **Me** page to export data.

Maintenance

Regularly clean the scale with a dry or slightly dampened soft cloth.

- Never wash the scale with water or immerse it in water.
- Never use detergents to clean the scale as this may cause it to malfunction or discolor.

Troubleshooting

Error Display	Description	Solution
	Overload. The device will turn off.	Stop using the scale for measurement.
	Low battery. The device will turn off.	Replace all batteries at the same time.
	Wi-Fi setup failure.	• Open the EufyLife app and follow the on-screen instructions to reconnect the scale with your device. • Move the scale to a place near your router before Wi-Fi setup. • Press and hold the UNIT button for 10 seconds to restore factory settings, then restart the Wi-Fi connection.
	Upgrade failed.	Open the EufyLife app and tap Device Settings , and then select Firmware Upgrade .

When Measuring

Issue	Cause	Solution
Abnormal measurement results: • Too high • Too low • Large difference between two consecutive measurements	Incorrect posture.	Step onto the scale barefoot and stand still.
	Inaccurate user profile.	Make sure the user profile is accurate.
	Device is placed on soft carpeting or uneven surface.	Place the device on a hard and flat surface.
TFT display does not turn on after stepping onto the scale.	Batteries not installed.	Install 4 new AAA batteries.
	Batteries are worn.	Replace all batteries at the same time.

Issue	Cause	Solution
Scale does not proceed to body composition analysis.	Taking a measurement while wearing shoes or socks.	Stand barefoot on the scale.
	Currently using the weight-only measuring mode in the EufyLife app.	Change the measuring mode in the EufyLife app.
	Incorrect posture.	Step onto the scale barefoot and stand still.
	Inaccurate user profile.	Make sure the user profile is accurate.

Specifications

Model	T9150
Power	6V (4× AAA batteries)
TFT Display	3.5 inch (320x 480pixel)
Measurement Units	Pound / Kilogram / Stone

Unit Increments	0.1lb / 0.05kg
Body Fat Unit Increments	0.1%
Body Water Unit Increments	0.1%
Muscle Mass Unit Increments	0.1kg
Measurement Range	3kg - 200kg / 6.6lb - 440lb / 0st:6.6lb - 31st:7lb
Activation Weight	5kg
Working Environment	Temperature: 10°C - 40°C / 50°F - 104°F Humidity: 20% - 90% RH
Storage Environment	Temperature: -10°C - 50°C / 14°F - 122°F Humidity: ≤90% RH
Operating Frequency Range	2400MHz - 2500MHz
Transmitted Power	Wi-Fi: 11.5 dBm Bluetooth: 2.5 dBm
Bluetooth Version	5.1
Wi-Fi Version	2.4G

Customer Service

For tutorial videos, FAQs, manuals, and more information, please visit:
support.eufylife.com.

Contact Us

✉ support@eufylife.com

… www.eufylife.com

☎ (US/CA) +1 (800) 994 3056
(UK) +44 (0) 1604 936200
(DE) +49 (0) 69 9579 7960
(日本) +81 03 4455 7823

For more support phone numbers,
please visit:
[https://support.eufylife.com/s/
phonecontactus](https://support.eufylife.com/s/phonecontactus)
Or scan the QR code below:

For tutorial videos, FAQs, manuals, and
more information, please visit:
<https://support.eufylife.com>
Or scan the QR code below:

f @EufyOfficial t @EufyOfficial @EufyOfficial p @EufyOfficial

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