
## Specification

**With App** - Push the button once before starting brushing, then it is connected to the App and starts to vibrate  
**Without App** - Push the button twice, then it starts to vibrate (Duration is 2 Mins, and it alerts every 30 Secs)

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
Product Name	BRUSH MONSTER Smart Toothbrush
Model	BMT100
Size	177x24x21mm
Net Weight	32g
Sonic Vibration	16,000 times/min
Power	AAA Dry battery x 1ea
Bluetooth	4.2 2.4 GHz
Material	ABS, PP, Silicone
Brush Material	DuPont Tynex® Nylon
Water Proof	IPX7
Using Age	Over 3years

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## Getting started with product and App

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


### Step 7. Taking pictures and collecting stickers

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**FCC Warning Statement**

FCC Part 15.19

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC Part 15.21

Any changes or modifications (including the antennas) to this device that are not expressly approved by the manufacturer may void the user's authority to operate the equipment.