

Superfine: User Manual

Meet your toothbrush

Replaceable sonic brush head
Power button + pressure sensor
Rechargeable handle
Magnetic charging cable
Stand + base
Travel case

Sonic brush head	Snaps into place, soft bristles pulse and oscillate. Only use quip supplied brush heads. Discontinue use if damaged.
Power button & LED Bar	Power on/off with one click
Pressure sensor	Pulses signal harsh brushing, helping protect teeth and gums.
Stand + base	Holds your quip in an elegant upright stance on any countertop or flat surface.
Travel case	Covers bristles and brush handle to protect on-the-go.
Charging cable	Magnetic end snaps onto the back of your brush for fast, convenient charging. Only use quip supplied cable. USB power adapter not included.

Get started

1. Charge up	Recycle packaging, then plug the charging cable into a power adapter (not included). To charge, snap the magnetic charging cradle to the back of the brush. Fully charged brush will have the entire LED bar light up for 3 times, and then turn off. Afterwards, tapping the light bar will light up the LEDs to indicate the current state of charge.
2. Set your intensity	To set up the intensity of your brush, pick up the brush to light up the LED bar, place your finger on the LED bead and slide up or down to adjust. Your brush will remember this for future use :)
3. Ready, set, brush	Wet bristles, add paste, then push the button to turn on brush. Gently guide the bristles from tooth to tooth, spending a few secs on the front, back, and chewing surfaces. No scrubbing needed , let the brush do the work.
4. Clean thoroughly	Spend 30 secs in each quadrant of your mouth, and switch to the next when your brush pulses once. Continue brushing until the vibrations stop at 2 mins.

*#quipTip: If the **light pulses**, loosen your grip – that's the pressure sensor helping you brush better.*

Lights to know

Keep an eye on the light! It communicates battery life, brush intensity, harsh pressure, and more.

Brushing intensities

The LED bead's location indicates the intensity

Power + battery

Fully charged full light bar

Charging light bar animates upwards, length of the bar indicates state of charge

Low battery light bar at lowest position pulsing 3 times

Brushing too hard light bar pulsing and animate downwards

Error full light bar pulsing

⚠ WARNING To reduce the risk of burns, electrocution, fire or injury to persons: Use only the charger provided, and specifically, an ETL or UL certified USB power supply, class III SELV output, 5V, 1A max. Discontinue use if this product or the power supply appears damaged in any way.

How to mount

Set up your magnetic mirror mount

1. Clean

Thoroughly clean and dry any smooth flat surface (like mirror, glass, or tile) where you'd like to attach your mount.

2. Prepare

Carefully peel off the thin plastic layer from the single-use adhesive on your mount.

3. Mount

Firmly press the mount onto its new home for 10–30 secs. Your brush will magnetically snap on just above the charging area.

#quiptip Learn how to easily remove your mount at getquip.com/help

Ready to refill?

1. Remove

Grip handle in one hand, brush head in the other and gently push up to unlock then rotate to remove.

2. Refresh

Use a damp cloth to scrub dried paste away from all sides of the handle and the top plastic tip.

3. Replace

Once all parts are clean and dry, slide the new brush head onto the handle, rotate and push until it clicks.

Model: S
FCC ID: 2AT6D-S

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference. and (2) this device must accept any interference received, including interference that may cause undesired operation.

Pursuant to FCC 15.21 of the FCC rules, changes not expressly approved by quip NYC, Inc. might cause harmful interference and void the FCC authorization to operate this product.

This product complies with the FCC RF radiation exposure limits set forth for an uncontrolled environment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an output on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Stay fresh & healthy

Refills auto-delivered

Get brush heads and more shipped every 3 months, like dentists recommend, at getquip.com/refills

Lifetime warranty

Your brush is covered for life as long as your refill plan is active (1 year if not).

Share the love

Refer friends or dentists to earn refill credits at getquip.com/referafriend

Got a quip question?

We're here to help you take good habits for a spin.
Get answers at getquip.com/help

[QR CODE]

Scan for a how-to video

[FB] [TWIT] [IG]

Follow us @getquip

Pro-fresh-ional tips

Round out your routine	Dentists recommend brushing 2 mins (twice a day!), flossing daily, plus regular cleanings with a professional.
Know when to swish	Brushing right after acidic foods or bevs can hurt tooth enamel. Swish water instead, and wait 30 mins to brush (or pop gum or mints).
Don't rinse right away	After brushing, spit your paste, but don't rinse! Wait 30 mins to let the formula work its magic. Repeat AM and PM.

Meet the family

Complete your routine with more oral care faves and fresh supplies shipped to your door.

Refillable Floss String
Anticavity Toothpaste
Cordless Water Flosser

[quip]

Healthy smiles start here