

We value every guest, so we will solve the after-sales problems encountered for you. We are eager to get your respect. If you are satisfied, please give us five points of praise. This is affirmation and support for our work! If you are not satisfied , Please contact us to solve, we cherish and value each of your suggestions.

Bluetooth

BodyScale
B1

Multifunctional household electronic scale

Multiple tests

Tempered glass

Smart App

Induction switch

HD Display

Bearing panel: High strength tempered glass weighing range :0.2~180 kg

Display Screen: LED HD Display Screen

Boot mode: automatic

Thank you for your trust, now start the wonderful weighing experience!

Connection scheme 1: By

A. use WeChat to sweep the left

B. search: at the top of the search window, Search for Flash Health

Usage:

1 WeChat Open Mini Program "Flash Health"

2 Authorized Login - Allow - Set Personal Data

Connection scheme 2: through APP software

Apple phone support IOS system

Support Android 4.3. above system

Scan the above QR code, first free download APP, open a beautiful weighing experience!

Bluetooth online weighing

1 When you need to understand the body's ingredients, you need to use the phone APP connect with the multifunctional body fat scale via Bluetooth devices. At this point you need to log in to the app market, in the search bar, enter "Youwo" in the app market to find the app, you can download and install to the mobile phone for free

2 For first use, enter your personal information (name, age, gender, height and target weight) After use without input, APP will analyze and store your measurement data, generate change curve for your reference

Measureme

4 Turn on Bluetooth, then turn on "You Wo" APP(Figure 4) if the display is not yet bound device, click on the binding (Figure 5), at this time, please step on your

The Bluetooth fat scale measures the body weight (Figure 7) and then the measurement appears

When you're done, click to view the data. will appear our measurement report!

have multiple measurements, combined with your weight signs information, APP will give reasonable exercise, diet and sleep guidance programs to assist you in managing weight, active activity, healthy eating, better sleep, and all-weather records to manage your health.

⑥ A maximum of 20 users can be recorded simultaneously in APP

Points fo

① When using barefoot, pay attention to the scale and feet do not touch water, so as not to slip and fall. The fuselage must be placed Horizontal hard ground, weighing as far as possible to make the body smooth, do not shake the body, keep static measurement. When measuring body composition, keep the measurement on barefoot station.

Each Bluetooth scale can only be connected to one mobile APP at the same time. Do not connect to multiple APP devices at the same time If the connection fails, restart the APP or see if the Bluetooth device is on, while reinstalling the battery of the scale. Then connect as normal;

For better weight monitoring, it is recommended to measure at the same time every day and get up in the morning. At the same time with the APP, you can record up to 8 people at the same time, so that you and your family can control the health weight at the same time;

After sal

① 7 days, normal use of machine problems, provide the card, return or replace the new machine;

② Within 365 days, the normal use of machine problems, enjoy free maintenance services, back and forth freight self-care; life-long

Product	Product Model:	Date of	Months	Day
Warranty matters	☐ 7 days return	☐ 1 year warranty	☐ Lifetime maintenance	

★★★★★

We value every guest, so we will solve the after-sales problems encountered for you. We are eager to get your respect. If you are satisfied, please give us five points of praise. This is affirmation and support for our work! If you are not satisfied , Please contact us to solve, we cherish and value each of your suggestions.

We have the right of final

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.