

SMART BLUETOOTH BODY FAT SCALE WARRANTY POLICY

① From the date of purchase, you can enjoy a year of replacement but not repair service.

② Within 15 days of purchase, the product quality problems (non-man-made damage), the merchants bear the back and forth freight, product quality for a year, within a year to provide new services.

Name		
Cellphone number		
Purchasing date		
Business name		
Product model		

Maintenance record	The fault phenomenon	
	Repair content	
	The date of repair	
	Maintenance personnel	

Product certificate

Place of Origin : Guangdong , China

inspector : _____The access time : _____

Guangdong Simad Industrial Co., Ltd.
Address : No.50 jinxing street, xiaojinkou street office, huicheng district, huizhou city, guangdong province,China
service hotline : +86-0752-6142656

Product specification

Size	27cm*25cm*2.2cm 28cm*26cm*2.2cm 28cm*28cm*2.2cm 28cm*24cm*2.5cm 29cm*26cm*2.2cm 30cm*29cm*2.5cm 30cm*30cm*2.5cm 31.2cm*31.2cm*2.5cm
Weight	1.1KG 1.2KG 1.3KG 1.5KG 1.7KG
Battery	2*AAA battery , 3*AAA battery , USB
Division value	0.1KG 0.01KG 0.05KG
Measurement Range	0.2-180kg
Display	Highlight LCD or LED
Material	ABS , Safety tempered glass

PS : Specific information for each model,According to the outer packing

display instructions

1. Lo (low battery)

3. O-Ld (Overweight prompt)

2. Err (test error)

4. unable to boot or display (replace the battery)

safety precautions

⊘

This product should not be used by those equipped with pacemaker, artificial lung or wearing electrocardiograph.

When using this product to measure body fat percentage, there is a weak electricity Circulation through the body may affect internal medical devices Normal operation, leading to serious consequences.

⊘

Do not use this product on slippery surfaces.

If the scale slips, it may hurt you.

Instructions for use

1 This product adopts high-precision sensor. In order to ensure the accuracy of measurement, please try to measure at the same time and under the same conditions and avoid the following situations:

☐ After strenuous exercise

☐ After bath or sauna

☐ After overeating

☐ After drink lots of water and wine

☐ Have a cold and fever

2 The recommended measurement time , keep at the same time andcondition every day.

6

7

8

10

12

14

16

18

breakfast

lunch

dinner

Measure before breakfast

Measure before lunch

Measure before dinner

3 The measurement needs to be barefoot, and keep the right and left feet respectively on the conductive plate, while the body is straight. Make sure the soles of your feet are clean when measuring. If you have dirt on the bottom of your feet or socks on your feet, it will not ensure the veracity of measurement data. As possible as measure in underwear only.

4 For children aged 0-6 years old, only weight and BMI data are provided. For teenagers aged 10-17 years old, fat rate can be detected. Because the under 18 age group is in the rapid growth stage,other data are for reference only.

5 Pregnant women, human dialysis patients, patients with swelling symptoms and the elderly over 80 years old will have different body fat rates than the general population, so the measured values are recommended as a reference.


SMART BLUETOOTH BODY SCALE

SPECIFICATION


Directions for use

The correct way to measure

- Please put the product on the hard flat floor, do not put on the carpet or other soft ground, or it will cause inaccurate measurement data.
- Clear before each measurement.
Reset method : After lightly stepping on the scale, the screen number will automatically return to 0.0KG or 0.00KG



The correct way to stand



APP download and measure

Step 1 : Install battery and APP download

1.APP Download


Scan the qr code on the right to complete the download, or search for "scale up" in the app store to download. (Chinese : 乐轻)

IOS can be downloaded from the apple store

Android can be downloaded on Google play




Step 2 : How to sign up(Mobile phone, email, third party registration and login)



APP download and measure


Step 3 : proper usage

- Mobile phone on bluetooth, enable location permissions
(Scale up APP Permissions need to allow location)
- Make sure the scale is placed on a flat , hard floor
- Step on the scale before the display of 0.0kg, and then measure
- Stand up straight with your bare feet on the scales and your feet aligned with the electrodes





view data , Sharing system , historical data

- After the measurement is completed, the report can be directly clicked to view each data standard, and the PDF format report can be generated.
- Reports can be Shared to third party software, community mobile phone.
- View historical data: click the upper right trend icon, and then click the historical data to view all measurement records



APP download and measure

Support holding baby to measure weight



FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

RF Exposure Information

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.