

Model:TZ280


TENSWALL

Smart Body Fat Scale

User Manual

www.tenswall.com

English



Warning

If it is strictly forbidden to disassemble, squeeze, puncture, or damage the scale. Do not expose the scale to high temperatures or exposed to the environment with a temperature higher than 50°C. Thank you for choosing our smart body fat scale products. Please handle and use it carefully.

Import: Measurements should only be used for references, not for medical diagnosis.

Note:

Not intended for use by those people:

- Infants and children

- Pregnant women

- Elderly people

Measuring:


1. It is easy to use when using the body fat scale. In wet conditions, please dry the scale with a dry cloth. Turn the scale on and stand on the scale to be informed of the ground, and then use the scale to obtain accurate data.

Maintenance:

When cleaning, please wipe the surface of the scale with a clean, dry cloth.

- Do not use any cleaning agent for surface cleaning.

Product Information



1. Install HWIFI APP

Method 1: Search "Tenswall" on Google Play or Apple App Store to download and install the application.
Method 2: Enter the URL on the mobile browser, or scan the QR code to download the Google Play or Apple App Store download page.
URL: <https://www.tenswall.com/app/appdown.html>








2. Turn on the phone's Bluetooth

After HWIFI APP is installed on the mobile phone with iOS system, turn on the phone's Bluetooth. Turn on the scale and turn on the mobile phone's Bluetooth. Turn on the mobile phone's Bluetooth and connect to the scale. After connecting, the mobile phone's Bluetooth connection is established. Please turn on the mobile phone's Bluetooth.

Notes:

1. Please clean the body fat scale to be kept dry, please use a clean, dry cloth when cleaning. Please use a dry cloth to clean the scale. 2. Two batteries are included. Please install the batteries according to the position and direction of the "plus" and "minus" side.



TENSWALL

TENSWALL

FCC Requirement

changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20cm between the radiator & your body. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.