

▲ Statement
To reduce the risk of malfunction/electric shock/inju-ry/fire/death and equipment damage, please read all the instructions and warnings in this manual before using this product

I. Users who have received tooth orthodontic, filling or plantation should consult densits price to use: 2. This product is designed for booth, gum and longue cleaning. Do not use it for any other purpose, if you feel unwell or painful, stop using the product and consult a doctor

### \* CAUTION

- Do not immerse the charging dock in water or use it in any liquid: any liquid:

  2 Males sure that the plug on the charging dock is completely dry before connecting the dock to the power supply

  3. The product is equipped with rechargeable battery. Do not shrow it into fire or put, use or recharge it in a high-temperature environment.

  4. Do not hit or expt to product. If it is damaged because million or the product is an analysis of the product of the product of the product.

- 1.Mark sure that the product is not damaged prior to use. Stop using it if it is damaged or falls heavily. 2 bn ont use a damaged charging dock;
  3 be not use the product if the handle, charging dock or brush head is damaged.
  4 be not use the charging dock in high temperatures such as under direct suright.
- 5.Do not use the toothbrush during a bath or shower.

## \* REMINDER:

- \* REMINDER:

  1. Do not share one brush head with others;
  2. Do not way, the power last on the changing dock when
  2. Do not way, the power last on the changing dock when
  3. Do not put metal deployers such as coins and clips on the
  changing dock.

  2. Do not put metal deployers such as coins and clips on the
  changing dock.

  3. De not put metal deployers such as coins and clips on the
  change of the coins of the policies of all publicies of all publicies

  5. Pleases stop using debremed or novel brustless. Only use
  touch heads recommended by the name stop to the
  head every three months. If there are signs of weat,
  regions the bush head in advinces.

# MODE DESCRIPTION

Brush head Brush head x 2 protection cap

Eaznul

PRODUCT DESCRIPTION

Clearing .....

Watering ....

Sensitely ...

0

Normal operation, the bottom blue light indicator is always on, the mode LED is always on;

During charging, the bottom LED flashes once every second; After battery is fully charged, the bottom indicator is off, and it enters sleep state; The battery level is less than 30%, and the bottom red light indicator is always on;

# The charging process takes at least 12 hours.

Search EaznuF in the app store or scan the following OR code to download the APP.





Install the brush head on the handle and make sure that there is a space between the two parts. (The space can prevent the head from colliding with the handle when it vibrates quickly and can improve the motor's power. Because stairs carl easily hide in the apace, we suggest that you disassemble and clean the brush head regularly.)



Wet the brush head and apply toothpaste on the brush head other. STEP 3

1.Tap the APP, turn on the Bluetooth, and press the toothbrush switch button according to the screen display of the APP competition (flower APP guides to start brushing.

2. Guide the brush head slowly from tooth to both with a 45 degree angel between bristles and gums line.

3. The toothbrush may sivibate frequently during use to remind you that you should brush another place.









正面



Power OnliPower Offillionde Switch Power onli United Text Good statlant office of the Control sight offilipress and reliesable the bottom, will Power offilibring power of, short press to puese, it will be control of puese, it will be control to puese, it will be control to puese the control puese of the control pueses power about and hold for 1.5 seconds been release to dark throwing, and the puese and hold for 1.5 seconds it is existent and hold for 1.5 seconds it is existent plays water to see prevailing, just press and hold for 1.5 seconds it is existent puese and hold for 1.5 seconds it is existent puese and hold for 1.5 seconds it is existent puese and hold for 1.5 seconds it is existent puese and hold for 1.5 seconds it is existent pueses and hold for 1.5 s

Whitening mode:Solve teeth whitening problems, inhibit dental plaque, remove surface dirt, whitening your teeth.

Gum care mode: Easily absorb the force of brushing, optimize deep clean effect and gentle care your gums

Sensitive mode:High cleaning performance and gentle brushing for sensitive teeth and ours.

## CLEANING AND MAINTENANCE

1.After brushing, rinse your brush head under running water. 2. After cleaning , cover the protective cap to keep the tooth-brush

4.After cleaning, wipe the handle and charging dock with a dry cloth and place in a dry place to ensure that the electric toothbrush is kept in its best condition.



Download the Eaznuf application, please press the electric toothbrush button for 0.5 seconds, and then turn the Bluetooth on of the phone to complete data synchronization



# 0 - m 0

Main page, toothbrush pag

Product contains rechargeable lithium battery ,which may be containinating. Please take battery out before you dispose it or send it to a recycling center appointed by local government and send the battery to a designated waste collection station. If you cannot take the battery out, you can send the product to our after-sales service center. We will remove the battery and dispose it in a nervicommatal-filendly manner.

Handle
Name of Produc

Name of Product: Eazural Smart sonic electric toothbrush Model:T0315 Rated voltage: DC 3.7V Rated Power: 2 Watt Charging time: about 12 hours Waterproof grade:IPX7 Made in China

\* Please use the standard USB adapter(at least 5V/1A) Travel advice: different regions have different power supply voltages, please choose a USB adapter that suits you

# NOTICES

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits

Standard: CONFORMS TO UL STD. 1431 CERTIFIED TO CSA STD. C22.2 NO. 60335-1 and C22.2 NO. 60335-2-52



