

98mm

98mm

98mm

98mm

FCC compliance statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

-- Reorient or relocate the receiving antenna.

-- Increase the separation between the equipment and receiver.

-- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

-- Consult the dealer or an experienced radio/TV technician for help.

FCC Radiation Exposure statement

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

131mm

طريقة الإستعمال

- امسح الرمز و هم يتنزل التطبيق
- قم بإشياء حسب المسخدم
- قم بتزويد البطاريات بالميزان
- أدخل إلى إضافة جهاز جديد اختار اسم الجهاز ليتم التوصيل
- قم على الميزان
- توأمم الجهاز بالتطبيق بطريقة أوتوماتيك
- قم على الجهاز مرة لمرة حتي تترى الأرقام تومض بشكل ثابت
- استمتع بفراشة نتججه الفحص
- عندما يمشي الرقم على الميزان مرتين يعني إنتهاء عملية القياس
- عندما يثبث الرقم 0 علي شاشة الميزان يعني تم الإنتهاء من عملية

تحذير

- قد لا تكون قياسات تكوينك الجسدي دقيقة عند استخدام الميزان مباشرة بعد ممارسة الأنشطة البدنية الشاقة
- هذا الطراز من ميزان Qutek الذكي غير مخصص لاستخدامه من قبل الأطفال دون الثامنة عشرة من العمر
- يجب عدم استعمال ميزان Qutek الذكي لتشخيص أي أوضاع صحية أو علاجه بل استشر طبيبك أو اختصاصي الصحة دائماً لتحديد ما هو أنسب لك

موانع الإستخدام

- لا تستخدمه إذا كنت حاملاً
- لا تستخدمه على الأرضيات الزلقة، لا تستخدمه عندما يكون جسدك مبللاً
- لا تستخدمه إذا كان لديك منظم للقلب أو آلة تنطيط كهربية القلب أو أي جهاز طبي داخلي

- Measure weight correctly: place the scale flat on hard floor and reset it before measurement. The weight should display 0.0 lb. before stepping on the scale.



- Measure body fat and other body composition parameters correctly: Stand straight on scale, and have both feet uniformly touching the electrodes as the below figure.



Specifications

Dimensions: 12.2" x 12.2" x 1.0" (310mm x 310mm x 25mm)
Weight: 3.97 lb. (1.8 kg)
Batteries: 4 x 1.5V AAA
Measurement Range: 0.4 ~ 400 lb. (0.2 ~ 180 kg)
Measurement Resolution: 0.2 lb. (0.1 kg)
Model# QSP511
Freq. Band: 2402MHz-2480 MHz
Max output power: 2.63dBm (Declaration for EU Compliance)

Trouble shooting

6800 kg Normal working	No display Change battery	6800 kg No wireless connection
Lo Change battery	EEr Overload	Cal Calibrating



Steps Body Composition scale User Manual

Qutek Wellbeing Introduces STEPS Body Composition Scale, designed to add a huge insight for a real health experience right at your feet.

This scale uses bioelectrical impedance analysis (BIA) technology to estimate key body health stats including body fat percentage to deliver actionable insight. With STEPS Body Composition Scale, you can now weigh yourself and your body composition data will sync automatically with Qutek App.

STEPS gives precise results to help you set your goal to stay on track while connecting with your body. Thank you for choosing STEPS Body Composition Scale to be your health partner!

Guarantee

The scale has been inspected for accuracy prior to leaving the manufacturing line, and comes with a 1 year guarantee against any defects in either workmanship or materials. **Apply warranty terms and conditions.

Address

Qutek Technology Limited
71-75 Shelton Street, Covent Garden, London, WC2H 9JQ, United Kingdom

Re#R200715



- STEPS must not be used on slippery floor, and in case your body is wet.
- STEPS is not recommended for pregnant women, patients with a pacemaker or artificial lung or wear an ECG or any other device.
- In case of malfunction, check the batteries in the scale and replace if needed.
- In case the device is damaged, do not use in order to avoid any inconvenient results.
- In case you do not plan on using the product for a week or more, please remove the battery.



- Your body composition data measurements may not be accurate when using the scale immediately after heavy physical activities.
- This model of STEPS body scale should not be used to diagnose or treat any medical conditions. It is best to consult with your physician for any advice on medications, diet or exercise programs.
- STEPS is not intended to be used by children under the age of 18.
- STEPS is not waterproof, never immerse the device in water. To clean, use a damp cloth. Do not use soap or any detergents.

- Use STEPS on a hard surface during use.
- Gently step on the scale to wake up and wait at least 10 seconds before weighing yourself. It is recommended never to jump or stamp to avoid damage.
- When not in use, please remove the batteries and place your device in a dry room avoiding heat, sunlight or dampness.

Power On Scale

Your Steps Body Composition Scale is shipped with 4 AAA batteries. Install the batteries inside the battery compartment and power on your scale.

After power on, you can measure body weight. However, to measure other body composition parameters such as body fat percentage, you need to use your Steps Body Composition Scale together with your smart phones or tablets. You will need to install mobile applications programs on your mobile devices first. See next.

Apps Install and Measurement

The Qutek Wireless Body Composition Smart Scales can connect wirelessly to mobile phones or tablets. To do so mobile devices need to install mobile application programs. Easiest way to download the mobile app is to use your phone to scan the below QR code, and then follow the instructions given.

Step 1: Download the APP



Scan QR code and download 'FeelFit' APPs



Register



Install battery cells

Step 2: Pairing scale with mobile phone



Enter my device/add new device/pair device



Step on scale



Pairing scale to APPs

Step 3: Measure and check the body analysis report



Step on again and keep to digital number stable



Check the body report in APPs



Number flashes twice means weight test finishes. Number 0 rolling to stable means body fat test finishes.

What's Included

Your Steps Body Composition Scale Model box includes:

- ① Steps Body Composition Scale
- ② Batteries : 4x 1.5V AAA
- ③ User Manual

技术要求