

## Instructionmanual



- ①PowerButton
- ②HeatingButton
- ③VibrationButton
- ④ChargingIndicator
- ⑤Light
- ⑥ModeButton
- ⑦ChargingInterface

Name	OperationMethod	IndicatorLight Status
Mode Button	Shortpresstocyclemodes, longprestodisablemode function	Blue(Soothing)/Purple(Vitality)/Red(Deep)
Heating Button	Shortpresstocycle temperatures, longpresto disableheating	Blue(41°C)/Purple(43°C)/Red(45°C)
Vibration Button	Shortpresstocycle intensities	Blue(Low)/Purple(Medium)/Red(High)
Power Button	Longpress3stoppower on/off,shortpresto togglelight	Whenthedeviceis turnedon,allobutton indicatorlightsare on, andthedefault colorisblue

- o BatteryCapacity:2000mAh
- o BatteryLife:120minutes
- o HeatingPower:5W
- o ProductDimensions:240\*170\*107mm
- o ChargingTime:Approximately4hours
- o ProductMaterial:PUFoam
- o ChargingInterface:Type-C

### 1.PowerOn/Off

- LongpressthePowerButtonfor3secondstoturnon(lights, andheatingenabledbydefault).
- Shortpresstotoggelighton/off
- LongpressthePowerButtonfor3secondstoturnoff(indicator lightsfullyextinguish)

ShortpresstheVibrationButtontocyclethrough:

- Low(blueindicatorlight)
- Medium(purpleindicatorlight)
- High(redindicatorlight)

\*Applicability: Only affects Soothing/Vitality Modes

### 3.ModeSwitching

Lightly press the Mode Button to cycle through:

- Shortpressonthemodebuttonturnsonthevibrationfunctionand  
alongpressforthreecondsturnsoffthevibration.
  - SoothingMode(gentlecontinuousvibration|blueindicatorlight)
  - VitalityMode(intermittentpulsedvibration|purpleindicator  
light)
  - DeepMode(gradualintensityincreasefromweaktostrong|red  
indicatorlight)
- \*SpecialNote:DeepModeignoresvibrationlevelsettings**

#### 4.HeatingAdjustment

ShortpresstheHeatingButtontocyclethrough:

- Low(41°C|blueindicatorlight)
- Medium(43°C|purpleindicatorlight)
- High(45°C|redindicatorlight)

## 5.ChargingInstructions

Use the provided Type-C cable with a 5V/1A adapter.

ChargingStatus:

- Red charging indicator light stay on → Charging in progress
- Indicator light turn off → Fully charged

## 6.LowBatteryReminder

RedlightflashesActionrequired:Chargeimmediatel

## 7.RemoteControlOperation

After manual power-on, use the remote control to:

- o Adjust settings (modes, heating, vibration)
  - o Power off the device
- \*Note: Remote cannot be used for power-on

## Precautions

1. Recommended to turn off device during charging disconnect power immediately after full charge
2. Avoid liquid contact with charging interface Clean with dry cloth only (do not submerge in water)
3. Do not wipe the product body with alcohol.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The equipment has been tested and found to comply with the limits. The device can be used in portable exposure condition without restriction.