GRAVITY

Eye Massager

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
Soothe facial tension instantly with the Gravity Eye Massager.

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Primeros pasos

Operating Instruction

1. Long press the key for 3 seconds to start the Multi-Function mode (pressure + vibration + hot compress)
2. Press the key again to start Relax mode (air pressure + heat compress)
3. Continue to press the key to start Rest mode (air pressure)
4. Continue to press the key to start Revive mode (pressure + vibration)
5. Continue to press the key and start to Heat Mode (hot compress)
6. Press the button for 3 seconds to turn off the machine function
7. Bluetooth will automatically turn on when the device starts
Device ID: YBY

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Portable et sans fil

Portátil e inalámbrico

Type-c charging mode; button operation

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Remote heating micro computer chip control technology

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Air bag massage (massage the muscles around the eye area by proper matching between air filling and air releasing)

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Power
4.5W

Input
DC 5V

5V2A adapter with international plugs

Battery Capacity
1200mAh

Battery Life
1 hour

Modes
Multi-Function Mode
Relax Mode
Rest Mode
Revive Mode
Heat Mode

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Resolución de problemas

Unable to boot

1. If the battery power is too low, charge immediately
2. Connect the charger, charge, and allow the host to automatically remove the protection (charge for 5 seconds)
3. The battery is not charged and the charger is damaged
4. The pressure intensity is too weak

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Remove contacts before use

If there is any dirt or foreign matter in the mask, wipe it away with a wet cloth before using.

Do not use while charging

Do not use during outdoor walking or driving

Do not use in hot or humid places

Do not use if you have been diagnosed with glaucoma, cataract, retinal detachment or any major eye-related disease

For patients with deep stroke, hypertension and hypotension, please consult a medical professional before using

Check the massage area frequently

If you feel any discomfort during use, cease use immediately

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This product should not be used by children or incapacitated persons without adult supervision

Use this product only according to instructions

THIS PRODUCT DOES NOT CONSTITUTE MEDICAL TREATMENT AND IS NOT A SUBSTITUTE FOR A MEDICAL EXAMINATION OR DIAGNOSIS. IF YOU ARE DEALING WITH A SERIOUS HEALTH CONDITION CHECK WITH YOUR HEALTH CARE PROVIDER BEFORE USING.

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FCC Warning

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

1. This device may not cause harmful interference
2. This device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance

with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Note: The Grantee is not responsible for any changes or modifications not expressly approved by the party responsible for compliance. Such modifications could void the user's authority to operate the equipment. The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

IC Warning

Avertissement CI

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions:

1. This device may not cause interference

2. This device must accept any interference, including interference that may cause undesired operation of the device.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

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About Gravity

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Acerca de Gravity

At Gravity, we believe in wellbeing for everyone. We meet you where you are on your wellness journey, and help you optimize your life with the latest in wellness innovation. Build your toolkit of Gravity essentials, and get ready to meet your best self.

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