

140mm

108mm

108mm

108mm

Specifications:

Display:All datas displayed when power on
Capacity:3-180KGS
Division:1D=100G/0.2LB
Unit:KG/LB(ST/LB(For LCD Display) or KG/LB(For LED Display)
Low Battery Indication
Overload Indication:Overload more than 2.5KGS,total 182.5KGS displayed"Err"
Best Working Temperature:-10~60℃
Moisure Environment:≤90%RH

Instructions for LEAONE Smart Scale APP.

- Instantly weighing**
Open the LEAONE APP and your Bluetooth, Enter into "Dashboard" interface. Stand on the scale and find "No paired" device, then click"to pair". After weiging finished, datas will be synchronized to APP.
If the Bluetooth disconnected, the weight datas will be saved into memory of the scale, then synchronized to APP while the Bluetooth connected.
 - Type weight datas into LEAONE Smart Scale APP**
The user can type weight datas into APP to keep a track
 - Dashboard**
Dashboard interface displays weight information of the user and family members, and generates the report and give health advice based on the measurement. Including Dashboard interface, Family Member, FAQ, and Settings.
 - Weight Trend**
Including Weekly report, Monthly report and Quarterly report, indicating the changes of weight through curves.
- Notice

For Bluetooth Body Fat Scale, there are 12 datas analyzed: Fat ,Water, Visceral Fat, Muscle, Bone, BMR, BMI, Weight, Protein, Obesity, Body Age, LBM.

5

Management of scale

Management of scale including two parts: scale management and connection management

Scale management

Manage the process of pairing, bind and unbind

Connection management

Including Bluetooth status indication and Bluetooth connection

Account and user

Including account management and user management

Account Management

Including registration, login in, change password, forget password, links to Facebook, Twitter...etc

User management

Manage user information, including Nickname, Birthday, Height, Gender...etc,also managing family member account and information (under the same registration account).

Warnings

- Please keep the scale away from wet floor, don't use the scale if the glass is wet in case of slip.
- Please stand on the middle of the scale in case of turnover.
- The platform of the scale is made of glass, please avoid sharp or harsh things scratching the surface, excessive forces or sudden shocks, unexpected fallen-off will break the glass and injure the user.
- The Bluetooth scale is not waterproof. Keep the surrounding environment dry, do not submerge it in water.
- Use soft cloth with neutral cleaning liquids to clean the surface, keep away contacting with organic solvent nor chemicals.
- When the batteries are in low power status, please remove the batteries in case of leaking liquids corrode the scale.
- This product should not be used by individuals with a pacemaker or other implanted medical devices as it passes a low level electrical signal through the body which may interfere with its operation. This product is also not suitable for women during pregnancy!

Notices

- Please stand on the scale in the right direction.
- Please stand on the scale no less than 10 seconds if you measure body fat.
- Please confirm Gender,Height,Weight datas are correct.

6


Bluetooth Smart Scale



正面

Step 1 Check Required Specifications for Phone or Tablet

- Before installing
① Before installing Lealone APP, kindly check if the phone or tablet meets following requirement:



Compatible with IOS7.0 and above with Bluetooth 4.0 and above


Phone 4s, iPhone 5, iPhone 5c, iPhone 5s
iPhone 6, iPhone 6 Plus, iPod touch 5, iPad 3
iPad 4, iPad Air, iPad Air 2, iPad mini
iPad mini 2, iPad mini 3



Compatible with Android 4.3 and above with Bluetooth 4.0 and above

SAMSUNG: galaxy s5, galaxy s4, galaxy s3, galaxy note4, galaxy note3, galaxy note2, galaxy grand2
XIAOMI: XIAOMI Note, XIAOMI4, XIAOMI3, XIAOMI2/2s, XIAOMI tablet, HONGMI Note 4G, HONGMI 1s
HUAWEI: Ascend P8, Ascend P7, Honor 6Plus, Honor 6, mate7, Ascend G7
If you are using Android newer than 5.1, Please meet the request to allow location setting.

Step 2 Install the batteries



For Battery Scales

- Open the cover of battery
- Install the batteries in the right direction
- Close the cover



For USB Charger Scales

Or make the scale recharged with the USB charger before electricity used up

2

Step 3 Install APP

- Switch on Bluetooth
Use your Devices to scan QR Code or search Lealone from APP Store/Google Play to install APP.

- Phone or Tablet connects to Net (WiFi, 3G or 4G)



For iOS



For Android

Step 4 Pair the devices, register and complete user information



3

Step 5 Weighing



Step 6 View the test result and sharing the datas



4

FCCStatement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
--Reorient or relocate the receiving antenna.
--Increase the separation between the equipment and receiver.
--Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
--Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

portable device statements
The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

反面

2018-11-16

三折页