

<p>Specification</p> <p>Capacity: 0.44lbs-396lbs/ 0.2kg-180kg Unit: lb/kg Voltage: 3V Power: 3*AAA batteries Low battery indication: display "Lo" and power off Overload indication: display "OL" Best working temperature: 5-35°C Storage temperature: -10-60°C</p> <p>AiLink App</p> <p>Instant weighing Open the AiLink app and turn Bluetooth on. The auto-pair process will complete automatically. Begin weighing - the weight data will synchronise with the app. If Bluetooth is disconnected, the data will be saved into the scale's memory and then synced to the app when Bluetooth reconnects.</p> <p>Current weight The current weight interface displays the weight information of the user and generates reports with health tips based on the recorded measurements.</p> <p>Functions Records body weight, fat, hydration, muscle and bone percentages, BMI & more.</p>	<p>About After-sales Service</p> <p>We value your complete satisfaction and strive to offer a high quality product and excellent customer service. If you have any problems with the product, feel free to contact us via Amazon with your order number or the following email address: vipservice2021@126.com</p> <p>Individuals with implanted pacemakers, artificial lungs, or wearable ECG devices should not use this product! When using this product to measure body fat percentage, a weak electrical current will pass through the body, which may interfere with the normal operation of medical devices and cause serious consequences.</p> <p>FCC Warning Statement Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference</p>	<p>will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:</p> <ul style="list-style-type: none">- Reorient or relocate the receiving antenna.- Increase the separation between the equipment and receiver.- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.- Consult the dealer or an experienced radio/TV technician for help. <p>This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.</p> <p>NOTE: This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter.</p> <p>The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.</p>	<div><div>Smart Body Scale</div><div><div><div>Bluetooth</div><div>Download on the App Store</div><div>GET IT ON Google Play</div></div></div></div>	<p>Step 1 Install AiLink App</p> <p>1. Switch on Bluetooth. 2. Connect to WiFi or 3G/4G/5G network. 3. Use your device (tablet or phone) to scan the QR code or search AiLink from AppStore/Google Play and install the app.</p> <div></div> <p>Step 2 Install batteries</p> <p>1. Open the battery cover. 2. Install the batteries as indicated. 3. Close the battery cover.</p> <div></div>	<p>Step 3 Register & Complete user information</p> <div><div>Open AiLink app</div><div>Register</div><div>Complete personal information</div></div> <div></div> <p>Step 4 Pair the Smart Body Scale</p> <div><div>Add device</div><div>Devices nearby</div><div>Stand on the scale barefoot & synchronize your data with the app</div></div> <div></div>	<p>Step 5 View the results</p> <div><div>Measuring</div><div>Data history</div></div> <div></div> <p>Warning</p> <p>1. Please keep the scale away from wet floor areas. Don't use the scale if the glass on top of the scale is wet to avoid injury. 2. Always stand in the middle of the scale. 3. The top of the scale is made of glass, be careful not to scratch the surface. 4. The scale is not waterproof. Keep it dry and never submerge it in water or any other liquid. 5. Use a soft cloth to clean the surface and never use harsh chemicals. 6. Data shown by the Smart Body Scale is for reference only and should not be used for medical purposes.</p>
---	--	---	--	---	--	--