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Instructions for use of Bluetooth body fat scale

Please refer to the instructions in detail before use. Bluetooth body fat scale supports two systems, IOS and Android versions.

▲ Precautions

- The after-sales package in weighing scale is one month (except for human factors such as being affected with damp, broken and self-disassembled).
- Just wait for the number to return to zero, and you can start measuring.
- Don't be close to the water source and fire source, and don't use chemical cleaner for cleaning, so as to avoid product failure or damage.
- Please keep it out of the reach of children.

▲ APP link description

- Download the [OKOK] APP on the mobile phone and install and register (scan the code to download or download it in the download application software that comes with the mobile phone).
- Turn on Bluetooth, then enter the APP, search for devices in my device, and stand on the scale to keep the screen on, so you can link and bind.

▲ operating instructions

- Stand on the scale with your feet steady.
- Let it stand for about 5 seconds and wait for the number to stabilize.
- Measurement data is for reference only.
- The weighing range is 0.2–180kg.

▲ Charging instructions (please ignore this item for battery models)

- Please charge it for 1–2 hours for the first time.
- The charging time is limited to 2 hours.
- Please use an ordinary mobile phone plug. It is not recommended to use the fast charging head of 2A. The battery symbol will jump when charging, and it will stop beating when it is full.

▲ Weigh items above 0.2kg

(Please press the button on the back to switch the required units first)

- Place weighing scale on a hard and flat ground.
- Put pressure on weighing scale to make its screen light up.
- Wait for the number to return to zero, then put the item to display the number, and the light item will automatically return to zero after 2 flashes.

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FCC Warning:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution: Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 0cm between the radiator and your body.