

Product Features



Automatic Timing



One Key Operation

Temperature-sensitive
heat pack

Air Pressure Massage



Vibration Relaxation




Intelligent BT

Product Diagram

1. On/Off Button/Volume Off Button/Mode Switch Button


2. Operating Status Indicator

3. Charging Interface



Product Configuration

The Host, Charging Line, Instruction



Function Introduction

1. long press

- Please long press “⏻” this key for 2 seconds to turn on/off the machine. The local music is turned on by default when the machine is turned on, and the local BT is turned on.
- The BT name of this product is: "Eye Massager".
- For the first use, you need to open the BT of your phone, search for "Eye Massager" and click to pair the device for the first time.
- For each subsequent use, the BT device successfully matched with this product will automatically connect when the mobile phone BT is turned on and no other devices are connected. After the BT connection, the voice prompts "BT has been connected".

2.short press

- Please short press ⏻ this key to switch modes. The default mode is "Automatic Mode".
- When pressing the key briefly, it can switch to the next mode automatically (with corresponding voice prompts).

There are 5 massage modes:

Mode 1: Automatic Mode	Airbag+Heating (strong)+Vibration (variable frequency)+15-minute Music Cycle
Mode 2: Sleep Mode	Heating (strong)+15-minute Default Mute
Mode 3: Clear Mode	Airbag+Heating (middle)+15-minute Music Cycle
Mode 4: Validity Mode	Airbag (strong)+Heating (strong)+Vibration (constant frequency)+15-minute Music Cycle
Mode 5: Eye Production Mode for Colored Contacts	Airbag+Heating (strong)+10-minute Music Cycle

3.Double click

Double click ⏻ this key to turn off the music and mute it directly.

FCC Warning Statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- - Reorient or relocate the receiving antenna.
- - Increase the separation between the equipment and receiver.
- - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- - Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.