



STATEMENT

Thank you for purchasing our Nano Resin Microneedling Pen. Before use, please carefully read this instruction and follow the safety precautions herein to prevent damage to the equipment or prevent the apparatus from being used normally.

LU&BA owns the copyright of all contents herein. Without our written permission, no one is allowed to copy, duplicate, or translate this instruction into other languages.

Proprietary information protected by copyright law contained herein shall be kept confidential by users and shall not be disclosed to unrelated third parties.

The Company is entitled to revise the instruction and change the product description at any time, and the final interpretation right is owned by the Company.

FOREWORD

Thank you for purchasing this product. This product is carefully made by us, combining unique technologies and designs, and meeting high-quality standards.

We are proud to meet the demands of our customers as much as possible. We welcome your comments and suggestions and believe that keeping in touch with our customers is very important for our future product line. We encourage you to read through the instructions to fully understand the operating characteristics of the instrument.

INSTRUCTION

1. Please read this instruction carefully before use.
2. Before using this apparatus, please clean the skin, wash hands, and disinfect.
3. During charging, DO NOT power on for operation or touch the charger with wet hands.
4. When placing the apparatus, it is required to remove the needles.
5. After using the apparatus, please wipe it and place it in a dry and ventilated area.
6. DO NOT open the device yourself; otherwise, LU&BA will not be responsible for subsequent operation or any warranty claims.

IoT STARTUP PROCESS

1. Download the LU&BA app on your mobile phone.
2. Enter your email to log in.
3. Add the customer's name in the upper left corner.
4. Enter the operation page and select the corresponding LU&BA instrument.
5. Press the Start button to scan the product consumable code.
6. One customer, one code.

CAUTIONS

The nano resin needles are disposable; secondary use is not allowed to prevent skin problems caused by cross-infection.

Please place this apparatus in a location out of the reach of children.

OPERATING INSTRUCTIONS

1. Long press the power switch for three seconds to turn on the indicator light.
2. Clean the skin and insert the collagen solution into the mouth of the instrument.
3. Open the APP and connect Bluetooth successfully.
4. Scan the code to boot.
5. Long press to power on; short press to adjust the gears, the equipment can be adjusted in 6 levels.
6. After finishing the operation, long press to power off and remove the needles.

APPARATUS PARAMETERS

Model No.:JY2024-10-01
Input:5V/0.5A
Energy Levels:6energy levels
Mode of Action: mechanical vibration and introduction
Power Supply Mode: USB charging

WARRANTY

The warranty is valid for two years from the date of purchase if the product has material, manufacturing, or functional defects not caused by the user.

The warranty is invalid in the following cases:

- The device is damaged due to improper operation, modification, or maintenance by unauthorized third parties.
- The device is damaged during transportation, such as from impact or improper handling.
- Damage caused by moisture or other liquids, foreign objects entering the device, or the use of unapproved and discouraged accessories contrary to the instructions.
- Use of accessories that do not comply with or match the specifications during the use of the device.
- Damage or loss of part or all of the device due to force majeure or accidents caused by natural disasters (such as earthquakes, floods, fires, and lightning, etc.).
- Surface damage, including but not limited to scratches and dents.

For the above situations, if maintenance is needed, users are required to bear the corresponding expenses.

DJM Medical Instrument Germany GmbH

Siemensring 98, 47877 Willich, Germany

Tel: 0049 2154 8864488

E-mail: info@djmiot.com

www.djmiot.com

Hong Kong Office of DJM:

Room 803, 8/F Block A, Hong Kong Industrial Centre, No 489-491 Castle Peak Road, Cheung Sha Wan, HK

Tel: (852)28136070

Email: operation@djmhk.com

wechat : DJM_HongKong


facebook: DJM HK

website: www.djm-group.com

Mainland China Office of DJM

Add: No.9 Yuexing Five Road, Yuehai Sireet, High-tech Zone Community, Nanshan District, Shenzhen City.

E-mail.: market@dimiot.com



This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Radiation Exposure statement

The device has been evaluatec to meel general RF exposure requirement. The device can be used in portable exposure condition without restriction.