# **PUSH™**

**Quick Start Guide**

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
**HM2411WB  
PUSH2.0**

**InBody**

---

## 1. Appearance and Part Names

## 2. How to Download the App

Easily track and manage your child's growth using your smartphone.

Download the InBody hi App from Google Play or the App Store.


< Android OS >


< iOS >


## 3. Precautions When Measuring

Adjust the length of the device according to the user's needs before measuring. This can be used if the sensor signal is disrupted due to high Percent Body Fat in the abdominal area.

Ultrasonic waves travel through air and tend to spread sideways. Clear the surrounding area to prevent interference from ultrasonic waves caused by objects like toys, and maintain an upright posture.


Align the Head Position with the end of the stadiometer.  
Measure the height based on the Head Position, even if slightly tilted.

## 4. Measuring with the App


- 1) Install InBody hi App and open it.
- 2) Press 'hi' button in the bottom center of the home in InBody hi App.
- 3) Perform the measurement according to the guidance of the InBody hi App screen.  
\* Refer to '5.Measuring Height' for the height measurement method.

## 5. Measuring Height

- 1) Press the button to turn on the device.  
The indoor temperature will be displayed on the screen.
- 2) Lean against the wall, place the device on your head, and attach the support vertically to the wall. Ensure the bottom of the device is facing downward.

- 3) Press the button to measure the distance from the ultrasonic sensor to the floor.
- 4) Press the button again to remeasure.

## 6. Error Messages

The following error messages may appear on the device screen during use.

### No Ultrasonic Signal Is Detected

Keep the area free of carpets, blankets, or other obstacles.

### Battery Is Low

The device powers off automatically after a certain time.

### Measured Length $\leq 5.0$ cm (1.97 in) or Sensor Obstructed

Measurement cannot be performed.

### Measured Length $\geq 230.0$ cm (90.55 in)

Measurement cannot be performed.

## 7. Setting the Unit

- 1) Press and hold the button for 10 seconds when the device is powered off.
- 2) When the unit setting screen appears, press the button to select the desired unit.

(°C / cm)

(°F / cm)

(°C / ft in)

(°F / ft in)

- 3) Press and hold the button for 3 seconds to confirm the selected unit.

## 8. Precautions When Replacing Batteries

- 1) Use AAA alkaline or manganese batteries.
- 2) Check if the positive (+) and negative (-) terminals of the battery are aligned correctly. Batteries may leak or generate heat if inserted incorrectly.
- 3) Do not mix old and new batteries. Replace all 3 AAA batteries with those of the same brand and model.
- 4) Do not handle batteries with wet hands.

## 9. Precautions for Use

- 1) Do not throw or drop the device, as it may cause injury. Store it in a safe place.
- 2) Dropping the device or applying excessive force to the ultrasonic sensor may cause damage.
- 3) Avoid measuring height in suddenly hot environments or near an air conditioner. Measure height in an area with a stable indoor temperature.
- 4) Measuring on carpets or blankets may absorb ultrasonic waves and cause inaccurate measurements.
- 5) Individuals with a high Percent Body Fat may experience inaccurate measurements due to abdominal interference.
- 6) The ultrasonic sensor used in this device has a spread angle of approximately 30 degrees. Height may not be measured accurately if there are objects nearby reflecting ultrasonic waves. Tidy up the surroundings.
- 7) Batteries are hazardous if they come into contact with the mouth.
- 8) Keep the device away from water.
- 9) InBody's free warranty service does not apply if the device is water-damaged or disassembled without authorization.
- 10) Remove the batteries when not using the device for an extended period.
- 11) This device includes one year of free service from the date of purchase. Damage and malfunctions caused by customer misuse are not covered under the warranty.

## 10. Specifications

Device Name	Portable Stadiometer
Measurement Items	Height, Temperature
Height Measurement Accuracy	±1 cm (±0.39 in)
Measurement Time	1 - 3 seconds
Current Consumption	Below 100 mA
Display Screen	LED
Device Size	65(W) × 65(L) × 300(H) mm 2.56(W) × 2.56(L) × 11.81(H) in
Operation Environment	41°F - 104°C (5 - 40°F), Up to 2,000 m, Max. 80 % RH up to 88°C (31°F), decreasing linearly to 50 % RH at 104°C (40°F)
Model Name	PUSH2.0
Height Measurement Range	5.0 - 230.0 cm (1.97 - 90.55 in)
Temperature Measurement Range	Indoor Temperature 5 - 40 °C (±1 °C), (41 - 104 °F (±1.8 °F))
Power Supply	DC 4.5 V (3 x 1.5 V AAA batteries)
Sensor	Ultrasonic Sensor, Temperature Sensor
Weight	203 g (0.45 lbs)
Storage Environment	-20 - 60 °C (-4 - 140 °F)

- \* The appearance and specifications of the device are subject to change without prior notice to improve performance.

## 11. Device Warranty

### Device Warranty Standard

Free warranty service will be provided only for a manufacturer's defect or a spontaneous failure within the warranty period.

In the Republic of Korea, any matters not specified in this document will follow the customer dispute resolution guidelines set by the Korea Fair Trade Commission. However, for users in other countries, the applicable local regulations will apply. For malfunctions or damage, please contact our customer service for assistance.

- \* For detailed warranty information, please refer to the InBody hi App.

Operating Instructions

Manufacturer

Authorized  
representative in the  
EUROPEAN COMMUNITY

European Conformity

Serial Number

WEEE mark

---

The devices have been tested against the relevant harmonized/designated standards and are in conformity with the essential requirements and other relevant requirements of the Radio Equipment Directive (RED) (2014/53/EU) and of the Radio Equipment Regulations (RER) (S.I. 2017/1206).

The devices are entitled to carry the CE and UKCA Marks, and a formal Declaration of Conformity (DoC) is available at the product manufacture web page.


This device was tested in compliance with CE RF Exposure limit.

(0 cm distance from human body (limit 4 W/kg @ 10g SAR)

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---

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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