

Me Health Updated Dielines

Info Card

- Card stock
- Double-sided print.

PMS 7464C + 40% white

PMS 7464C


PMS 178C

PMS 2965C

~ 42mm

~ 73 mm

Front




Power

welcome guide
see reverse

~ 42mm

~ 102 mm

Front




Hub

welcome guide
see reverse

Back

~ 73 mm

~ 42mm




Set Up:

- 1 Remove components from packaging.
- 2 Connect power cord to the back of the Hub.
- 3 Place Hub in a location with good cell reception and plug in.

Back

~ 102 mm

~ 42mm



Network Connection:
Pulse White = Start Up / Searching

Flashing Green = Communicating

Solid Green = Connected / Success

Red = Error (move to new location)

Device Connection:
Pulse White = Start Up / Searching

Flashing Green = Communicating


Solid Green = Connected / Success

Red = Error (move to new location)

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Hub Outer sticker

This will wrap around the back. Use clear sticker to seal front opening.



*Specs on back UPC and QR Codes/Numbers are all placeholders.