

Company Name: Jiangsu Yuyue Medical Equipment & Supply Co., Ltd.
Add: No.1 Baisheng Road Development Zone, Danyang, Jiangsu province, 212300
China
TEL: +86 18252963550
FAX: +86-511-86900827
E-mail: qa.yuwell@yuyue.com.cn

DECLARATION

Date: 2025-03-03

We, Jiangsu Yuyue Medical Equipment & Supply Co., Ltd. hereby declare that product name: Pulse Oximeter (FCC ID: 2A2JJ-YX500; models: YX500M, YX500S, YX500L, YX501S, YX501M, YX501L) are electrically identical with the same electromagnetic emissions and electromagnetic compatibility characteristics. And the YX500M has been tested by Bay Area Compliance Laboratory Corp. (Kunshan).

The differences of them are as follows:

Model	Differences	Power on Mode
YX500S		
YX500M		
YX500L	The hardware structure is entirely identical, with the only variation being the size differences in the ring dimensions. The software for these three models is exactly the same, except for the displayed model names.	Touch button for power on.
YX501S		
YX501M		
YX501L	The hardware structure is entirely identical, with the only variation being the size differences in the ring dimensions. The software for these three models is exactly the same, except for the displayed model names.	Automatic power on.
The YX500 and YX501 series share an identical hardware structure. The key distinction between them lies in their boot-up mechanisms, which are software-implemented. Specifically, the YX500 series initiates the boot		

process by detecting interrupts from touch buttons in its software. In contrast, the YX501 series triggers startup through software recognition of interrupts from a three-axis accelerometer.

Please contact me if there is need for any additional clarification or information.

Best Regards.

Signature:

A handwritten signature in black ink that reads "Caixia Yang". The signature is written in a cursive style with "Caixia" on the top line and "Yang" on the bottom line.

Printed Name: Caixia Yang

Title: RA Engineer